
# Kurzanleitung · Quick Start Guide Mode d'emploi condensé · Istruzioni in breve Instrucciones básicas VITA Linearguide 3D-MASTER®

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E

### Instructions for use

Dear Customer,  
Congratulations on your purchase of the new  
VITA Linearguide 3D-MASTER!

With the VITA Linearguide 3D-MASTER you  
can determine the correct tooth shade swiftly  
and accurately. It includes the 29 shade tabs  
of the VITA SYSTEM 3D-MASTER.

### Assortment structure:

The VITA Linearguide 3D-MASTER consists of  
an acrylic box (protective housing) containing  
a VITA Valueguide 3D-MASTER and five VITA  
Chroma/Hueguides 3D-MASTER. The VITA  
Valueguide 3D-MASTER, which is located  
on top, contains 6 shade samples of the  
lightness levels 0 – 5, with an intermediate  
(neutral) hue M and the chroma level 2: 0M2,  
1M2, 2M2, 3M2, 4M2, 5M2.

The 5 VITA Chroma/Hueguides 3D-MASTER  
positioned vertically in the lower part of the  
acrylic box contain the 29 shade samples as  
found in the VITA SYSTEM 3D-MASTER.

They are arranged for the 5 VITA  
Chroma/Hueguides 3D-MASTER in the  
following way:

1. Lightness level 0/1: 0M1, 0M2, 0M3, 1M1,  
1M2
2. Lightness level 2: 2M1, 2L1.5, 2R1.5, 2M2,  
2L2.5, 2R2.5, 2M3
3. Lightness level 3: 3M1, 3L1.5, 3R1.5, 3M2,  
3L2.5, 3R2.5, 3M3
4. Lightness level 4: 4M1, 4L1.5, 4R1.5, 4M2,  
4L2.5, 4R2.5, 4M3
5. Lightness level 5: 5M1, 5M2, 5M3

### Procedure:

- Open the acrylic box by pulling apart each  
end.
- First take the VITA Valueguide 3D-MASTER  
out of the opened acrylic box
- Make a first preselection with the VITA  
Valueguide 3D-MASTER by comparing the  
shade tabs with the natural tooth. Here  
you determine which lightness level group  
commencing from 0 to 5 best corresponds  
to the tooth shade.
- Based to your first preselection take the  
corresponding VITA Chroma/Hueguide  
3D-MASTER (0/1, 2, 3, 4 or 5) out of the  
acrylic box and determine the chroma and  
hue.
- You have now determined the correct tooth  
shade quickly and reliably.

### Tips for shade-taking:

- Determine tooth shades, if possible, in  
daylight conditions or under standardized  
daylight lamps or daylight corrected

4500K – 5500K and not under conventional  
room lighting.

- All shade selections should be determined  
swiftly, typically accepting your first  
decision since the eye tires after  
approximately 5 – 7 seconds.
- Office surroundings are preferably  
furnished in neutral colors. Patients should  
be advised to remove influencing factors,  
e.g. lipstick or cosmetics and cover up  
brightly colored clothing with a neutral  
grey bib.

### Hygiene tips and instructions for care:

The acrylic parts of the 5 VITA  
Chroma/Hueguides 3D-MASTER and the  
interior part of the box are made of  
a high-quality, easy-care material that  
is temperature-resistant (up to 180 °C).  
The entire VITA Valueguide 3D-MASTER,  
VITA Chroma/Hueguides 3D-MASTER and  
the interior part of the acrylic box can be  
autoclaved at a temperature of up to max.  
140 °C.

### Important:

The two exterior parts of the acrylic box  
cannot be autoclaved! For the sterilization  
of the shade guides and the interior part,  
the acrylic box must be taken apart as follows:

- Pull the two halves of the acrylic box past  
the retaining mechanism and remove the  
interior part.
- In order to reassemble, simply push the  
two halves of the acrylic box back together  
over the interior part and past the retaining  
mechanism.

### Disinfection:

The acrylic parts of the VITA Linearguide  
3D-MASTER can be surface disinfected.  
We recommend a customary, mild spray  
disinfectant or disinfectant wipes. According  
to the manufacturers' product information  
these disinfectants are effective against  
HIV and HBV, and are available from most  
dental dealers. Please heed the respective  
manufacturer's instructions for use.  
Disinfectants containing phenol, compounds  
with phenylphenol groups or methyl ethyl  
ketone will damage the acrylic parts.

### Note:

The acrylic parts of the VITA Linearguide  
3D-MASTER are neither monomer- nor  
UV-resistant, and must be protected from  
intense sunlight.

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