

## Bond strength of self-adhesive luting composites after chemical curing

### a) Materials and methods

The respective restorative material was embedded in a temporary crown and bridge material and ground flat. According to the respective working instructions, the surface of the restorative material was etched with VITA ADIVA CERA-ETCH (5 % hydrofluoric acid gel, HF) or sandblasted with Al<sub>2</sub>O<sub>3</sub> (50 µm, 2 bar) and then the appropriate primer was applied.

Light-curing composite Harvard Restore (Harvard Dental International) was used to produce cylindrical test specimen (Ø approx. 4 mm, h max. 5 mm). The respective luting composite was applied to the base of the cylinder and the cylinder was placed on the prepared surface of the restorative material (12 mm x 14 mm).


The samples were stored in a dry place at room temperature (23°C) for 10 minutes and then immersed in water at 37°C for 24 hours.

In the shear bond strength test with a mechanical testing machine and a feed rate of 0.5 mm/min, the maximum force was determined until the test cylinder sheared off and the bond strength relating to the bonding area was also calculated. Each value is based on five measurements.

### b) Source

Harvard Dental International GmbH, Dr. Dierk Lübbers, Scientific Affairs Manager, Hoppegarten, Report: 11/2016

### c) Result



### d) Conclusion

The self-adhesive luting composite VITA ADIVA S-CEM demonstrated excellent adhesion to all restorative materials tested.