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VITA VM.LC PRIMER I+II

Haftvermittler

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VITA VM.LC PRIMER I+II

Bonding agent

Please read the instructions for use before using the product. VITA VM.LC PRIMER I+II are bonding agents for reliable bonding of veneering composites to dental substructure materials.

Intended purpose

VITA VM.LC PRIMER is an accessory of the VITA VM.LC veneering resins. VITA VM.LC products are veneering resins for dental restorations.

Patient target group

No restrictions

Intended user

Dental professionals only: dentists and dental technicians (Rx only)

Indications

- Adhesive bonding of metal substructure surfaces (non-precious metal alloys, PM alloys, titanium) with (meth-)acrylate-based resins/composites
- Adhesive bonding of ZrO₂ substructure surfaces with (meth-)acrylate-based resins/composites
- Adhesive bonding of cross-linked composites/resins such as VITA CAD Temp or of high-performance polymers (PEEK, PEKK) with (meth-)acrylate-based resins/composites

Contraindications

- Applications that are not listed under indications

Product use / application

VITA VM.LC PRIMER I:

Add a sufficient quantity of VITA VM.LC PRIMER I to a mixing dish and apply to the corresponding surface with a disposable brush or instrument and allow to air dry (approx. 60s). If not all of the intended surface has been coated, repeat the above step.

VITA VM.LC PRIMER II:

Use a disposable brush to apply a uniform thin layer of VITA VM.LC PRIMER II. Then polymerize with an approved light-curing unit according to VITA VM.LC polymerization instructions (see www.vita-zahnfabrik.com/vmlc). After light curing, the reactive (sticky) surface must be coated with opaque material and veneering composite without any contamination or long interruptions.

Storage/disposal

They can be disposed of with household waste. The products labelled with a pictogram for hazardous substances are to be disposed of as hazardous waste. Recyclable waste (such as attachments, paper and plastics) must be disposed of using appropriate recycling systems. If necessary, contaminated product residues should be pretreated in accordance with regional regulations and disposed of separately. Storage temperature: Not above 28°C / 82°F

Composition

Primer I:

Acetone 93-99%, adhesive monomers < 5%

Primer II

Methyl methacrylate 40-80%, bifunctional methacrylates 20-40% Photoinitiator < 3%

Safety information

- While work is in progress, wear suitable safety goggles/face protection, gloves and protective clothing.
- In case of formation of dust, an extraction system or a dust protection mask must be used.

Product reliability

Information on reporting serious incidents in connection with medical devices, general risks associated with dental treatments, residual risks and (if applicable) short clinical safety and performance reports (SSCPs) can be found at https://www.vita-zahnfabrik.com/product_safety The safety data sheets can be downloaded at www.vita-zahnfabrik.com or requested by fax at (+49) 7761-562-233.

Recommended combination product

Sub-structure material	Pretreatment (unless otherwise described by the framework manufacturer)	Use		
		only VMLC Primer I	First step: VMLC Primer I Second step: VMLC Primer II	only VMLC Primer II

Non-precious metal alloys	Thoroughly sandblast with 110-250 µm Al ₂ O ₃ disposable abrasive material at 2.5 -3.5 bar. After sandblasting, clean with oil-free compressed air.	+	++	-
Precious metal alloys	For high gold alloys, use retentions and prepare with milling tools. Thoroughly sandblast with 110-250 µm Al ₂ O ₃ disposable abrasive material at 2.5 -3.5 bar. After sandblasting, clean with oil-free compressed air.	-	++	-
Titanium alloys	Thoroughly sandblast with 50 µm Al ₂ O ₃ disposable abrasive material at max. 2 bar. After sandblasting, clean with alcohol or steam and then dry with oil-free compressed air.	+	++	-
Zirconia (e.g., VITA YZ...)	Thoroughly sandblast with 50 µm Al ₂ O ₃ disposable abrasive material at max. 2 bar. After sandblasting, clean in an ultrasonic bath and then dry with oil-free compressed air.	+	++	-
PMMA (e.g., VITA CAD-Temp...)	Thoroughly sandblast with 50 µm Al ₂ O ₃ disposable abrasive material at max. 2 bar. After sandblasting, clean and dry if necessary.	-	++	+
High performance polymers (e.g., PEEK, PEKK...)	Thoroughly sandblast with 50-110 µm Al ₂ O ₃ disposable abrasive material at 2-3 bar. After sandblasting, clean with alcohol or oil-free compressed air.	-	++	+

++ highly recommended

+ recommended

- not recommended / not indicated

Disclaimer

Please note: Our products must be used in accordance with the instructions for use. We do not accept any liability for damage resulting from incorrect handling or usage. The user is furthermore obliged to check the product before use with regard to its suitability for the intended area of applications. We cannot accept any liability if the product is used in conjunction with materials or equipment from other manufacturers that are not compatible or not authorized for use with our product and this results in damage. The VITA Modulbox is not necessarily a component of the product.

Date of issue of these instructions for use: 10-2022

After the publication of these instructions for use, any previous versions become obsolete. The current version can be found at www.vita-zahnfabrik.com

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Highly flammable liquid and vapour. / Causes skin irritation / Causes serious eye irritation / May cause drowsiness or dizziness

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Highly flammable liquid and vapour / Causes skin irritation / May cause an allergic skin reaction / May cause respiratory irritation / Harmful to aquatic life with long lasting effects

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VITA VM.LC PRIMER I+II

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Instructions for use | short version

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