## VITA POWDER SCAN SPRAY

DE

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EN

### Statement of intended use

VITA Powder Scan Spray is a spray-on, titanium dioxide free and blue coloured pigment suspension with mint aroma for intraoral application (tooth surface) as well as for extraoral use (plaster die / plaster model) for the fabrication of the opto-electronic impression for CAD/CAM restorations. With the VITA Powder Scan Spray, the tooth or model surface is matted in order to obtain an optimum image with the intraoral camera (e. g. with Sirona CEREC 3 or CEREC AC) and/or with the scanner.

VITA Powder Scan Spray makes it possible to achieve a homogeneous masking layer with precise reproduction of the contours. The blue shading facilitates the checking of the application. After scanning, VITA Powder Scan Spray can be easily removed from the tooth surface or the model using a water jet.

### Composition

Titanium dioxide free pigment suspension in ethanol. Isobutane as a propellant.

### Area of applications

#### • Intraoral use

As a contrast powder in order to obtain an optimum, diffuse distribution on the surface of the tooth substance (dentine and enamel) for the opto-electronic impression taking, e.g. with the CEREC intraoral camera.

In order to highlight interferences (e.g. premature contacts, occlusion contacts, approximal contacts) during the try-in of restorations.

#### • Extraoral use

As a contrast powder for model materials, e.g. dies, in order to capture these with a scanner.

### Note

Spray onto clean and dry surfaces only. Failure to heed this advice can lead to failure of the restoration.

### Interactions

There are no known interactions to date. In seldom, individual cases allergic reactions are possible.

### Danger

#### Hazard warnings / Safety advice

- Pressurized container: May burst if heated.
- Avoid breathing spray
- Use only outdoors or in a well-ventilated area
- Do not pierce or burm, even after use.
- Protect from sunlight. Do not expose to temperatures exceeding 50 °C/122 °F.
- Keep out of the reach of children.

### How to use

#### • Preparing for spraying

1. Attach spray head with spraying needle to the spray can. Before each application, pull gently to check that the needle/nozzle is securely seated in the spray head.

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2. The Spraying needle can be rotated by 360° in order to reach an ideal position with regard to the areas to be sprayed.
3. The spraying needle can be bent slightly. Never bend to a right angle!
4. Thoroughly clean and dry areas to be sprayed.
5. Shake the spray can well before spraying! The mixing ball must be clearly audible.

#### • Spraying procedure

Before spraying on the tooth for the first time, spray briefly outside the patient's mouth!

1. One brief squirt is sufficient for each tooth.
2. Make sure the container is held upright when spraying.
3. While spraying, press the spray head briefly to release a small amount of the material each time. Use sparingly.
4. The distance between the spray can and the surface to be sprayed should – if possible – be approx. 10 – 15 mm. The surface must be completely coated with a thin, homogeneous layer.
5. Capture the image with the scanner according to the manufacturer's instructions.

After applying the pigment film in the mouth, let it dry briefly.

For reasons of hygiene and to always obtain optimal spray distribution, it is recommended that the spray application tip be cleaned after each use.

#### • Cleaning the surfaces after successful image capture

1. After image capture, the surface can be cleaned using a water or air jet. In order to clean the surface completely, repeat the procedure if necessary.
2. For intraoral use, rotating prophylaxis brushes can be additionally used for cleaning.

### Storage

Store at room temperature. Protect from direct sunlight and temperatures over 50 °C /122 °F.

### Disposal

Dispose of only after removal of all residue.

### Labelling

The shelf life (hour glass symbol) and the batch number (LOT) are displayed on the container.

### Supply schedule

VITA Powder Scan Spray, container 75 ml, incl. 2 irrigation cannulas.

Prod. No. ECSCAN75V2

### Description of symbols

Expiration date Lot number

Observe instructions for use Date of manufacture Manufacturer

The product was developed for dental use and must be used according to the instructions for use. Information and illustrations contained therein do not constitute a guarantee of properties and are not binding. Date of issue of this information for use: 05.2019.

VITA Powder Scan Spray is certified according to the guideline for medical devices and bears the certification

RX only

CEREC® and inLab® are registered Trademarks of Sirona Dental GmbH, Bensheim, Germany.

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CEREC® e inLab® sono marchi registrati della Sirona Dental GmbH, Bensheim, Germany.
PRISMAN GmbH
VITA Zahnfabrik H. Rauter GmbH & Co. KG