

VITA VIONIC VIGO®

Instructions for use for the digital workflow with exocad® I Full version



VITA shade determination

VITA shade communication

VITA shade reproduction

VITA shade control

Date of issue: 2022-09

VITA – perfect match.

VITA

The material system for digital dentures



Dear Customers,
Thank you for choosing VITA VIONIC VIGO as your digital tooth solution!

VITA VIONIC VIGO is a new generation of digital denture teeth, which are optimized to support digital workflows. The VITA VIONIC open material system optimizes each process step of CAD/CAM denture fabrication. Natural tooth esthetics and time-saving production steps enable functional and patient-specific restorations.

To ensure the effective and safe use of this product, please read this brochure carefully and completely.

For detailed information on the CAD/CAM fabrication of full dentures, please refer to the manual provided by your CAD/CAM system manufacturer.

We hope you enjoy VITA VIONIC SOLUTIONS and wish you great results!

Your VITA Product Management Team

Explanation of symbols:



System/technology info



Process



Please note



Links/Tutorials



Information

➤ 1. The Material System	4
➤ 2. System Components	5
2.1 VITA VIONIC WAX	5
2.2 VITA VIONIC BASE	6
2.3 VITA VIONIC VIGO	7
2.4 VITA VIONIC BOND	8
➤ 3. Overall Workflow	9

➤ 4. Scan-Import Process	10
➤ 5. The CAD Process	11

➤ 6. (CAM) Processing	13
6.1 Production of full-sized wax try-ins (process 1)	13
6.2 Fabrication of the wax try-ins with denture teeth (process 2) ...	14
6.3 Fabrication of the final denture bases	15
6.4 Unpacking of the VITA VIONIC VIGO teeth	16

➤ 7. Bonding	17
7.1 Bonding of the VITA VIONIC VIGO teeth in the base	17
➤ 8. Finalization	18
8.1 Finalization of the final denture bases	18

➤ 9. Moulds, technical data, and information ..	20
9.1 Overview of available tooth moulds	20
9.2 Intended purpose	21
9.3 Patient target group	21
9.4 Intended users	21
9.5 Indications	21
9.6 Contraindications	21
9.7 Storage/disposal	21
9.8 Chemical composition	21
9.9 Physical properties	22

9.10 Product reliability	22
9.11 Information and explanations of symbols.	22

1. The Material System



Information:

- What? VITA VIONIC is a compatible material system for process-reliable CAD/CAM denture manufacturing at the touch of a button.
- With what? The VITA VIONIC material system includes:
 - VITA VIONIC VIGO: Denture tooth for CAD/CAM dentures.
 - VITA VIONIC WAX: Wax discs for the fabrication of full-sized wax try-ins and wax setups.
 - VITA VIONIC BASE: PMMA discs for the fabrication of final denture bases.
 - VITA VIONIC BOND: Bonding solution for fixing denture teeth in the base.

*) For the fabrication of the (wax) try-in and the denture base, milling and 3D printing methods from VITA technology partners can be used.

2. System Components

2.1 VITA VIONIC® WAX



Information:

- What? VITA VIONIC WAX is a millable blank made of a high-quality, millable and dimensionally stable synthetic wax with a high melting temperature.
- What for?
 - VITA VIONIC WAX White is used for the economical fabrication of full-sized try-ins. They are milled completely from wax (base, incl. teeth) and are used to test the midline, occlusion plane and phonetics.
 - VITA VIONIC WAX, pink, is used for the fabrication of denture bases for try-in. For this purpose, denture teeth are fixed in the milled cavities. All parameters can be checked, and if necessary, corrections can be implemented.
- With what? VITA VIONIC WAX is available in the colors pink and white.

Please note:

- Not suitable for the direct fabrication of final prostheses.
- Wet machining is recommended for the milling process.
- Must be processed at room temperature.
- Must be protected from direct sunlight.

2.2 VITA VIONIC® BASE



VITA VIONIC Discs



Milled denture base

Information:

- What? VITA VIONIC BASE is a millable blank made of a high-quality, industrially polymerized acrylic polymer that features shade stability (PMMA) for the CAD/CAM fabrication of denture bases for full dentures.
- What for? For the CAD/CAM fabrication of final denture bases in conjunction with VITA denture teeth (VITA VIONIC VIGO).
- With what? VITA VIONIC BASE is available in three colors (Deep Pink, Light Pink, Orange Pink) and in two different heights (26 and 30 mm).

Please note:

- Must not be used in patients who are allergic to PMMA.
- Suitable only for the manufacture of full dentures with VITA denture teeth (VITA VIONIC VIGO).
- Can be lined and repaired with a commercially available cold polymer:
Manufacturing recommendation: Cold polymer FuturaGen (Schütz Dental GmbH) in the colors orange, transparent pink and opaque pink.



2.3 VITA VIONIC VIGO®



Information:

- What? VITA VIONIC VIGO is a VITA tooth specially designed for digital denture workflows.
- What for? For manufacturing dentures using milled or printed denture bases*.
- With what? VITA VIONIC VIGO is available in the following moulds and shades:
 - VITA VIONIC VIGO (anterior): 8 x maxilla and 4 x mandible anterior tooth moulds.
 - VITA VIONIC VIGO (posterior): 4 x maxilla and 4 x mandible posterior tooth moulds.
 - Available in: VITA classical A1–D4/Bleach shades (0M1, A1, A2, A3, A3.5, B3, D3)



Please note:

- Store and process at room temperature.
- Protect from direct sunlight.

*) For compatible materials and systems please visit: www.vita-zahnfabrik.com/vionic_compatibility

2.4 VITA VIONIC® BOND



VITA VIONIC BOND bonding system

Information:

- What? VITA VIONIC BOND is a self-curing, two-component bonding system (BOND I + II) based on methyl methacrylate (MMA).
- What for? It is used for the final bonding of VITA VIONIC VIGO teeth in the cavities of CAD/CAM fabricated denture bases made of VITA VIONIC BASE*.
- With what? The VITA VIONIC BOND KIT consists of VITA VIONIC BOND I (glass vial), VITA VIONIC BOND II (glass bottle) and an applicator (microbrush).

Please note:

- The packaging unit is a single dose for full dentures, consisting of upper and lower jaw dentures.
- Store in the refrigerator between 5 and 10°C where it is dark and dry; observe the expiration date and protect from sunlight.
- VITA VIONIC BOND contains methyl methacrylate (MMA). MMA is a hazardous substance that is highly flammable and has a sensitizing effect. Avoid contact with skin and inhalation of the fumes.
- You can find detailed instructions in the Safety Data Sheets at www.vita-zahnfabrik.com.

*) For compatible materials and systems please visit: www.vita-zahnfabrik.com/VIONIC_compatibility

3. The Overall Workflow



Information:

1. Scan and model analysis.
2. Selection of VITA VIONIC VIGO denture teeth and CAD-design.
3. Manufacturing of a wax setup or a full-sized try-in with VITA VIONIC WAX or a recommended printing material.
4. Fabrication of a denture base with VITA VIONIC BASE or a recommended printing material.
5. Bonding of the denture teeth with VITA VIONIC BOND and finishing.

Please note:

System requirements and compatible materials for VITA VIONIC SOLUTIONS:
Please visit the following URL: www.vita-zahnfabrik.com/vionic_compatibility.

4. Scan Import Process




1 Create the patient case, incl. master data.



2 After the scan, open the patient case for the scan upload.




3 Define the indication "full dentures" via the material configuration.




4 Select the production method: "Mill base" or "Print base."




5 Start the scan selection via the "CAM" button.




6 Import the scan data for the upper jaw.



7 Import the scan data for the lower jaw.



8 Import the scan data for the relation determination (bite registration).



9 To complete the import, confirm the alignment of the data in the space, if necessary.


Information:

- During the scanning process, prepare the models so that no large shadows are visible in the area of the vestibular fold.
- Spray the bite rim with scan spray (e.g., VITA Scan Spray) to prevent reflections and record all necessary information.
- Please follow the instructions of your scanner manufacturer.

Please note:

- For detailed information on the scan process, please see the instructions for use of the scanner manufacturer.
- When selecting the production method for VITA VIONIC VIGO, only select the options "Mill base" or "Print base."
- When importing the scan data, please always pay attention to the prompt in the dialogue box to specify which object must be loaded in this step.

5. The CAD Process




1 Determine the occlusion plane.




2 Conduct the model analysis according to TiF (UJ and LJ).




3 Design upper jaw boundary dimension.



4 Design lower jaw boundary dimension.




5 Define the insertion direction. Pay attention to undercuts.



6 Block out upper jaw undercuts.



7 Define the insertion direction. Pay attention to undercuts.



8 Block out lower jaw undercuts.



9 Selection of VITA VIONIC VIGO teeth according to the setup concept.



! Note: The different setup concepts are selected via the designation of the tooth set. For example here: VITA VIONIC VIGO 23L Lingualized. A bilateral and/or unilateral crossbite set up is available.


Information:

- Use bite rim landmarks, like canine position or lip line, to set accurate markings for model analysis.
- The setup concept is selected by quadrant. The corresponding upper and lower quadrants always have the same functional setup.
- The VITA VIONIC VIGO tooth library offers different setup concepts, and all angle classes can be depicted.
- Support of bilateral and unilateral crossbite selection.
- Select the suitable anterior and posterior teeth and the desired setup concept, according to the TiF model analysis.
- Keep functional relation of teeth by using group function or symmetric movement function. The single tooth movement function should be limited to anterior esthetic modifications.

Please note:

- For detailed information on the CAD process, please see the instructions for use of the software provider.

5. The CAD Process




10 Modify the setup using Tooth Placement.



11 Modify the gingival morphology, if necessary



12 Adhesive gap: 0.02 mm.
Drill compensation: 1 mm.



13 Nesting of the denture bases (maxilla and mandible).

Information:

- Modify the setup suggested by the CAD software. For this, use the bite rim for alignment.
 - If required, design the gingiva individually according to your preferences in the denture base design step.
 - Set the adhesive gap to 0.02 mm for milling denture bases.
 - Set drill compensation to 1 mm. The smallest recommended tool for VITA VIONIC VIGO is 2 mm. Smaller tools are not required for milling the tooth pockets.
 - Printed (3D printing) denture bases may require different settings.
- For further information, please contact your printing material provider and / or visit the following URL:
www.vita-zahnfabrik.com/vionic_compatibility.

Please note:

- For detailed information on the CAD process, please see the instructions for use of the software provider.

6. (CAM) Processing

6.1 Fabrication of full-sized wax try-ins (Alternative Process 1)



1 Clamp the white wax disc in the holder system.



2 Mill the wax disc for the full-sized wax try-in (maxilla).



3 Mill the wax disc for the full-sized wax try-in (mandible).



4 Separate out the try-in with a hot wax knife.



5 The full-sized try-ins on the master models.



6 Conduct the intraoral try-in.

Information:

- This process is used to manufacture full-sized wax try-ins (note: there are two processes for the fabrication of try-ins; see 2.1 VITA VIONIC WAX).
- For the full-sized try-in, blend in or smooth out excess (protrusions) flush with the denture base.
- Central parameters, such as the midline, course of the occlusion plane and phonetics, must be examined during the try-in.
- If all parameters are in order, you can mill the final denture base and conduct the CAM modification of the denture teeth.
- Any necessary adjustments after the try-in are performed digitally with the CAD software.
- The final prosthesis is then fabricated using the CAM process (based on the modified and finalized setup).


Please note:

- Wet machining is recommended for the milling process
- Based on CAM technology, a final prosthesis can only be fabricated after checking the wax try-in.
- In the case of comprehensive adjustments after the wax try-in, conduct another try-in as a control, if necessary.
- The full-sized try-ins must be disinfected before and after the try-in.

▶ 6.2 Fabrication of the wax try-in with denture teeth (Alternative process 2)



1 Clamp the pink wax disc in the holder system.



2 Mill the wax disc (maxilla and mandible) for the try-in.



3 Separate out the try-in with a hot wax knife.



4 Preparation of VITA VIONIC VIGO (see 6.3).



5 Fix the denture teeth in the cavities with wax.



6 The finished wax base with denture teeth secured with wax.


💡 Information:

- Wet machining is recommended for the milling process.
- This process is used to manufacture denture bases for wax try-ins. The VITA VIONIC VIGO teeth are fixed in the milled cavities (note: there are two processes for the fabrication of try-ins, see 2.1 VITA VIONIC WAX).
- Blend in and smooth out the excess (protrusions) flush with the denture base and fix the denture teeth for the cavities try-in.
- Central parameters, such as the midline, course of the occlusion plane and phonetics, must be examined during the try-in.
- If all parameters are in order, you can mill the final denture base or perform the conventional fabrication.
- Any necessary adjustments after the try-in can be conducted digitally (CAD software) or manually (by hand).
- If significant changes have been made to the anterior setting, the try-in is scanned as a new bite rim, and the setup of the teeth is modified using the CAD software.
- The wax dentures must be disinfected before and after the try-in.

6.3 Fabrication of the final denture bases



1 Clamp the PMMA disc in the holder system.




2 Mill the final maxilla denture base.



3 Mill the final mandible denture base.



4 Separate the denture bases from the disc.



5 Final milled denture bases (here, maxilla).



6 Sandblast the milled cavities of the denture bases.

Information:

- Separate the denture bases from the discs with a mill suitable for PMMA and grind the excess (protrusions) flush with the base.
- Carefully sandblast the milled cavities in the denture bases with Al_2O_3 (50 μm , 2–3 bar) and remove the mill chippings with compressed air (with water separator). Observe the operating instructions for the blasting unit used.

Please note:

- Please use appropriate safety glasses and wear face protection/respiratory protection.

Links/Tutorials:

- Learn more in tutorial videos: www.vita-zahnfabrik.com/tutorial/vionic/all/vigo/denturebase

▶ 6.4 Unpacking the denture teeth



1 Check the labeling to ensure correct tooth selection.



2 Open the blister carefully to avoid teeth falling out.



3 The tooth numbers are given in the product information.



4 Please use forceps to remove teeth from the blister.

Explanation of VITA VIONIC VIGO labeling



💡 Information:

- To avoid any contamination of the VITA VIONIC VIGO teeth, please wear talcum powder free examination gloves.
- Please open the single blister pockets carefully to avoid teeth falling out.
- Please read the instructions for use
- The product tooth number is documented in the instructions for use.

⚠ Please note:

- Please use appropriate safety glasses and wear face protection/respiratory protection.

7. Bonding

7.1 Bonding of the VITA VIONIC VIGO denture teeth in the base



1 Check the fit of the teeth for any slight proximal interferences.



2 Mix VITA VIONIC BOND I and II with each other.



3 Mix the two components for 30 seconds.



4 Apply the adhesive system to the adhesive surfaces of the teeth.



5 Moisten the milled cavities generously with the adhesive system.



6 Position the denture teeth in the cavities.



7 Attach the denture teeth using light pressure.



8 Final bonded denture teeth ready for curing.

Information:

- Position all teeth one after the other in the cavities of the denture base to remove any occurring proximal interference.
- One portion of adhesive supports one full denture (one upper and one lower denture).
- Mix the two components using the applicator for 30 seconds to prevent air bubbles.
- After that, immediately apply and process the bonding systems at room temperature (> 20 °C). The adhesive begins curing after 10 minutes.
- Adhesive with a low viscosity may evaporate if it is handled for too long. To ensure that the adhesive surface is completely moistened, a generous application is recommended.
- To create a secure bond, the denture must rest for at least 30 minutes without load after the teeth have been inserted.
- After 20 minutes in the water bath of a pressure pot (55 °C, 2 bar), curing is complete. Alternatively, curing can be achieved by storing the denture for 12 hours at room temperature without any load.
- Any excess can be removed by careful blast polishing (1 – 2 bar) after bonding the teeth.

Please note:

- The mixing ratio of the two components is ideally matched to one another and may not be changed.
- If the low-viscosity consistency of the bonding system changes and threads begin to appear, for example, VITA VIONIC BOND can no longer be used.
- Please use appropriate safety glasses and wear face protection/respiratory protection.

Links/Tutorials:

- Learn more in tutorial videos: www.vita-zahnfabrik.com/tutorial/vionic/all/vigo/bond

8. Finalization


8.1 Finalization of the final denture bases



1 Use VITA VM LC MODELLING LIQUID for moistening.



2 Moisten the sandblasted areas with VITA VM LC MODELLING LIQUID.



3 VITA VM LC flow materials for interdental spaces.



4 Close the interdental spaces with VITA VM LC flow.



5 Finalize and polish the denture surfaces.



6 Final digitally-fabricated full dentures.



7 Check occlusion in the articulator.

Information:

- Sandblast the interdental spaces with Al_2O_3 (50 μm , 2–3 bar) to ensure a good bonding of VITA VM LC flow to the base material and the denture teeth.
- Clean the sandblasted surfaces with compressed air (water separator).
- Moisten the sandblasted areas with VITA VM LC MODELLING LIQUID for the bond between the denture teeth and VITA VM LC flow materials.
- Then use VITA VM LC flow (five gingiva shades and/or Window) to close the interdental spaces.
- Please conduct the surface processing and polishing according to the analogous full denture.
- Check the occlusion in the articulator. After adjustment of the muscle balance, a reocclusion by the practitioner is recommended after a wearing time of about two days.

! Please note:

- When using the light-curing microparticle composite VITA VM LC flow, please observe the instructions for use of this product.
- The curing parameters may vary depending on the device used. Please observe the manufacturer's recommendations (detailed instructions can be found at www.vita-zahnfabrik.com).
- Please use appropriate safety glasses and wear face protection/respiratory protection.





Links/Tutorials:

- Learn more in tutorial videos: www.vita-zahnfabrik.com/tutorial/vionic/all/vigo/finalization



9. Moulds, technical data, and information

9.1 Overview of available tooth moulds

Upper anteriors VITA VIONIC VIGO®						
 ovoid	O43  10.7 7.9	O45  11.3 8.3				
 triangular	T44  11.7 8.1	T46  11.6 8.2	T48  12.5 8.7	T50  12.9 9.2		
 rectangular	R45  10.7 8.6	R47  12.1 8.9				
Lower anteriors VITA VIONIC VIGO®						
	L33  8.9 4.9	L35  9.5 5.2	L37  10.0 5.5	L39  10.5 5.8		
Posterior VITA VIONIC VIGO®						
21L 9.2  10.0	22L 9.6  10.4	23L 10.3  11.3	24L 10.9  11.9			

 **Note:** The VITA VIONIC VIGO prefabricated tooth is a specially developed product with a design perfected for digital production processes. The teeth are shortened in the cervical area and the contact surfaces to the denture base are preconditioned for bonding. As a result, the ready-to-use teeth can be directly taken from the packaging and bonded into the cavities of the digitally fabricated denture base.

VITA VIONIC VIGO Anterior	UJ / No.	R47	R49	O45	T50	T46	O43	T44	T48
	LJ / No.	L37	L39	L35	L39	L35	L33	L33	L37
	mm	47.0	49.0	44.5	49.9	46.2	43.3	44.0	48.3
	mm	37.2	39.2	35.2	39.2	35.2	33.2	33.2	37.2
VITA VIONIC VIGO Posterior	UJ / No.	22L	23L	22L	24L	22L	21L	21L	22L
	LJ / No.	22L	23L	22L	24L	22L	21L	21L	22L
	mm	31.8	34.3	31.8	35.9	31.8	30.3	30.3	31.8
	mm	32.6	35	32.6	36.9	32.6	31	31	32.6

 9.2. Intended purpose

VITA VIONIC VIGO teeth are acrylic denture teeth used for dental treatments.

 9.3. Patient target group

No restrictions

 9.4. Intended user

Dental professionals only: dentists and dental technicians

 9.5. Indications

The VITA VIONIC VIGO denture tooth is designed for the fabrication of permanent and/or temporary partial and full dentures. For use by or upon instruction of a dentist or dental surgeon. For professional use only – not intended for the general public.

 9.6. Contraindications

None known

 9.7. Storage/disposal

Store in original packaging at room temperature. Store in a dry location. Protect from direct sunlight.

The products labelled with a pictogram for hazardous substances are to be disposed of as hazardous waste. Recyclable waste (such as attachments, paper and plastics) must be disposed of using appropriate recycling systems. If necessary, contaminated product residues should be pretreated in accordance with regional regulations and disposed of separately.

 9.8. Chemical composition

Chemical composition	Wt%
Polymethyl methacrylate (PMMA) CAS no. 9011-14-7	90–98
Silicon dioxide CAS no. 14808-60-7	2–9
Pigments	0–1



9.9. Physical properties

Physical properties	Value
Appearance according to ISO 22112	Complies with standard
Tooth shape according to ISO 22112	Complies with standard
Shade and blending of shades according to ISO 22112	Complies with standard
Free of porosity and other defects according to ISO 22112	Complies with standard
Surface treatment according to ISO 22112	Complies with standard
Adhesive strength to denture base according to ISO 22112	Complies with standard
Resistance to fading, deformation, and cracking according to ISO 22112	Complies with standard
Shade stability according to ISO 22112	Complies with standard



9.10. Product reliability

Information on reporting serious incidents in connection with medical devices, general risks associated with dental treatments, residual risks, and (if applicable) short clinical safety and performance reports (SSCPs) can be found at https://www.vita-zahnfabrik.com/product_safety.

The safety data sheets can be downloaded at www.vita-zahnfabrik.com or requested by fax at (+49) 7761-562-233.



9.11. Information on and explanations of symbols

Medical device		Manufacturer	
For professionals only	Rx only	Date of manufacture	
Observe instructions for use		Expiration date	
Do not use if the packaging is damaged. Follow the instructions for use		Product number	
Store in a dry location		Batch description	
Protect from sunlight			

VITA VIONIC BOND I		
	Danger	H225 Highly flammable liquid and vapor.
	Important	H315 Causes skin irritation. H317 May cause an allergic skin reaction. H335 May cause respiratory irritation.
VITA VIONIC BOND II		
	Danger	H225 Highly flammable liquid and vapor.
	Important	H315 Causes skin irritation. H317 May cause an allergic skin reaction. H335 May cause respiratory irritation.
VITA VM LC MODELLING LIQUID		
	Important	H315 Causes skin irritation. H317 May cause an allergic skin reaction. H319 Causes serious eye irritation. H335 May cause respiratory irritation.
VITA VM LC flow GINGIVA		
	Important	H315 Causes skin irritation. H317 May cause an allergic skin reaction. H319 Causes serious eye irritation. H412 Harmful to aquatic life with long lasting effects.
VITA VM LC flow WINDOW		
	Important	H315 Causes skin irritation. H317 May cause an allergic skin reaction. H319 Causes serious eye irritation. H412 Harmful to aquatic life with long lasting effects.
 	Safety at work and health protection	When working with the product, wear suitable safety goggles/face protection and light respiratory protection.

WE ARE HAPPY TO HELP

More information about the products and processing is also available at www.vita-zahnfabrik.com



Hotline Sales Support

Mrs. Carmen Holsten and her team (Internal Sales Department) will be glad to assist you with orders or questions about the delivery, product data and marketing materials.

► **Phone +49 (0) 7761 / 56 28 84**
Fax +49 (0) 7761 / 56 22 99
8:00 a.m. to 5:00 p.m. CET
E-mail info@vita-zahnfabrik.com

For additional international contact information please visit www.vita-zahnfabrik.com/contacts



Technical Hotline

If you have technical questions concerning VITA product solutions, you can contact Dr. Michael Tholey and his Technical Service team.


► **Phone +49 (0) 7761/56 22 22**
Fax +49 (0) 7761/56 24 46
8:00 a.m. to 5:00 p.m. CET
e-mail info@vita-zahnfabrik.com

For additional product informations please visit www.vita-zahnfabrik.com/vionivigo



VITA PROSTHETIC SOLUTIONS –

For ideal dental prosthetics: natural, reliable, rich in variation.




> Efficient CAD/CAM dentures for greater productivity.

Increasing productivity with digital technology

The VITA DIGITAL PROSTHETICS combines innovation, efficiency and precision with precisely matched material and CAD/CAM system solutions for full dentures.

More information about VITA VIONIC VIGO®
is available at www.vita-zahnfabrik.com/vioniccvgio



Please note: Our products must be used in accordance with the instructions for use. We accept no liability for any damage resulting from incorrect handling or usage. The user is furthermore obliged to check the product before use with regard to its suitability for the intended area of applications. We cannot accept any liability if the product is used in conjunction with materials and equipment from other manufacturers that are not compatible or not authorized for use with our product and this results in damage. The VITA Modulbox is not necessarily a component of the product. Date of issue of this information: 2022-09

After the publication of this information for use any previous versions become obsolete. The current version can be found at www.vita-zahnfabrik.com

VITA Zahnfabrik has been certified and the following products bear the mark

CE 0124

VITA VIONIC VIGO®, VITA VIONIC® BOND, VITAVM®LC flow, VITAVM®LC MODELLING LIQUID

The company GDF Gesellschaft für dentale Forschung und Innovationen GmbH is certified in accordance with the Medical Device Directive, and the following product bears the mark

CE 0297

VITA VIONIC® BASE

VITA Zahnfabrik is certified, and the following products bear the mark

CE

VITA VIONIC® WAX

DentalCAD® and exocad® are registered trademarks of exocad GmbH, Julius-Reiber-Str. 37, D-64293 Darmstadt, Germany. Other products/systems of other manufacturers mentioned in this document are registered trademarks of the respective manufacturers.

Rx only    

CH REP

VITA Zahnfabrik H. Rauter GmbH & Co.KG, Bad Säckingen (Germany)
Zweigniederlassung Basel c/o Perrig AG, Max Kämpf-Platz 1, 4058 Basel

VITA



VITA Zahnfabrik H. Rauter GmbH & Co.KG

Spitalgasse 3 · D-79713 Bad Säckingen · Germany

Tel. +49(0)7761/562-0 · Fax +49(0)7761/562-299

Hotline: Tel. +49(0)7761/562-222 · Fax +49(0)7761/562-446

www.vita-zahnfabrik.com · info@vita-zahnfabrik.com

 facebook.com/vita.zahnfabrik