VITA – perfect match.

VITA

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VITA LUMEX® AC Verblendkeramik

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
Instructions for use

VITA LUMEX® AC veneering ceramic

Please read the instructions for use before using the product.

Intended use

VITA LUMEX AC is a leucite-reinforced, glass-ceramic veneering system for veneering all common, full-ceramic substructure materials (zirconia, lithium disilicate and feldspar ceramics) and titanium substructures. For the production of restorations without a substructure (e.g., veneers).

Patient target group

No restrictions

Indication

Indication range:

- Full and partial veneering of zirconia
- Full and partial veneering of lithium disilicate
- Full and partial veneering of titanium grades 4 and 5
- Partial veneering of feldspar ceramic
- Reconstruction without a substructure

Materials:

- Zirconia substructures (CTE approx. 10.0 to 10.5 x 10⁻⁶ K⁻¹)
- Glass-ceramic substructures (CTE approx. 9.0 to 10.5 x 10⁻⁶ K⁻¹)
- Titanium substructures (CTE approx. 9.0 to 10.5 x 10⁻⁶ K⁻¹)

Contraindication

- Substructures with unsuitable CTE values and material properties
- In patients with allergies or sensitivities to the ingredients
- In cases of insufficient space

Product application/control

Workflow:

Please adhere to the following workflow when working with this product:

1. Substructure preparation
2. Wash application (optional for lithium disilicate)
3. Opaque firing (ZrO₂ (optional) and titanium)
4. Application of DENTINE
5. Application of ENAMEL
6. Application of corrective material (optional)
7. Adjustments of the shape
8. Characterization/glaze firing

Firing:

Predry. °C	→ min.	↑ °C/min.	approx. temp. °C**	↕ °C	↘ °C	Vac.
400	4:00	50	Washbake/ Opaque firing ZrO ₂ + titanium 800	1.0	500*	on
	6:00		First dentine firing 760			on
	08:00	50	Glaze firing - VITA AKZENT PLUS GLAZE LT Paste 750			–

* Long-term cooling down (without holding time) to the corresponding temperature is recommended for the respective last planned firing cycle of the veneering ceramic.

** The firing temperatures indicated represent the highest and lowest recommended firing cycle of the veneering ceramic firing as examples. If necessary, additional firing programs can be found in the detailed instructions for use.

The specified firing temperatures are only reference values. The firing results are influenced by various parameters such as model, performance and age of the firing unit used. Detailed information on processing and firing can be found in the document "Instructions for Use I Full Version" at www.vita-zahnfabrik.com/lumex

Recommended product combinations

VITA AKZENT PLUS stains can be used for characterizing VITA LUMEX AC restorations.

Storage/disposal

The veneering ceramic must be stored in a dry place. The products labelled with a pictogram for hazardous substances are to be disposed of as hazardous waste. Recyclable waste (such as attachments, paper and plastics) must be disposed of using appropriate recycling systems.

Technical data

Physical properties	Unit of measure	Value
3-point flexural strength	MPa	approx. 100
Solubility in acids	µg/cm ²	approx. 10
Coefficient of thermal expansion (CTE)	10 ⁻⁶ K ⁻¹ (25-400°C)	approx. 8.8

Chemical composition	Wt%
SiO ₂	60–75
Al ₂ O ₃	3–10
K ₂ O	5–12
Na ₂ O	4–11
B ₂ O ₃	5–12
CaO	< 3
Li ₂ O	< 3
Pigments	< 10

General risks/reporting incidents

Information on the general risks associated with dental treatments and on reporting serious incidents in connection with medical devices can be found at www.vita-zahnfabrik.com/service_risks
Safety data sheets can be downloaded at www.vita-zahnfabrik.com or requested by fax from (+49) 7761-562-233.

Disclaimer

Please note: Our products must be used in accordance with the instructions for use. We do not accept any liability for damage resulting from incorrect handling or usage. The user is furthermore obliged to check the product before use with regard to its suitability for the intended area of applications. We cannot accept any liability if the product is used in conjunction with materials or equipment from other manufacturers that are not compatible or not authorized for use with our product and this results in damage. The VITA Modulbox is not a required component of the product. Date of issue of this information: 2022-08
After the publication of this information for use, any previous versions become obsolete. The current version can be found at www.vita-zahnfabrik.com

VITA Zahnfabrik has been certified and the following products bear the mark: