VITA – perfect match.

VITA

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VITA VM[®] 11 Verblendkeramik

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Instructions for use

VITA VM[®] 11 veneering material

Please read the instructions for use before using the product.

Intended use

VITA VM 11 is a low fusing fine-structure feldspar ceramic that has been developed especially for individualizing restorations made of zirconia reinforced lithium silicate ceramic (ZLS).

Patient target group

No restrictions

Indication

- For individualizing restorations made of zirconia reinforced lithium silicate glass ceramic (ZLS), especially for VITA SUPRINITY PC

Materials:

- VITA SUPRINITY PC zirconia reinforced lithium silicate glass ceramic, CTE (25-500°C), approx. 11.9 – 12.3 x 10⁻⁶ x K⁻¹

Contraindication

- For substructures not within the recommended CTE range
- For parafunctions (e.g., bruxism)
- If minimum layer thicknesses of the ceramic can not be adhered to
- In cases of inadequate oral hygiene

Product application/control

Workflow:

Please adhere to the following workflow when using this product:

- Substructure preparation (e.g., processing with rotary instruments, crystallization),
- Application of DENTINE,
- Application of ENAMEL,
- First dentine firing
- Adjustments of the shape,
- Characterization/glaze firing

Firing:

Predry. °C	→ min.	↗ °C/min.	approx. temp. °C*	→ min.	Vac.
400	6:00	55	1. dentine firing 800	1:00	on
	4:00	80	Glaze firing with AKZENT PLUS 800		–

* As examples, the firing temperatures given represent the highest and lowest recommended firing cycle of the veneering ceramic firing. If necessary, additional firing programs can be found in the detailed instructions for use.

The specified firing temperatures are only reference values. The firing results are influenced by various parameters such as model, performance and age of the firing unit used.

Detailed information on processing and firing can be found in the document "Working Instructions VITA VM11" at www.vita-zahnfabrik.com

Recommended product combinations

VITA AKZENT PLUS stains can be used for characterizing VITA VM11 restorations.

Storage/disposal

The veneering ceramic must be stored in a dry place. The products labelled with a pictogram for hazardous substances are to be disposed of as hazardous waste.

Recyclable waste (such as attachments, paper and plastics) must be disposed of using appropriate recycling systems.

Technical data

Physical properties	Unit of measure	Value
3-point flexural strength	MPa	approx. 100
Solubility in acids	µg/cm ²	approx. 8
Coefficient of thermal (25-500°C) expansion (CTE)	10 ⁻⁶ K ⁻¹	11.2 – 11.6

Chemical composition	Wt%
SiO ₂	62 – 65
Al ₂ O ₃	8.5 – 12
K ₂ O	9 – 12
Na ₂ O	5 – 7.5
B ₂ O ₃	4 – 6
CaO	1 – 2
ZrO ₂	< 1

General risks/reporting incidents

Information on the general risks associated with dental treatments and on reporting serious incidents in connection with medical devices can be found at www.vita-zahnfabrik.com/service_risks Safety data sheets can be downloaded at www.vita-zahnfabrik.com or requested by fax at (+49) 7761-562-233.

Disclaimer

Please note: Our products must be used in accordance with the instructions for use. We do not accept any liability for damage resulting from incorrect handling or usage. The user is furthermore obliged to check the product before use with regard to its suitability for the intended area of applications. We cannot accept any liability if the product is used in conjunction with materials or equipment from other manufacturers that are not compatible or not authorized for use with our product and this results in damage. The VITA Modulbox is not a required component of the product.

Date of issue of this information: 2021-07

After the publication of this information for use, any previous versions become obsolete. The current version can be found at www.vita-zahnfabrik.com

VITA Zahnfabrik has been certified, and the following products bear the mark:

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