# VITA AKZENT® LC GLAZE

Gebrauchsanweisung | Kurzversion  
Instructions for use | short version  
Mode d'emploi | Version courte  
Instrucciones de uso | Versión abreviada  
Istruzioni d'uso | -versione breve

## VITA – perfect match. VITA

# VITA

VITA Zahnfabrik · H. Rauter GmbH & Co. KG  
Spitalgasse 3 · 79713 Bad Säckingen  
Germany · Allemagne · Alemania · Germania  
Tel. +49 (0) 7761 / 562-0  
Fax +49 (0) 7761 / 562-299  
www.vita-zahnfabrik.com  
info@vita-zahnfabrik.com

922-12770 – 0621 (X.)TP –  
Version (01)

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### Lichthärtender Glasurlack

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## Instructions for use (EN)

### Light-curing glaze varnish

Please read the instructions for use before using the product.

### Intended use

VITA AKZENT LC GLAZE is a transparent, light-curing methacrylate-based varnish for extraoral surface sealing of dental restorations made of hybrid ceramic, resin veneering materials, CAD/CAM composites, prefabricated teeth and denture-base resins.

### Patient target group

No restrictions

### Indication


Indication range:

- Restorations made of hybrid ceramic (VITA ENAMIC)
- Restorations made of light-curing veneering material (e.g., VITA VM LC)
- Restorations made of CAD/CAM composites (e.g., VITA CAD-Temp)
- Prefabricated teeth (e.g., VITAPAN)
- Denture bases (e.g., VITA VIONIC BASE)

### Contraindication

- In patients with allergies or sensitivities to the ingredients
- Do not use on occlusal contact points of restorations
- Not approved for intraoral use

### Product application/control Pretreatment steps according to material

\* These are only general recommendations for pretreatment. Depending on the type of material, the specific manufacturer's instructions must be observed.

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## Applying VITA AKZENT LC GLAZE

- Use microbrush disposable applicators or brushes for the application of glaze.
- Apply GLAZE swiftly and without streaks.
- Close bottle again after use.

Please note the recommended processing time of approx. three minutes and clean the applicator brush with VITA AKZENT LC CLEANER after use.

### Polymerization

The varnish can be polymerized with many standard dental light-curing units with a spectral range of 350 nm – 500 nm.

### Recommended polymerization devices/lamps

Dental technical polymerization devices			
Company	Polymerization device	Polymerization time*	Note
Shofu	Solidite V	3 min	The object must be placed in the cone of light in the center of the chamber. The object must not be placed on the floor!
Shofu	Solidite EX	4 min	see Solidite V
DequDent/Dentoply	Eclipse® junior VLC Curing Unit	Basic 3 – 3 min	The object must be placed in the center of the chamber.
DequDent/Dentoply	Triad® 2000	6 min	rotating on the rotary table
Brodent	Iris Lux Power Unit 2	1 x program FI (90 s)	-
3M ESPE	Vivoio Beta vario	2 x program; 7 min incl. 10 x vacuum	-
GC	LabLight Duo	3 min	The object must be placed in the cone of light high in the center of the chamber. The object must not be placed on the floor!
Hager & Werken	Speed LabLight®	3 min	The object must be placed in the cone of light high in the center of the chamber. The object must not be placed on the floor!
Kulzer	HLite® power	90 s	The object must be placed high in the chamber. The object must not be placed on the floor!
Kulzer	Herflash	90 s	see HLite Power
Kulzer	UnixSP®	90 s	see HLite Power
Kulzer	Dentastar® XS	90 s	see HLite Power
Ivoclar Vivadent	Lumamat® 100	Heating level 0 V0 + 0 V0 + 10:00 min (BP = 10:00 min.)	The times and heating levels must be programmed by the user! V0 = Preheating V0 + Heating BP = Resulting curing time

### Note:

- Dental polymerization devices are always preferable to dental hand lamps, as they ensure much more even exposure of the object.

Dental polymerization lamps			
Company	Polymerization device	Polymerization time*	Note
Ivoclar Vivadent	BioPhase G2	4 x 20 s	"High" mode
Ultradent	VALO LED	4 x 20 s	"Standard" mode The light cones of the two wavelengths must completely overlap -> very small distance to the restoration

\* The specified polymerization times are minimum times.

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Urethan(meth-)acrylate 40 – 60

Siliziumdioxid 8 – 11

Ethyl-phenyl(2,4,6-trimethylbenzoyl)phosphinat 2 – 6

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## Please note:

- The necessary light intensity and wavelength of the polymerization device must be observed to ensure complete curing.
- To avoid irritation of the mucous membrane, appropriate polymerization of VITA AKZENT LC GLAZE must be ensured.
- You can find the latest information on the recommended polymerization devices at: [www.vita-zahnfabrik.com/akzentlc](http://www.vita-zahnfabrik.com/akzentlc)

### Important information:

- **All surfaces must be completely polymerized and hard. There must be no sticky surfaces.**
- **After final polymerization, polish the glazed surface with a suitable polishing paste (such as Renfert Polish hybrid materials, for example) and a soft goat hair brush. Then use a dry cotton buff to achieve the final gloss.**

### Recommended product combinations

VITA AKZENT LC GLAZE can be used to seal restorations that have been characterized using VITA ENAMIC STAINS.

### Storage/disposal

Store at 4 °C - 25 °C (39 °F - 77 °F). Storage in a refrigerator is recommended.

Do not use after the expiry date!

The products labelled with a pictogram for hazardous substances are to be disposed of as hazardous waste. All other products are to be disposed of as mixed industrial waste.

### Technical data

#### Chemical composition (% by weight)

Methyl methacrylate and multifunctional

methacrylates 30 – 40

Urethane(meth-)acrylates 40 – 60

Silicon dioxide 8 – 11

Ethyl phenyl(2,4,6-trimethylbenzoyl)phosphinate 2 – 6

Other < 1

Pigments < 2

#### General risks/reporting incidents

Information on the general risks associated with dental treatments and on reporting serious incidents in connection with medical devices can be found at [www.vita-zahnfabrik.com/service\\_risks](http://www.vita-zahnfabrik.com/service_risks)

Safety data sheets can be downloaded at

[www.vita-zahnfabrik.com/vita-akzent-lc](http://www.vita-zahnfabrik.com/vita-akzent-lc).

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## Hazardous substances

Danger

Contains methyl methacrylate, two propenoic acid reaction products

with pentaerythritol, various acrylates

H317 May cause an allergic skin reaction

H 318 Causes serious eye damage

H 412 Harmful to aquatic life with long-lasting effects

### Disclaimer

Please note: Our products must be used in accordance with the instructions for use. We do not accept any liability for damage resulting from incorrect handling or usage. The user is furthermore obliged to check the product before use with regard to its suitability for the intended area of applications. We cannot accept any liability if the product is used in conjunction with materials or equipment from other manufacturers that are not compatible or not authorized for use with our product and this results in damage. The VITA Modulbox is not necessarily a component of the product. Publication of these working instructions: 2021-07

After the publication of these working instructions, any previous versions become obsolete. The current version in each case can be found at [www.vita-zahnfabrik.com](http://www.vita-zahnfabrik.com)

VITA Zahnfabrik has been certified and the following products bear the marks:

VITA AKZENT® LC

CE 0124 Rx only MD

CERASMART® is a registered trademark of GC Dental Products Corp, Aichi, Japan

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