

VITA AMBRIA® pressing parameters

VITA VACUMAT 6000 MP (VITA Zahnfabrik)

Pressing programs VITA AMBRIA										
Programs	Pre-drying temperature	Rise time	Heating rate / temperature gradient	Pressing temperature	Holding time	Pressing time 1	Pressing pressure 1	Pressing time 2	Pressing pressure 2	Cooling temperature
	°C	min	°C/min	°C	min	min	bar	min	bar	°C
100 g investment ring	700	3:36	50	880	25:00	3:00	3.0	1:00	3.0	–
200 g investment ring	700	3:48	50	890	30:00	3:00	3.0	1:00	3.0	–

We recommend calibrating the furnace prior to the first pressing (see VITA AMBRIA calibration information #10642).

Firing programs VITA AKZENT Plus / VITA LUMEX AC								
Programs	Pre-drying temperature	Pre-drying time	Temperature gradient	Temperature gradient	Firing temperature	Firing time	Vacuum holding time	Long-term cooling
	°C	min	°C/min	min	°C	min	min	°C
Power Glaze*	400	4:00	55	7:16	800	2:00	–	500**
Stains fixation	400	4:00	80	3:45	700	1:00	–	–
1. Dentin	400	6:00	50	7:12	760	1:00	7:12	500**
2. Dentin	400	6:00	50	7:06	755	1:00	7:06	500**
glazing	400	00:00	80	4:23	750	1:00	–	500**
Glaze material GLAZE LT, FLUOGLAZE LT Spray/Powder	400	4:00	50	7:00	750	1:00	–	500**
Glaze material GLAZE LT Paste	400	8:00	50	7:00	750	1:00	–	500**
Correction	400	4:00	50	6:30	725	1:00	6:30	500**

*) Power Glaze is a firing program which is carried out before the respective first firing cycle. Power Glaze increases the flexural strength of the VITA AMBRIA restoration to > 550 MPa.

***) Long-term cooling to 500 °C is recommended for the respective last firing cycle. The lift position for VITA VACUMAT 6000 M furnaces should be > 75%. The fired items must be protected from drafts after opening the furnace.

The values specified here are guide values and are provided solely as a point of reference. Variations in the firing result can occur. The firing result depends on the performance of the furnace in question, and on the manufacturer and age of the furnace. Due to this, the guide values must be individually adapted to the corresponding furnace.

This information has been carefully prepared and verified by VITA, but does not constitute a guarantee.

VITA AMBRIA® pressing parameters

Programat EP5010 (Ivoclar Vivadent)

Pressing programs VITA AMBRIA					
Programs	Standby	Heating rate	End temperature	Holding time	Stop speed
	B (°C)	t↗ (°C/min)	T (°C)	H (min)	E (µm/min)
100 g investment ring	700	50	880	25:00	300
200 g investment ring	700	50	890	30:00	300

The pressing parameters specified are guide values and apply to Ivoclar Vivadent Programat EP5010.

With older generation furnaces, the temperature in the firing chamber may vary by approx. ± 10 °C, depending on the age of the firing muffle. Temperature adjustments can not be excluded.

We recommend calibrating the furnace prior to the first pressing (see VITA AMBRIA calibration information #10642).

Firing programs VITA AKZENT Plus / VITA LUMEX AC									
Programs	Standby temperature	Closing time	Heating rate	Firing temperature	Holding time	Vacuum 1	Vacuum 2	Long-term cooling	Cooling rate
	B (°C)	S (min)	t↗ (°C/min)	T ₁ (°C)	H ₁ (min)	V1 (°C)	V2 (°C)	L (°C)	tl (°C/min)
Power Glaze*	403	4:00	55	800	2:00	–	–	500	0
Setting stains	403	4:00	80	700	1:00	–	–	–	–
1. Dentin	403	6:00	50	760	1:00	400	759	500**	0
2. Dentin	403	6:00	50	755	1:00	400	754	500**	0
glazing	403	00:00	80	750	1:00	–	–	500**	0
Glaze material GLAZE LT, FLUOGLAZE LT Spray/Powder	403	4:00	50	750	1:00	–	–	500**	0
Glaze material GLAZE LT Paste	403	8:00	80	750	1:00	–	–	500**	0
Correction	403	4:00	50	725	1:00	400	724	500**	0

*) Power Glaze is a firing program which is carried out before the respective first firing cycle. Power Glaze increases the flexural strength of the VITA AMBRIA restoration to > 550 MPa.

***) Long-term cooling to 500 °C is recommended for the respective last firing cycle. The lift position for VITA VACUMAT 6000 M furnaces should be > 75%. The fired items must be protected from drafts after opening the furnace.

The values specified here are guide values and are provided solely as a point of reference. Variations in the firing result can occur. The firing result depends on the performance of the furnace in question, and on the manufacturer and age of the furnace. Due to this, the guide values must be individually adapted to the corresponding furnace.

This information has been carefully prepared and verified by VITA, but does not constitute a guarantee.

VITA AMBRIA® pressing parameters

Vario Press 300 (Zubler GmbH)

Pressing programs VITA AMBRIA; Program type: i-Press									
Programs	Start temperature	Heating rate	Holding temperature	Holding time	Extra pressing time	Maximum pressing time	Pressure	Vacuum level	Opening time
	°C	°C/min	°C	min	min	min	low	mm	min
100 g investment ring	700	50	880	20:00	00:00	5:00	Low	720	00:00
200 g investment ring	700	50	890	25:00	00:00	5:00	Low	720	00:00

We recommend calibrating the furnace prior to the first pressing (see VITA AMBRIA calibration information #10642).

Firing programs VITA AKZENT Plus / VITA LUMEX AC																
Programs	Start temperature	Predrying	Predrying time	Closing time	Homogenization temperature	Homogenization time	Heating rate	Final temperature	Holding time	Predrying	Predrying time	Homogenization temperature	Homogenization time	Heating rate	Final temperature	Holding time
	°C	-	min	min	°C	min	°C/min	°C	min	/	/	°C	min	-	-	°C
Power Glaze*	400	No	/	4	400	0.00	55	800	2	/	/	800	/	No	/	/
Setting stains	400	Yes	2	2	400	0.30	80	700	1	/	/	700	/	No	/	/
1. Dentin	400	Yes	4	2	400	0.30	50	760	1	/	/	760	/	Yes	Heat.	760
2. Dentin	400	Yes	3	3	400	0.30	50	755	1	/	/	755	/	Yes	Heat.	755
glazing	400	No	/	4	400	0.00	80	750	1	/	/	750	2	No	/	/
Glaze material GLAZE LT, FLUOGLAZE LT Spray/Powder	400	Yes	4	2	400	0.00	50	750	1	/	/	750	2	No	/	/
Glazing with GLAZE LT Paste	400	Yes	8	2	400	0.00	50	750	1	/	/	750	2	No	/	/
Correction	400	Yes	2	2	400	0.30	50	725	1	/	/	725	2	Yes	Heat.	725

*) Power Glaze is a firing program which is carried out before the respective first firing cycle. Power Glaze increases the flexural strength of the VITA AMBRIA restoration to > 550 MPa.

**) Long-term cooling to 500 °C is recommended for the respective last firing cycle.

The values specified here are guide values and are provided solely as a point of reference. Variations in the firing result can occur. Due to this, the guide values must be individually adapted to the corresponding furnace. This information has been carefully prepared and verified. However, it is provided without guarantee. Special advanced pressing programs for various furnace types are provided free of charge by Zubler GmbH, Ulm.

VITA – perfect match. **VITA**