CAD/CAM Disc PINK

Gebrauchsinformation

Instructions for use

Mode d'emploi

Istruzioni d'uso

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EN

Instructions CAD/CAM Disc PINK

For making denture bases. Acrylic discs made out of PMMA (Polymethylmethacrylate) for use in dental milling systems.

Indications

PMMA Blanks are made of a PMMA-acrylic milling material. These Blanks are for the making of denture bases. This material can be individualized with any commercially available denture staining product.

Development / Polishing

The development is carried out with cross-cut, carbide burs.

To polish, use a commercially available, high-gloss polishing paste and goat hair brushes.

Important notice: Never let the acrylic overheat when using burs mills or using steam jets etc.

Process

Optimally, while milling, the temperature conditions should lie between 10–25°C.

Please always check the precautions found in the milling machine instructions.

The milling parameter strongly depends on the type of machine, the CAD CAM software and the drill-geometry used. If you are unsure, please go by these factors when choosing the parameter for your PMMA-acrylics.

The holes are drilled according to the desired geometrical properties of the acrylic teeth. Apply a thin coat of chemical roughening agent to the teeth which have already been roughened and teeth which make up the undercut. Ensure that there are no remaining pools of liquid. Allow the prepared teeth to dry for 2–3 minutes. Then position the teeth on the palatal plate fix with a thin layer of PMMA-denture acrylic. Use a pressure chamber to cure the material (2–4 bar, 45°C, 15 min).

Safety notice

Dust is usually generated when milling the discs. This could damage the lungs or cause the eyes and skin to become irritated. Please ensure that you always work with a properly functioning suction.

As a result of this, please follow the following safety precautions:

- Always work under a nearby suction
- Avoid letting the product make contact with the eyes or mucous membranes.
- Do not inhale the dust.
- Protective masks and safety goggles must be worn at all times.

Storage

The material should be kept dry and stored at 10–30°C. Keep out of direct sunlight and damp conditions.

Disposal

Dispose of any leftover material along with household waste.

Warranty

Any technical advice concerning our recommendations about the application method, whether written, spoken or given during practical instruction, has been based on our personal experience. Therefore, all advice should be thought of as guidelines only.

Durability

The maximum shelf life is printed on the label. Do not use after the expiration date.

Side effects

With proper use of this medical device, unwanted side effects are extremely rare. Reactions of the immune system (allergies) or local discomfort, however, cannot be ruled out completely. Should you learn about unwanted side effects – even if it is doubtful that the side effect has been caused by our product – please kindly contact us.

Contra-indications/interactions

If a patient has known allergies against or hypersensitivities towards a component of this product, we recommend not to use it or to do so only under strict medical supervision. In such cases, we will supply the composition of our medical device upon request. The dentist should consider known interactions and cross-reactions of the product with other materials already in the patient's mouth before using the product.

Note

Please supply the dentist with the above informations, if this medical device is used to produce a special model. Please also consider safety data sheets.

Physical Properties

Composition: > 98% PMMA,
< 2 % inorganic pigments

Vickers hardness HV50: 14,2
Flexural strength: 91 MPa
E- Module: 2650 MPa
Water absorption: 26,8 µg/mm³
Water solubility: 1,6 µg/mm³
Residual monomer: < 1,0 %
Expires after: 10 years

The technical data are median values and may vary slightly for the delivery batches.

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