VITA SUPRINITY PC®

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VITA SUPRINITY PC®

Please read the instructions for use before using the product.

VITA SUPRINITY PC is a zirconia reinforced lithium silicate glass ceramic for dental CAD/CAM applications for the fabrication of inlays, onlays, partial crowns, veneers, anterior and posterior crowns and anterior and posterior single tooth restorations on implant abutments.

Intended purpose

VITA SUPRINITY PC products are ceramic materials for dental treatments.

Patient target group

No restrictions.

Intended user

Dental professionals only: dentists and dental technicians

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Contraindication

- General
 - inadequate preparation results
 - insufficient remaining natural tooth substance
 - Insufficient space available
- Parafunction
Restorations made of VITA SUPRINITY PC are contraindicated for patients diagnosed with excessive masticatory functions, in particular, teeth grinders and clenchers. Restoring devitalized teeth of patients with hyperfunctions is absolutely contraindicated.
- Bridges
- Veneering: Full veneers on molar crowns using veneering ceramic

Product application/control

Please adhere to the following workflow when using this product:

- 1) Preparation/impression taking
- 2) CAD/CAM processing of the block
- 3) Finishing
- 4) Manual polishing / staining technique / veneering technique
- 5) Bonding

Recommended parameters for crystallization of VITA SUPRINITY PC restorations.

Predry. °C	→ ↕↔↕	↗ min.	↗ °C/min.	T °C	→ ↕↔↕	VAC min.	↘ °C
400	4.00	7.49	55	830	8.00	8.00	600

* The firing chamber must not be opened during long-term cooling.

The specified firing temperatures are only reference values. The firing results are influenced by various parameters such as model, performance and age of the firing unit used. Detailed information on processing and firing can be found in the document "Instructions for Use | Full Version" at www.vita-zahnfabrik.com/suprinity.

Bonding recommendation

VITA SUPRINITY PC	Adhesive bonding	Self-adhesive bonding	Conventional bonding
Veneers	●	X	X
Inlays, onlays, partial crowns	●	X	X
Anterior crowns	●	●	○
Posterior crowns	●	●	○

● recommended ○ possible X not possible

Recommended product combinations

For recommended product combinations with the VITA system components and corresponding processing recommendations and restrictions, please refer to Instructions for Use | Full Version.

Storage/disposal

Store the VITA SUPRINITY PC blocks in the original packaging and in a dry place. They can be disposed of with household waste. The products labelled with a pictogram for hazardous substances are to be disposed of as hazardous waste. Recyclable waste (such as attachments, paper and plastics) must be disposed of using appropriate recycling systems. The materials must not be contaminated with foreign substances (e.g., during the grinding process). If necessary, contaminated product residues should be pretreated in accordance with regional regulations and disposed of separately.

Technical data

Physical properties	Unit of measure	Value**
CTE	10 ⁻⁶ · K ⁻¹	11.9 – 12.3
3-point flexural strength	MPa	approx. 420
Modulus of elasticity	GPa	approx. 70
Hardness according to Vickers (HV)	MPa	approx. 7000
Chemical solubility	µg/cm ²	approx. 40

Chemical composition	Wt%**
ZrO ₂ (zirconia)	8-12
SiO ₂ (silicon dioxide)	56-64
Li ₂ O (lithium oxide)	15-21
La ₂ O ₃ (lanthanum oxide)	0.1
Pigments	< 10
Various	> 10

* Information according to ISO 6872 / ** Source: Internal study, VITA

Safety information

When working with the product, wear suitable safety goggles/face protection and light respiratory protection.

Product reliability

Information on reporting serious incidents in connection with medical devices, general risks associated with dental treatments, residual risks and (if applicable) short clinical safety and performance reports (SSCPs) can be found at www.vita-zahnfabrik.com/product_safety. Safety data sheets can be downloaded at www.vita-zahnfabrik.com or requested by fax at (+49) 7761-562-233.

Disclaimer

Our products must be used according to the instructions for use. We accept no liability for any damage resulting from incorrect handling or usage. The user is furthermore obliged to check the product before use with regard to its suitability for the intended area of applications. We cannot accept any liability if the product is used in conjunction with materials and equipment from other manufacturers that are not compatible or not authorized for use with our product and this results in damage. The VITA Modulbox is not necessarily a component of the product. Date of issue of this information for use: 2022-03.

After the publication of this information for use, any previous versions become obsolete. The current version can be found at www.vita-zahnfabrik.com.

VITA Zahnfabrik has been certified and the following products bear the mark:

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CH REP VITA Zahnfabrik H. Rauter GmbH & Co.KG, Bad Säckingen (Germany)
Zweigniederlassung Basel c/o Perrig AG, Max Kämpf-Platz 1, 4058 Basel
VITA Zahnfabrik · H. Rauter GmbH & Co. KG