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### VITA VIONIC VIGO®

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### VITA VIONIC VIGO®

#### denture teeth

Please read the instructions for use before using the product.

VITA VIONIC VIGO is a VITA tooth especially developed for the digital denture workflow. VITA VIONIC VIGO is suitable for manufacturing dentures using a CAD/CAM-produced, milled or printed denture base.

#### Intended purpose

VITA VIONIC VIGO teeth are resin artificial teeth for dental treatments.

#### Patient target group

No restrictions

#### Intended user

Dental professionals only: dentists and dental technicians

#### Indication

The VITA VIONIC VIGO denture tooth is intended for the fabrication of permanent and/or temporary partial and full dentures. For use by or on the order of a dental professional, such as DDS or DMD. Not for use by the general public or OTG.

#### Contraindication

None known

#### Product application/control

Please adhere to the following workflow when using VITA VIONIC VIGO denture teeth:

1. Check label for correct tooth selection
2. Carefully open the blister cap (for the product-related tooth number, see the product information in the blister packaging)
3. Use tweezers to remove teeth from blister
4. Optional wax try-in with VITA VIONIC VIGO teeth:
  - 4a. Wax into cavities for milled wax try-in teeth
  - 4b. Clean the teeth after try-in and remove any residual waxAlternatively, the try-in can be printed as a monoblock or milled from wax without using the physical VITA VIONIC VIGO teeth.
5. For the fit check, insert teeth into CAD/CAM-fabricated denture base
6. Use VITA VIONIC BOND to bond teeth into the denture base

Detailed information on working with this product can be found in the document Instructions for Use I full version at [www.vita-zahnfabrik.com/vionicvigo](http://www.vita-zahnfabrik.com/vionicvigo).

#### Recommended product combinations

Recommended product combinations with the VITA system components, as well as corresponding processing recommendations and restrictions regarding the combinability of the product, can be found in the instructions for use I full version. For an updated overview of compatible materials and systems, see [www.vita-zahnfabrik.com/vionic\\_compatibility](http://www.vita-zahnfabrik.com/vionic_compatibility).

#### Storage/disposal

Store in original packaging at room temperature. Store in a dry location. Protect from direct sunlight.

The products labelled with a pictogram for hazardous substances are to be disposed of as hazardous waste. Recyclable waste (such as attachments, paper and plastics) must be disposed of using appropriate recycling systems. If necessary, contaminated product residues should be pretreated in accordance with regional regulations and disposed of separately.

## Technical data

Chemical composition	Wt%
Polymethyl methacrylate (PMMA) CAS no. 9011-14-7	90-98
Silicon dioxide CAS no. 14808-60-7	2-9
Pigments	0-1

Physical properties	Value
Appearance according to ISO 22112	Complies with standard
Tooth shape according to ISO 22112	Complies with standard
Shade and blending of shades according to ISO 22112	Complies with standard
Free of porosity and other defects according to ISO 22112	Complies with standard
Surface treatment according to ISO 22112	Complies with standard
Adhesive strength to denture base according to ISO 22112	Complies with standard
Resistance to fading, deformation, and cracking according to ISO 22112	Complies with standard
Shade stability according to ISO 22112	Complies with standard

#### Safety note

When working with the product, wear suitable safety goggles/face protection and light respiratory protection.

#### Product reliability

Information on reporting serious incidents in connection with medical devices, general risks associated with dental treatments, residual risks, and (if applicable) short clinical safety and performance reports (SSCPs) can be found at [https://www.vita-zahnfabrik.com/product\\_safety](https://www.vita-zahnfabrik.com/product_safety).

The safety data sheets can be downloaded at [www.vita-zahnfabrik.com](http://www.vita-zahnfabrik.com) or requested by fax at (+49) 7761-562-233.

#### Disclaimer

Please note: Our products must be used in accordance with the instructions for use. We do not accept any liability for damage resulting from incorrect handling or usage. The user is furthermore obliged to check the product before use with regard to its suitability for the intended area of applications. We cannot accept any liability if the product is used in conjunction with materials or equipment from other manufacturers that are not compatible or not authorized for use with our product and this results in damage. The VITA Modulbox is not necessarily a component of the product. Date of issue of this instruction for use: 2022-02.

After the publication of this information for use, any previous versions become obsolete. The current version can be found at [www.vita-zahnfabrik.com](http://www.vita-zahnfabrik.com).

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