VITA PHYSIODENS®

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VITA PHYSIODENS®

Denture teeth

Please read the instructions for use before using the product.

Anterior and posterior denture teeth for the conventional prosthetics workflow. Suitable for all common occlusal concepts. Especially suited for: BLP (BIO-Logical Prosthetic) according to Dr. E. End.

Intended purpose

VITA PHYSIODENS teeth are resin artificial teeth for dental treatments.

Patient target group

No restrictions

Intended user

Dental professionals only: dentists and dental technicians (Rx only)

Indication

- Complete dentures
- Implant constructions
- Telescopic and conus constructions
- Attachments
- Model casting prostheses

Contraindication

None known

Product application/control

For the conventional production of a denture, please follow the common prosthetic concepts and principles, considering the requirements and needs of the individual patient case.

Recommended product combinations

VITAFOL® H: Protects denture and ceramic teeth during denture fabrication

VITACOLL®: Bonding agent for denture teeth

VITAVM®LC: Light-curing, microparticle composite for extraoral use with fixed and removable restorations

VITA AKZENT® LC: Reliable reproduction of individual shade effects

Storage/disposal

The denture teeth can be disposed of in domestic waste. The products labelled with a pictogram for hazardous substances are to be disposed of as hazardous waste. Recyclable waste (such as attachments, paper and plastics) must be disposed of using appropriate recycling systems. If necessary, contaminated product residues should be pretreated in accordance with regional regulations and disposed of separately.

Technical Data

Chemical composition	Wt%
Polymethyl methacrylate CAS no. 9011-14-7	84 - 86
Silicon dioxide CAS no. 14808-60-7	14 - 15
pigments	< 1

Physical properties	value
Appearance according to ISO 22112	Complies with standard
Dimension of teeth according to ISO 22112	Complies with standard
Shade and blending of shades according to ISO 22112	Complies with standard
Free of porosity and other defects according to ISO 22112	Complies with standard
Surface treatment according to ISO 22112	Complies with standard
Bonding to denture base polymers according to ISO 22112	Complies with standard
Resistance to blanching, distortion, and crazing according to ISO 22112	Complies with standard
Shade stability according to ISO 22112	Complies with standard
Dimensional stability according to ISO 22112	Complies with standard

Product reliability

Information on reporting serious incidents in connection with medical devices, general risks associated with dental treatments, residual risks, and (if applicable) short clinical safety and performance reports (SSCPs) can be found at https://www.vita-zahnfabrik.com/product_safety.

The safety data sheets can be downloaded at www.vita-zahnfabrik.com or requested by fax at (+49) 7761-562-233.

Disclaimer

Please note: Our products must be used in accordance with the instructions for use. We do not accept any liability for damage resulting from incorrect handling or usage. The user is furthermore obliged to check the product before use with regard to its suitability for the intended area of applications. We cannot accept any liability if the product is used in conjunction with materials or equipment from other manufacturers that are not compatible or not authorized for use with our product and this results in damage. The VITA Modulbox is not necessarily a component of the product. Date of issue of this information for use: 2022-09.

After the publication of this information for use, any previous versions become obsolete. The current version can be found at www.vita-zahnfabrik.com

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CH REP VITA Zahnfabrik H. Rauter GmbH & Co.KG, Bad Säckingen (Germany)
Zweigniederlassung Basel c/o Perrig AG, Max Kämpf-Platz 1, 4058 Basel
VITA Zahnfabrik · H. Rauter GmbH & Co. KG