VITA CAD-Temp® IS

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VITA CAD-Temp® IS

VITA CAD-Temp IS (Implant Solutions) is a composite material made of acrylate polymer for the fabrication of long-term temporary restorations on implants.

Intended purpose

VITA CAD-Temp products are polymer-based temporary crown and bridge materials for dental treatments.

Patient target group

No restrictions.

Intended user

Dental professionals only: dentists and dental technicians, Rx only

Indication

VITA CAD-Temp IMPLANT SOLUTIONS is indicated for CAD/CAM fabrication of temporary single-tooth abutment crowns for anterior and posterior restorations on titanium bases for a clinical wearing period of up to one year.

Requirements for this indication:

- Reliable adhesive bonding between the VITA CAD-Temp abutment crown and the titanium base:
 - Sufficient geometry (diameter, height) of the titanium base.
 - Avoidance of sharp margins and edges.
- Cervical support of the abutment crown on the titanium base:
 - Chamfer or rectangular step with rounded inner angle and at minimum width of 0.6 mm.

Note: The processing instructions for the titanium bases are general recommendations that are not subject to a guarantee. In the event of questions regarding which titanium base is suitable for the indication in question, please contact your implant manufacturer. Strict compliance with the processing instructions provided by the manufacturer for the recommended bonding materials is crucial to clinical success.

Overview of indications - VITA CAD-Temp IMPLANT SOLUTIONS

Temporary anterior and posterior tooth abutment crowns

Contraindication

- Superstructures with a highly asymmetric design and elongated extensions are contraindicated for reasons of structural stability.
- Permanent restorations
- Not for patients with PMMA intolerance


Product application/control

Geometric requirements for the titanium base

Note: Please observe the information provided by the respective manufacturers with regard to using the titanium bases.

The following points need to be observed:


- The diameter or the size must be adapted to the clinical situation and to the implant system selected.
- The titanium base may only be processed if the specifications of the respective manufacturer are observed.
- The geometric requirements for safe structural stability of the crown restoration must be fulfilled.

Dimensions of the titanium base

- Height of adhesive surface: at least 4.0 mm
- Marginal shoulder width: at least 0.6 mm

Geometric requirements - VITA CAD-Temp IS

- Occlusal wall thickness: at least 1.0 mm
- Circumferential wall thickness around the titanium base: at least 0.8 mm, max. 6.0 mm.
- Marginal shoulder width: at least 0.6 mm

* Note: Maximum height of the abutment crown* = twice the height of the TiBase** + max. 2 mm.

Compatibility of titanium bases

The current country-specific lists of titanium bases approved for VITA IMPLANT SOLUTIONS blanks can be found on the homepage under the following links:

www.vita-zahnfabrik.com/VITA_CADTemp_IS#titan_compatibility

Processing

Designing the restoration

- The single-element VITA CAD-Temp crown is constructed using a VITA CAD-Temp monoColor IS-16 block with CEREC 4.4 software or inLab 15.0 software or higher.
- A temporary VITA CAD-Temp crown is then milled in the desired shade from a VITA CAD-Temp IS-16 block (S or L), depending on the corresponding connection size.

Retworking of the VITA CAD-Temp abutment crown (extraoral)

- After the milling / grinding process, separate the sprue from the block using a diamond cutting disc and then remove using a fine-cut tungsten carbide bur while taking the shape of the emergence profile and the proximal contact points into consideration. Please ensure that at least the minimum layer thickness is observed.

Verification of the fit between the titanium base and the VITA CAD-Temp abutment crown

- Fix the TiBase on the laboratory analog and carefully place the abutment crown on the titanium base, while taking the anti-rotation lock (groove) into account, and checking the fit closely.

Polishing of the abutment crown

Prepolishing:

- Use a suitable silicone polisher and a small goat-hair brush.
- Use a low rotational speed: (< 5,000 rpm).

High-gloss polishing:

- High-gloss polishing is done using commercially available acrylic polishing pastes, such as VITA Polish Hybrid.

Note: Careful polymerization and polishing are essential requirements in order to obtain an optimum result and prevent the formation of deposits and gingivitis.

Cleaning

- In the ultrasonic unit

Extraoral bonding of the titanium base with the VITA CAD-Temp abutment crown

- Precise and careful preparation of the bonding surfaces is required in order to achieve an optimum adhesive bond between the titanium base and the VITA CAD-Temp abutment crown.

Note: The diameter of the titanium base may not be reduced, for example, by grinding.

- To protect the connection surface between the implant and the titanium base, fixation of the titanium base inside a laboratory analog should be carried out.
- The contact surfaces between the titanium base and the implant may not be sandblasted or processed in any other way!
- Shortening the titanium base is not recommended!
- Prior to bonding, verify whether the VITA CAD-Temp abutment crown can be easily attached to the titanium base without any gaps.
- Place the VITA CAD-Temp abutment crown on the titanium base and mark the positional relation using a waterproof pen.
- This supports correct subsequent positioning, when bonding the crown to the titanium base.
- The bonding surfaces of the VITA CAD-Temp abutment crown and the titanium base must be free from dust and grease.
- For extraoral bonding of the titanium base and the VITA CAD-Temp monoColor abutment crown, use a suitable opaque adhesive composite: Multilink Hybrid Abutment (Ivoclar Vivadent).

Conditioning the titanium base

- Protect the connection surface between the implant and the titanium base with wax or silicone. Seal the screw channel with Teflon tape. Leave some Teflon tape sticking out at the top. This allows it to be removed more easily later on.
- Use high-grade aluminium oxide (Al₂O₃) and carefully sandblast only the bonding surface of the titanium base
 - Particle size: 50 µm,
 - Pressure max. 2.0 bar until a matte surface is obtained.
- Then remove wax or silicone. Remove the Teflon tape from the screw channel. Clean the titanium base using an ultrasonic bath, with alcohol or using a steam jet, and dry using oil-free air.
- The surface to be bonded must not be touched after cleaning to avoid any contamination, which may have a negative effect on the subsequent adhesive bond.
- Application of a suitable bonding agent, such as Monobond Plus (Ivoclar Vivadent), using a disposable brush or microbrush.
- Apply Monobond Plus and allow to act for 60 seconds.
- Then dry with oil-free air.

Note: Please observe the instructions for use of the manufacturers of the corresponding products!

Then seal the screw channel again using a foam pellet or Teflon tape, prior to bonding with the abutment crown.

Cleaning of the abutment crown

- In an ultrasonic bath, remove grease using ethanol. Then dry with oil-free air.
- After cleaning, use a microbrush or disposable brush to thinly coat the adhesive surface with a primer, such as SR Connect (Ivoclar Vivadent) that contains MMA (methyl methacrylate), allow to act for 30 seconds and then cure with a polymerization device. Please observe the manufacturer's instructions!

Permanent extraoral bonding of the titanium base with the VITA CAD-Temp abutment crown

Note: For bonding the titanium base to the abutment crown, please use a suitable opaque adhesive composite:

Multilink Hybrid Abutment (Ivoclar Vivadent).

- Application of Multilink Hybrid Abutment (Ivoclar Vivadent) to the titanium base and to the adhesive surface of the VITA CAD-Temp abutment crown. Slide the VITA CAD-Temp abutment crown slightly, rotating back and forth onto the titanium base to approx. two-thirds. This way, uniform wetting of the two adhesive surfaces is ensured. Align both components in such a way that the positional markings are matched.
- Carefully slide the VITA CAD-Temp abutment crown on to the marginal final position so that the rotation and position lock of the titanium base engages in the groove of the abutment crown interface.
- Press the parts firmly together for five seconds.
- Check correct positional relationship in the final position:
- Smooth transition (no gaps) between the crown and titanium base!
- Remove the foam pellet from the screw channel.
- Remove any excess material in the screw channel using a microbrush.

Important: Do not remove excess adhesive in the circumferential area before the hardening phase two to three minutes after mixing. During this step, fix the parts by applying slight pressure.

For final curing of the adhesive composite, apply glycerin gel (e.g., VITA ADIVA OXY-PREVENT) at the joint gap between VITA CAD-Temp and the titanium in order to prevent an oxygen inhibited layer.

Curing time (autopolymerization): at least seven minutes.

Important: Fix the restoration to be bonded on the titanium base until autopolymerization is completed and do not move it.

If bonding material residue is present in the screw channel, remove it using suitable rotary instruments. Do not damage the titanium base!

Note: Please observe the instructions for use provided by the manufacturers of the respective products!

Polishing of the adhesive joint

- Prepolish the adhesive joint carefully using a suitable silicone polisher and a small goat-hair brush.
- Use a low rotational speed: (< 5,000 rpm).
- Subsequent high-gloss polishing using acrylic polishing pastes, such as VITA Polish Hybrid.
- Avoid generation of excessive heat!

Optional: Individualization

Individualization of the emergence profile of VITA CAD-Temp abutment crowns with VITA VM LC flow, veneering composite or restorative composites is possible.

Please see Instructions for use / full version.

Polishing

• Prior to bonding to the abutment crown, polish proximal areas outside the mouth. Prepolishing:

- Use a suitable silicone polisher and a small goat-hair brush.
- Use a low rotational speed: (< 5,000 rpm).

High-gloss polishing:

- Using commercially available acrylic polishing pastes, such as VITA Polish Hybrid.
- Avoid generation of excessive heat!

Note: Careful polymerization and polishing are essential requirements in order to obtain an optimum result and prevent the formation of deposits and gingivitis.

Overview of the process steps for extraoral bonding of the VITA CAD-Temp crown with a titanium base

| Process steps | VITA CAD-Temp inter- face Abutment crown | Titanium base |
|---|--|--|
| 1. Sandblasting with Al ₂ O ₃ | 50 µm, max. 2.0 bar | 50 µm, max. 2.0 bar |
| 2. Cleaning the surface | Ethanol, allow to evaporate | Ultrasonic unit, ethanol, allow to evaporate |
| 3. Surface conditioning | Apply SR Connect, allow to act for 30 seconds and polymerize subsequently | Apply Monobond Plus, allow to act for 60 seconds and blow with air |
| 4. Adhesive bonding | Multilink Hybrid Abutment Curing time (autopolymerization): at least seven minutes* | |
| 5. Coverage of the adhesive joint | Glycerin gel (VITA ADIVA OXY-PREVENT) | |
| 6. Polishing of the adhesive joint | Standard commercially-available polisher for composites / acrylics | |

* Important: Fix the restoration to be bonded on the titanium base until autopolymerization is completed and do not move it.

Disinfection/sterilization

Prior to seating, it is recommended to disinfect VITA CAD-Temp crowns with an ethanol-containing surface disinfectant or to sterilize them. The local legal regulations and hygiene standards applicable to dental practices need to be observed.

Steam sterilization can be carried out with a triple fractionated prevacuum in compliance with the following parameters:

- Sterilization time: 10 min
- Steam temperature 134° C/ 273° F

The abutment crown needs to be seated immediately after sterilization and may not be stored temporarily!

Please note:

The user/dentist is responsible for the sterility of the VITA CAD-Temp abutment crown. It must be ensured that only suitable equipment and materials, as well as product-specific validated methods, are used for sterilization. The equipment that is used should be maintained properly and serviced regularly.

Intraoral seating of the abutment crown on the implant

Note: For screw attachment with the implant, please use the tools provided by the implant manufacturer and observe the specified torques.

- Fixation of the crown on the implant.
- Manual insertion of the corresponding implant screw.
- Tightening of the implant screw using a torque wrench.
- Observe the specifications provided by the manufacturer!
- Ensure that the treatment area and screw channel remain dry.
- Insert a sterile cotton pellet, foam pellet or Teflon tape into the screw channel using a spherical plugger.
- Apply a bonding agent compatible with the restorative composite onto the inner surfaces of the screw channel.
- Seal the screw channel with temporary restorative material or restorative composite, depending on the planned time in situ, in the appropriate shade.
- Verification of the proximal and occlusal contact points.

Overview of the process steps for intraoral sealing of the VITA CAD-Temp screw channel with filling composite

| Process steps | | VITA CAD-Temp Abutment crown |
|---------------|---|---|
| 1. | Protect the implant screw head | Foam pellet, Teflon tape |
| 2. | Application of bonding agent in screw channel | Bonding agent |
| 3. | Adhesive filling | Temp. filling material or filling composite* |
| 4. | Adhesive filling | Standard commercially-available polisher, polishing paste for composites / acrylics |

*Depending on the clinical time in situ of the crown

Fine morphological adjustments

- The occlusion must be completely free of interferences. Remove unwanted occlusal contacts using diamond burs (40 µm).
- Finishing and polishing (intraoral)
- Prepolish the VITA CAD-Temp abutment crown with suitable silicone polishers.
- Use a low rotational speed: (< 5,000 rpm).
- Then polish to a high-gloss finish using standard commercially-available acrylic polishing agents that are also suitable for intraoral use.
- Avoid generation of excessive heat!
- Pay attention to margins and contact points when finishing and polishing the restoration.

Detailed information on processing can be found in the document "Instructions for Use | full version" at www.vita-zahnfabrik.com/de/Dentist-Solutions/CAD/CAM.

Recommended product combinations

For recommended product combinations with the VITA system components and corresponding processing recommendations and restrictions regarding the combination possibilities of the product, please refer to the Instructions for Use | full version.

Storage/disposal

Disposal via household waste.

The products labelled with a pictogram for hazardous substances are to be disposed of as hazardous waste. Recyclable waste (such as attachments, paper and plastics) must be disposed of using appropriate recycling systems. If necessary, contaminated product residues should be pretreated in accordance with regional regulations and disposed of separately.

Technical data

Physical properties

| Properties | Unit of measure | Value |
|-----------------------------|---|--------------|
| Flexural Strength | MPa (Nmm ⁻²) | > 80 |
| Modulus of elasticity | MPa (Nmm ⁻²) | Approx. 2800 |
| Softening temperature (DSC) | °C | approx. 118 |
| Inorganic filler content | Wt% | approx. 14 |
| Water absorption | complies with EN ISO 10477 Polymer-based crown and bridge materials | |
| Solubility | complies with EN ISO 10477 Polymer-based crown and bridge materials | |
| Shade stability | complies with EN ISO 22112 Artificial teeth for dental prostheses | |

Chemical composition

| Component | by weight % |
|-------------------------|-------------|
| Polymethyl methacrylate | 82 - 86 |
| SiO ₂ | 14 - 17 |
| Pigments | < 1 |

Safety information

Please refer to the corresponding safety data sheets for detailed information. The corresponding safety data sheets can be downloaded at www.vita-zahnfabrik.com/sds.

When working with the product, wear suitable safety goggles/face protection, gloves and safety clothing. In case of formation of dust, use an extraction system or wear a face mask.

Product reliability

Information on reporting serious incidents in connection with medical devices, general risks associated with dental treatments, residual risks and (if applicable) short clinical safety and performance reports (SSCPs) can be found at https://www.vita-zahnfabrik.com/product_safety


Disclaimer

Please note: Our products must be used in accordance with the instructions for use. We do not accept liability for any damage resulting from incorrect handling or usage. The user is furthermore obliged to check the product before use with regard to its suitability for the intended area of applications. We cannot accept any liability if the product is used in conjunction with materials and equipment from other manufacturers that are not compatible or not authorized for use with our product and this results in damage. The VITA Modulbox is not necessarily a component of the product. Date of issue of this information for use: 2023-02

After the publication of this information for use, any previous versions become obsolete. The current version can be found at www.vita-zahnfabrik.com.

VITA Zahnfabrik has been certified and the following product bears the CE mark:

CE 0124**VITA CAD-Temp® IS****Rx only (only for professional use)**

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VITA Zahnfabrik H. Rauter GmbH & Co.KG, Bad Säckingen (Germany)
Zweigniederlassung Basel c/o Perrig AG, Max Kämpf-Platz 1, 4058 Basel
VITA Zahnfabrik · H. Rauter GmbH & Co. KG