VITA ENAMIC® IS

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VITA ENAMIC® IS

VITA ENAMIC IS (Implant Solutions) is a hybrid ceramic for the fabrication of abutment crowns and mesostructures on implants.

Intended purpose

VITA ENAMIC products are ceramic materials for dental treatments.

Patient target group

No restrictions.

Intended user

Dental professionals only: dentists and dental technicians, Rx only

Indication



VITA ENAMIC IMPLANT SOLUTIONS is indicated for permanent mesostructures and abutment crowns on titanium bases.

Requirements for this indication:

- Reliable adhesive bonding between the VITA ENAMIC abutment crown and the titanium base:
 - Sufficient geometry (diameter, height) of the titanium base.
 - Avoidance of sharp margins and edges.
- Cervical support of the abutment crown on the titanium base:
 - Chamfer or rectangular step with rounded inner angle and a minimum width of 0.6 mm.

Note: The processing instructions for the titanium bases are general recommendations that are not subject to a guarantee. In the event of questions regarding which titanium base is suitable for the indication in question, please contact your implant manufacturer. Strict compliance with the processing instructions provided by the manufacturer for the recommended bonding materials is crucial to clinical success.

Indication: VITA ENAMIC IMPLANT SOLUTIONS

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Contraindication

- Superstructures with a highly asymmetric design and elongated extensions are contraindicated for reasons of structural stability.
- Temporary bonding of a VITA ENAMIC crown on a VITA ENAMIC mesostructure
- Monolithic abutment bridges
- Free-end restorations
- Parafunction (for example bruxism)


Product application/control

Geometric requirements for the titanium base

Note: Please observe the information provided by the respective manufacturers with regard to using the titanium bases.

The following points need to be observed:


- The diameter or the size must be adapted to the clinical situation and to the implant system selected.
- The titanium base may only be processed if the specifications of the respective manufacturer are observed.
- The geometric requirements for safe structural stability of the crown restoration must be fulfilled.

Dimensions of the titanium base

- Height of adhesive surface: at least 4.0 mm
- Marginal shoulder width: at least 0.6 mm

Geometrical requirements - VITA ENAMIC IS abutment crown


- Occlusal wall thickness: at least 1.0 mm
- Circumferential wall thickness around the titanium base: at least 0.8 mm, max. 6.0 mm.

- Marginal shoulder width: at least 0.6 mm

* Note: Maximum height of the abutment crown* = twice the height of the TiBase** + max. 2 mm.

Geometrical requirements - VITA ENAMIC IS mesostructure

- Maximum height of the mesostructure* = twice the height of the TiBase
- Wall thickness of the mesostructure around the screw channel: min. 0.8 mm
- Crown width with regard to the screw channel of the mesostructure - circumferential: max. 6.0 mm

Processing

Designing the restoration

- The VITA ENAMIC abutment crown or the ENAMIC mesostructure is constructed using a VITA ENAMIC IS-16 block with CEREC 4.4 software or inLab 15.0 software or higher.
- An abutment crown or a mesostructure is then milled in the desired shade from a VITA ENAMIC IS-16 block (abutment crown) or IS-14 block (mesostructure), (S or L), depending on the corresponding connection size.

Reworking of the VITA ENAMIC abutment crown or the mesostructure (extraoral)

- After the milling process, separate the sprue from the block using a diamond-coated cutting disc, and then grind with a diamond instrument. Please ensure that at least the minimum layer thickness is observed.

Verification of the fit between the titanium base and the restoration

- Fix the titanium base on the laboratory analog and carefully place the restoration (abutment crown, mesostructure) on the titanium base, while taking the anti-rotation lock (groove) into account, and checking the fit closely.

Polishing the restoration

Prepolishing:

- Using a suitable polisher from the VITA ENAMIC Polishing Set.
- Use a low rotational speed: (< 7,000 rpm).

High-gloss polishing:

- High-gloss polishing is done using commercially available acrylic polishing pastes, such as VITA Polish Hybrid.

Note: Careful polishing is an essential requirement in order to obtain an optimum result and prevent the formation of deposits and gingivitis.

Cleaning

- In the ultrasonic unit

Extraoral bonding of the titanium base with the restoration

- Precise and careful preparation of the bonding surfaces is a prerequisite in order to achieve an optimum adhesive bond between the titanium base and the VITA ENAMIC restoration.

Note: The diameter of the titanium base may not be reduced, for example, by grinding.

- To protect the connection surface between the implant and the titanium base, fixation of the titanium base inside a laboratory analog should be carried out.
- The contact surfaces between the titanium base and the implant may not be sandblasted or processed in any other way!
- Shortening the titanium base is not recommended!
- Prior to bonding, verify whether the VITA ENAMIC restoration can be easily attached to the titanium base without any gaps.
- Place the VITA ENAMIC restoration on the titanium base and mark the positional relationship using a waterproof pen.
- This supports correct subsequent positioning, when bonding the crown to the titanium base.
- The bonding surfaces of the VITA ENAMIC restoration and the titanium base must be free from dust and grease.
- For extraoral bonding of the titanium base and the VITA ENAMIC restoration, use a suitable opaque adhesive composite: Multilink Hybrid Abutment (Ivoclar Vivadent).

Conditioning the titanium base

- Protect the connection surface between the implant and the titanium base with wax or silicone. Seal the screw channel with Teflon tape. Leave some Teflon tape sticking out at the top. This allows it to be removed more easily later on.
- Use high-grade aluminium oxide (Al₂O₃) and carefully sandblast only the bonding surface of the titanium base
 - Particle size: 50 µm,
 - Pressure max. 2.0 bar until a matte surface is obtained.
- Then remove wax or silicone. Remove the Teflon tape from the screw channel. Clean the titanium base using an ultrasonic bath, alcohol or steam jet, and dry with oil-free air.
- The surface to be bonded must not be touched after cleaning to avoid any contamination, which may have a negative effect on the subsequent adhesive bond.
- Application of a suitable bonding agent, such as Monobond Plus (Ivoclar Vivadent), using a disposable brush or microbrush.
- Apply Monobond Plus and allow to act for 60 seconds.
- Then dry with oil-free air.

Note: Please observe the instructions for use of the manufacturers of the corresponding products!

Then seal the screw channel again using a foam pellet or Teflon tape, prior to bonding with the restoration.

Conditioning of the VITA ENAMIC abutment crown/mesostructure

- Cover any polished outer surface, as needed, in order to avoid accidental etching.
- Apply VITA ADIVA CERA-ETCH (hydrofluoric acid gel, 5%) to the surfaces to be bonded.

Etching time: 60 seconds
 - Completely remove any remaining acid by using H₂O spray or clean in the ultrasonic bath.
 - Then dry with oil-free air. After drying, the etched surfaces have a whitish opaque appearance.
 - Apply silane bonding agent (such as Monobond Plus from Ivoclar Vivadent) to the etched surface.
 - Allow the silane bonding agent to act and let it dry completely.
 - After this pretreatment, it is essential to avoid all contamination of the bonding surface in order to prevent any negative impact on the adhesive bonding.

Note: Please observe the instructions for use of the manufacturers of the corresponding products!

Optional: Individualization

Individualization of the emergence profile of VITA ENAMIC abutment crowns with VITA VM LC flow, veneering composite or restorative composites is possible. Please see Instructions for use / full version.

Permanent extraoral bonding of the titanium base with the VITA ENAMIC abutment crown/mesostructure

- Precise and careful preparation of the bonding surfaces is a prerequisite in order to achieve an optimum adhesive bond between the titanium base and the VITA ENAMIC abutment crown or mesostructure.

Note: In order to bond the titanium base with the abutment crown, please use a suitable opaque adhesive composite based on methacrylate: Multilink Hybrid Abutment (Ivoclar Vivadent).

Note: Please observe the instructions for use of the manufacturers of the corresponding products!
- Application of the adhesive composite Multilink Hybrid Abutment to the titanium base and to the adhesive surface of the VITA ENAMIC restoration.

- Slide the VITA ENAMIC restoration slightly, rotating back and forth onto the titanium base to approx. two-thirds. This way, uniform wetting of the two adhesive surfaces is ensured. Align both components in such a way that the positional markings are matched.
- Carefully slide the VITA ENAMIC restoration on to the marginal end position so that the rotation and position lock of the titanium base engages in the groove of the restoration interface.
- Press the parts firmly together for five seconds.
- Check correct positional relationship in the final position:
 - Smooth transition (no gaps) between the restoration and titanium base!
 - Remove the foam pellet from the screw channel.
 - Remove any excess material in the screw channel using a microbrush.

Important: In the circumferential area, only remove excess adhesive in the hardening phase two to three minutes after mixing. During this step, fix the parts by applying slight pressure.

For final curing of the adhesive composite, apply glycerin gel (e.g., VITA ADIVA OXY-PREVENT) at the joint gap between VITA ENAMIC and the titanium in order to prevent an oxygen inhibited layer. Curing time (autopolymerization): at least seven minutes.

Important: Fix the restoration to be bonded on the titanium base until autopolymerization is completed and do not move it.

If bonding material residue is present in the screw channel, remove it using suitable rotary instruments. Do not damage the titanium base!

Note: Please observe the instructions for use provided by the manufacturers of the respective products!

Polishing of the adhesive joint

Using the instruments provided with the VITA ENAMIC Polishing Set and at a low speed (< 5,000 rpm), carefully prepolish the adhesive joint and polish to a high-gloss finish.

- Additionally, a goat hair brush and diamond polishing paste (VITA Polish Hybrid) can be used for final high-gloss polishing. Then repolish with a dry cotton buff.

Overview of the process steps for extraoral bonding of the VITA ENAMIC abutment crown/mesostructure with a titanium base

	Process steps	Interface Abutment crown/ Mesostructure	Titanium base
1.	Sandblasting with high-grade aluminum oxide (Al ₂ O ₃)	–	50 µm, max 2.0 bar
2.	Cleaning the surface	Ethanol, allow to evaporate	Ultrasonic unit, ethanol (allow to evaporate)
3.	Etching (extraoral)	VITA ADIVA CERA-ETCH, 5% HF, 60 s	–
4.	Cleaning the surface	Use water, either by spraying off or in an ultrasonic bath	–
5.	Conditioning/ silanizing	Apply Monobond Plus, allow to act for 60 seconds and blow with air	Apply Monobond Plus, allow to act for 60 seconds and blow with air
6.	Adhesive bonding	Multilink Hybrid Abutment Curing time (autopolymerization): at least seven minutes*	
7.	Coverage of the adhesive joint	VITA ADIVA OXY-PREVENT	
8.	Polishing of the adhesive joint	VITA ENAMIC POLISHING SET	

* Important: Fix the restoration to be bonded on the titanium base until autopolymerization is completed and do not move it.

Note: Please use the above adhesive products in accordance with the manufacturer's system requirements. For example, Monobond Plus (Ivoclar Vivadent) only in combination with Multilink Hybrid Abutment (Ivoclar Vivadent).

Disinfection/sterilization

Prior to seating, it is recommended to disinfect the VITA ENAMIC restoration with an ethanol-containing surface disinfectant, or to sterilize it. The local legal regulations and hygiene standards applicable to dental practices need to be observed. Steam sterilization can be carried out with a triple fractionated prevacuum in compliance with the following parameters:

- Sterilization time: 10 min
- Steam temperature 134° C/ 273° F

The restoration needs to be seated immediately after sterilization and may not be stored temporarily!

Please note:

The user/dentist is responsible for the sterility of the VITA ENAMIC restoration. It must be ensured that only suitable equipment and materials, as well as product-specific validated methods, are used for sterilization. The equipment that is used should be maintained properly and serviced regularly.

Intraoral seating of the restoration on the implant

Note: For screw attachment with the implant, please use the tools provided by the implant manufacturer and observe the specified torques.

- Fixation of the restoration on the implant.
- Manual insertion of the corresponding implant screw.
- Tightening of the implant screw using a torque wrench.
Observe the specifications provided by the manufacturer!
- Ensure that the treatment area and screw channel remain dry.
- Insert a sterile cotton pellet, foam pellet or Teflon tape into the screw channel using a spherical plugger.
- Apply a bonding agent compatible with the restorative composite onto the inner surfaces of the screw channel.
- Seal the screw channel with temporary restorative composite.

Fine morphological adjustments

- The occlusion must be completely free of interferences. Remove unwanted occlusal contacts using diamond burs (40 µm).

Finishing and polishing (intraoral)

- Prepolish the VITA ENAMIC abutment crown with suitable silicone polishers.
- Use a low rotational speed: (< 5,000 rpm).
- Then polish to a high-gloss finish using standard commercially-available acrylic polishing agents that are also suitable for intraoral use.
- Avoid generation of excessive heat!
- Pay attention to margins and contact points when finishing and polishing the restoration.

Recommended product combinations

For recommended product combinations with the VITA system components and corresponding processing recommendations and restrictions regarding the combination possibilities of the product, please refer to the Instructions for Use | full version.

Storage/disposal

Disposal via household waste.

The products labelled with a pictogram for hazardous substances are to be disposed of as hazardous waste. Recyclable waste (such as attachments, paper and plastics) must be disposed of using appropriate recycling systems. If necessary, contaminated product residues should be pretreated in accordance with regional regulations and disposed of separately.

Technical data

Physical properties

Property	Unit	Value
Flexural Strength (ISO 6872)	MPa	150–160
Fracture toughness	MPa√m	1.5
Modulus	GPa	30
Hardness	GPa	2.5

Chemical composition of the fine-structure feldspar ceramic network*

Oxides	by weight %
SiO ₂	58–63
Al ₂ O ₃	20–23
Na ₂ O	6–11
K ₂ O	4–6
B ₂ O ₃	0.5–2
CaO	< 1
TiO ₂	< 1

* The values of the chemical composition listed above are dependent on the lot. Chemical elements (oxides) that are contained in very low concentrations and required (e.g., for coloring) are not listed.

Chemical composition of the polymer network

The polymer network consists of methacrylate polymer.

Safety information

Please refer to the corresponding safety data sheets for detailed information. Safety data sheets can be downloaded at www.vita-zahnfabrik.com/sds.

When working with the product, wear suitable safety goggles/face protection, gloves and safety clothing. In case of formation of dust, use an extraction system or wear a face mask.

Product reliability

Information on reporting serious incidents in connection with medical devices, general risks associated with dental treatments, residual risks and (if applicable) short clinical safety and performance reports (SSCPs) can be found at https://www.vita-zahnfabrik.com/product_safety

Disclaimer


Please note: Our products must be used in accordance with the instructions for use. We do not accept liability for any damage resulting from incorrect handling or usage. The user is furthermore obliged to check the product before use with regard to its suitability for the intended area of applications. We cannot accept any liability if the product is used in conjunction with materials and equipment from other manufacturers that are not compatible or not authorized for use with our product and this results in damage. The VITA Modulbox is not necessarily a component of the product. Date of issue of this information for use: 2023-02

After the publication of this information for use, any previous versions become obsolete. The current version can be found at www.vita-zahnfabrik.com.

VITA Zahnfabrik has been certified and the following product bears the CE mark:

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VITA Zahnfabrik H. Rauter GmbH & Co.KG, Bad Säckingen (Germany)
Zweigniederlassung Basel c/o Perrig AG, Max Kämpf-Platz 1, 4058 Basel
VITA Zahnfabrik · H. Rauter GmbH & Co. KG