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4 °C / 39 °F 25 °C / 77 °F

Please read the instructions for use before using the product.

VITA AKZENT LC is a light-curing methacrylate-based stain/glaze system for extraoral surface characterization of dental restorations made of hybrid ceramic, resin veneering materials, CAD/CAM composites, prefabricated teeth, denture base resins and 3D printing acrylic polymers. It can also be used for internal characterization with the layering technique of veneering composites.

**Intended use**

VITA AKZENT LC products are unfilled resin sealant and coating materials for dental treatments.

**Patient target group**

No restrictions

**Intended user**

Dental professionals only:  
dentists and dental technicians. (Rx Only)

**Indication**

Indication range:


- Restorations made of hybrid ceramic (VITA ENAMIC)
- Restorations made of light-curing veneering material (e.g., VITA VM LC)
- Restorations made of CAD/CAM composites (e.g., VITA CAD-Temp)
- Prefabricated teeth (e.g., VITAPAN)
- Denture bases (e.g., VITA VIONIC BASE)
- Restorations and denture bases made of 3D printing acrylic polymers

**Contraindication**

- In patients with allergies or sensitivities to the ingredients
- Do not use on occlusal contact points of restorations
- Not approved for intraoral use

**Product application/control**

**Pretreatment steps according to material**

\* These are only general recommendations for pretreatment. Depending on the type of material, the specific manufacturer's instructions must be observed.

\*\*If no inhibition layer is present during internal characterization, apply VITA VM LC MODELLING LIQUID, allow to act for 30 to max. 60 seconds and blow off. Characterize subsequently.

**Use:****Characterization / Glazing with VITA AKZENT LC****Note:**

- Shake bottle well before use (for 10s). The mixing ball must be clearly heard!
- Close bottles immediately after use.
- The brushes included in the assortment are recommended for applying the stains.
- Apply the stains in a thin layer.
- Subsequent coating with GLAZE is not absolutely necessary.

**Tip:**

- Before characterization with the EFFECT STAINS, the rough surfaces of the restoration to be characterized can first be wetted with a thin layer of GLAZE to identify its base shade. Then polymerize GLAZE together with the EFFECT STAINS.

- The chroma of the stains can be reduced by adding GLAZE.
- Polymerize individual stain layers in steps.
- Use microbrush disposable applicators or brushes for the application.
- Apply STAINS to the intended area.
- Apply GLAZE swiftly and without streaks.
- Close bottle immediately after use.

To avoid fast abrasion of the stain, VITA AKZENT LC should not be used at the occlusal contact points of restorations.

Please note the recommended processing time of approx. three minutes and clean the paint brushes with VITA AKZENT LC CLEANER after use.

**Polymerization**

**To cure VITA AKZENT LC, a light-curing device is required whose light sources emit rays in the wavelength range of < 430 nm.**

This is true of common light curing units that emit light in the wavelength range of 350 to 500 nm.

**Recommended polymerization devices/lamps**

Dental-technical polymerization devices			
Company	Polymerization device	Polymerization time*	Note
Shofu	Solidilite V	3 min	The object must be placed in the cone of light in the center of the chamber. The object must not be placed on the floor!
Shofu	Solidilite EX	4 min	see Solidilite V
DeguDent/Dentsply	Eclipse® junior VLC Curing Unit	Basic 3 = 3 min	Material group: in-joy The object must be placed in the center of the chamber.
DeguDent/Dentsply	Triad® 2000	6 min	rotating on the rotary table
Bredent	bre.Lux Power Unit 2	1 x program F1 (90 s)	–
3M ESPE	Visio Beta vario	2 x program: 7 min incl. 10 s vacuum	–
GC	Labolight Duo	3 min	The object must be placed in the cone of light high in the center of the chamber. The object must not be placed on the floor!
Hager & Werken	Speed Labolight®	3 min	The object must be placed 3 cm higher in the light cone of the halogen lamp. Caution with denture acrylics: Do not overheat!
Kulzer	HiLite® power	90 s	The object must be placed high in the chamber. The object must not be placed on the floor!
Kulzer	Heraflash	90 s	see HiLite Power
Kulzer	UniXS®	90 s	see HiLite Power
Kulzer	Dentacolor® XS	90 s	see HiLite Power
Ivoclar Vivadent	Lumatmat® 100	Heating level 0 VB = 0 VG = 10:00 min (BP = 10:00 min)	The times and heating levels must be programmed by the user! VB = Pre curing VG = Tempering BP = Resulting curing time

\* The specified polymerization times are minimum times.

**Note:**

- Dental-technical polymerization devices are preferable to dental hand lamps, as they ensure a much more even exposure of the object.

Dental polymerization lamps			
Company	Polymerization device	Polymerization time*	Note
Ivoclar Vivadent	Bluephase G2	4 x 20 s	"High" mode
Ultradent	VALO LED	4 x 20 s	"Standard" mode. The light cones of the two LEDs must completely overlap -- > very small distance to the restoration required

**Please note:**

- To ensure the complete curing of VITA AKZENT LC, the polymerization units must feature the required intensity and wavelength.
- Many dental hand lamps do not have the wavelength range required for VITA AKZENT LC. Please observe the instructions for use of the manufacturer of your device.
- The necessary light intensity and wavelength of the polymerization device must be observed to ensure complete curing.
- You can find the latest information on the recommended polymerization devices at: [www.vita-zahnfabrik.com/akzentlc](http://www.vita-zahnfabrik.com/akzentlc)

**Important information:**

- All surfaces must be completely polymerized and hard. There must be no sticky surfaces.**
- After final polymerization, polish the glazed surface with a suitable polishing paste (such as Renfert Polish hybrid materials) and a soft goat hair brush. Then use a dry cotton buff to achieve the final gloss.**

**Recommended product combinations**

For recommended product combinations with the VITA system components and corresponding processing recommendations and restrictions regarding the combination possibilities of the product, please refer to the full version of the Instructions for Use.

**Storage/disposal**

Store at 4 °C - 25 °C (39 °F - 77 °F). Storage in a refrigerator is recommended. Do not use after the expiry date!

The products labelled with a pictogram for hazardous substances are to be disposed of as hazardous waste. Recyclable waste (such as attachments, paper and plastics) must be disposed of using appropriate recycling systems. If necessary, contaminated product residues should be pretreated in accordance with regional regulations and disposed of separately.

**Technical data**

Chemical composition	Wt%
Methyl methacrylate and multifunctional methacrylates	30 – 40
Urethane(meth-)acrylates	40 – 60
Silicon dioxide	8 – 11
Ethyl-phenyl(2,4,6-trimethylbenzoyl)phosphinate	2 – 6
Other	< 1
Pigments	< 2

**Product reliability**

Information on reporting serious incidents in connection with medical devices, general risks associated with dental treatments, residual risks and (if applicable) short clinical safety and performance reports (SSCPs) can be found at

[www.vita-zahnfabrik.com/product\\_safety](http://www.vita-zahnfabrik.com/product_safety).

Safety data sheets can be downloaded at [www.vita-zahnfabrik.com/vita-akzent-lc](http://www.vita-zahnfabrik.com/vita-akzent-lc).

## Hazardous substances

### VITA AKZENT LC EFFECT STAINS/CHROMA STAINS/GLAZE

Contains: methyl methacrylate; methyl 2-methylprop-2-enoate; methyl 2-methylpropenoate  
7,7,9-Trimethyl-4,13-dioxo-3,14-dioxa-5,12-diazahexadecan-1,16-diylobismethacrylat  
2-Propenoic acid, reaction products with pentaerythritol  
Ethyl phenyl(2,4,6-trimethylbenzoyl)phosphinat  
pentaerythritol tetraacrylate  
pentaerythritol triacrylate  
2,2'-ethylenedioxydiethyl dimethacrylate

Danger

H225 Highly flammable liquid and vapour.  
H315 Causes skin irritation.  
H317 May cause an allergic skin reaction.  
H318 Causes serious eye damage.  
H335 May cause respiratory irritation.  
H412 Harmful to aquatic life with long lasting effects.

### VITA AKZENT LC CLEANER

Gefahr

H225 Highly flammable liquid and vapour.  
H319 Causes serious eye irritation.

## Precautionary statements

P210 Keep away from heat/sparks/open flames/hot surfaces.  
— No smoking.  
P233 Keep container tightly closed.  
P280 Wear protective gloves/protective clothing/eye protection/face protection.

## Disclaimer

Please note: Our products must be used in accordance with the instructions for use. We accept no liability for any damage resulting from incorrect handling or usage. The user is furthermore obliged to check the product before use with regard to its suitability for the intended area of applications. We cannot accept any liability if the product is used in conjunction with materials and equipment from other manufacturers that are not compatible or not authorized for use with our product and this results in damage. The VITA Modulbox is not necessarily a component of the product. Date of issue of this information: 2024-02

After the publication of this information for use any previous versions become obsolete. The current version can be found at [www.vita-zahnfabrik.com](http://www.vita-zahnfabrik.com).

VITA Zahnfabrik has been certified and the following products bear the CE mark:

**VITA AKZENT® LC**

**CE0124** Rx only MD

CERASMART® is a registered trademark of GC Dental Products Corp, Aichi, Japan.

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CERASMART® es una marca registrada de GC Dental Products Corp, Aichi (Japón).
CH REP VITA Zahnfabrik H. Rauter GmbH & Co.KG, Bad Säckingen (Germany)
Zweigniederlassung Basel c/o Perrig AG, Max Kämpf-Platz 1, 4058 Basel
VITA Zahnfabrik H. Rauter GmbH & Co.KG