# VITA AMBRIA® PRESS SOLUTIONS – Presstechnik

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VITA

# VITA AMBRIA® PRESS SOLUTIONS - Press technique

Please read the instructions for use before using the product.

With VITA AMBRIA PRESS SOLUTIONS, you can expect an ideally matched material system made of zirconia-reinforced lithium disilicate glass ceramic pellets and various system components for the fabrication of restorations using the press technique.

## Intended purpose

VITA AMBRIA products are ceramic materials for dental treatments.

## Patient target group

No restrictions.

## Intended user

Dental professionals only: dentists and dental technicians (Rx only).

## Indication

- Occlusal veneers (table tops)\*, onlays\*, partial crowns\*
- Inlays\*, onlays\*, partial crowns\*
- Crowns in the anterior and posterior areas
- Three-unit bridges in the anterior area up to the second premolar as a terminal abutment
- Single tooth restorations as implant suprastructures for anterior and posterior teeth
- Three-unit bridges as implant suprastructures up to the second premolar on implant abutments
- Single tooth mesostructures in the anterior and posterior areas
- Abutment crowns in the anterior and posterior areas

\*) For adhesive bonding only

## Contraindication

- in cases of inadequate preparation results (such as tangential preparation, for example)
- in cases of insufficient hard tooth substance
- in cases of insufficient space
- hyperfunction: for patients diagnosed with excessive occlusal function, in particular teeth grinders and clenchers
- restoring devitalized teeth of patients with hyperfunctions
- endodontic crowns
- posterior bridges in the area of molars
- in cases of bridges with more than three units
- inlay-retained bridges / Maryland bridges
- cantilever bridges
- in patients with allergies or sensitivities to the ingredients
- conventional or self-adhesive insertion of inlays, onlays, veneers, partial crowns and occlusal veneers (table top)
- temporary seating of restorations

## Product application/control

Workflow:

- 1) Preparation/impression taking
- 2) Model preparation (manual or CAD/CAM model preparation)
- 3) Investing
- 4) Pressing
- 5) Divesting / finishing
- 6) Power Glaze
- 7) Optional: Veneering with suitable ceramic veneering materials (with VITA LUMEX AC, for example, reduced substructure and tooth-shape supporting contouring)
- 8) Manual polishing / staining technique (with VITA AKZENT PLUS, for example)
- 9) Bonding

## Firing:

### Press parameters T / HT

Pre-dry. °C	↗ °C/min.	→ Press muffle / temp. °C, min.	↓ Pressing time 1	↓ Pressing time 2	↘ °C	Vac.
700	50	100 g, 880°C	3:00 min.	1:00 min.	-	on
		200 g, 890°C				

### Press parameters MO

Pre-dry. °C	↗ °C/min.	→ Press muffle / temp. °C, min.	↓ Pressing time 1	↓ Pressing time 2	↘ °C	Vac.
700	50	100 g, 870°C	3:00 min.	1:00 min.	-	on
		200 g, 875°C				

### Firing parameters Power Glaze

Pre-dry. °C	→ min.	↗ °C/min.	Temp. approx. °C	→ min.	↘ °C	Vac.
400	04:00	55	800°C	02:00	500	-

The press results are influenced by various parameters such as model, performance and age of the press unit used. Detailed information on processing and press process can be found in the document "Instructions for Use | Full Version" at [www.vita-zahnfabrik.com/ambria](http://www.vita-zahnfabrik.com/ambria).

## Recommended product combinations

For recommended product combinations with the VITA system components and corresponding processing recommendations and restrictions regarding the combination possibilities of the VITA AMBRIA products, please refer to the Instructions for Use | full version.

## Storage/disposal

The products labelled with a pictogram for hazardous substances are to be disposed of as hazardous waste. Recyclable waste (such as attachments, paper and plastics) must be disposed of using appropriate recycling systems. If necessary, contaminated product residues should be pretreated in accordance with regional regulations and disposed of separately. The materials must not be contaminated with foreign substances (e.g., during the press process).

## Technical data

Physical properties	Unit of measure	Value
CTE	10 <sup>-6</sup> · K <sup>-1</sup>	approx. 9.4
Biaxial strength (after pressing)	MPa	approx. 400
Biaxial strength (after Power Glaze)	MPa	approx. 550
Solubility	µg/cm <sup>2</sup>	approx. 30

Chemical composition	Wt%
ZrO <sub>2</sub> (zirconia)	8-14
SiO <sub>2</sub> (silicon dioxide)	58-66
Li <sub>2</sub> O (lithium oxide)	12-16
Pigments	< 10
Various	> 10

## Safety note

- While work is in progress, wear suitable safety goggles/face protection.
- Only perform work under an extraction unit.
- Wear light face mask when working.

## Product safety

Information on reporting serious incidents in connection with medical devices, general risks associated with dental treatments, residual risks and (if applicable) short clinical safety and performance reports (SSCPs) can be found at [www.vita-zahnfabrik.com/product\\_safety](http://www.vita-zahnfabrik.com/product_safety). The safety data sheets can be downloaded at [www.vita-zahnfabrik.com](http://www.vita-zahnfabrik.com) or requested by fax at (+49) 7761-562-299.

## Disclaimer

Please note: Our products must be used in accordance with the instructions for use. We accept no liability for any damage resulting from incorrect handling or usage. The user is furthermore obliged to check the product before use with regard to its suitability for the intended area of applications. We cannot accept any liability if the product is used in conjunction with materials and equipment from other manufacturers that are not compatible or not authorized for use with our product and this results in damage. The VITA Modulbox is not necessarily a component of the product. Date of issue of this information: 2023-06.

After the publication of these instructions for use, any previous versions become obsolete. The current version can be found at [www.vita-zahnfabrik.com](http://www.vita-zahnfabrik.com).

VITA Zahnfabrik has been certified and the following products bear the mark:

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CH REP VITA Zahnfabrik H. Rauter GmbH & Co. KG, Bad Säckingen (Germany)
Zweigniederlassung Basel c/o Perrig AG, Max Kämpf-Platz 1, 4058 Basel
VITA Zahnfabrik · H. Rauter GmbH & Co. KG