VITA NP BOND Paste

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EN

VITA NP BOND Paste

Please read instructions for use before using the product. VITA NP BOND paste is a ceramic material that balances the differing coefficients of thermal expansion between non-precious-metal alloys and ceramics (such as VITA VMK MASTER or VITA VM 13, for example). It also protects the ceramic from metal oxides emitted from the metal.

Intended purpose

VITA NP BOND products are ceramic materials for dental treatments.

Patient target group

No restrictions.

Intended user

Dental professionals only: dentists and dental technicians (Rx only).

Indication

VITA NP BOND can be used with all non-precious metal alloys such as cast iron, CAD/CAM or laser-sintered alloys. It is compatible with a CTE value (600°C) of 13.8-15.2 x 10-6 x K-1.

Contraindication

None known.

Product application/control

Framework preparation:

- Prepare the framework in accordance with the criteria specified by the alloy manufacturer.
- Depending on the manufacturer specifications, oxidation firing must be carried out.
- Once oxidation firing has been completed, the metal surface must be sandblasted with 110 µm and approx. 2-3 bar pressure.
- Afterwards, the metal surface may no longer be touched.
- Clean with pressurized air.

Mixing:

- Before using for the first time, VITA NP BOND Paste should be mixed thoroughly using the enclosed spatula.
- If VITA NP BOND Paste has dried out, it can be mixed to its original consistency again using VITA PASTE FLUID. Other liquids are not suitable.

Application:

- Using a brush, apply an opaque layer of VITA NP BOND Paste to the framework (similar to a washbake) and then perform firing in the ceramic furnace.

Firing chart:

Predry. temp. in °C	→ min.	↗ min.	↗ °C/ min.	Temp. approx. °C	→ min.	VAC min.
600	6.00	6.00	60	960	1.00	6.00

The values given are intended only as general guidelines. The firing results depend on the performance, make and age of the furnace.

A washbake is no longer necessary once VITA NP BOND Paste has been used. Once firing of the NP BOND Paste has been completed, apply OPAQUER as usual and then perform secondary firing of the ceramic.

More information on processing, as well as an application film, is available on our Internet page <https://www.vita-zahnfabrik.com/npbond>

Recommended product combinations

For recommended product combinations with the corresponding VITA system components, corresponding processing recommendations and restrictions regarding the combination possibilities of the product, please refer to the Instructions for Use | Full version.

Storage/disposal

They can be disposed of with household waste. The products labelled with a pictogram for hazardous substances are to be disposed of as hazardous waste. Recyclable waste (such as attachments, paper and plastics) must be disposed of using appropriate recycling systems. If necessary, contaminated product residues should be pretreated in accordance with regional regulations and disposed of separately. If necessary, contaminated product residues should be pretreated in accordance with regional regulations and disposed of separately.

Technical data

Chemical composition	Wt%
SiO ₂	27 – 34
Na ₂ O	2 - 3
Al ₂ O ₃	4 - 6
K ₂ O	3 - 4
CaO	1 - 2
B ₂ O ₃	2 - 3
BaO	< 1
Sn O ₂	4 - 5
ZrO ₂	3 - 4
TiO ₂	35 - 54
CeO ₂	<1
MgO	<1
Fe ₂ O ₃	<1
Others	<1

Product reliability

Information on reporting serious incidents in connection with medical devices, general risks associated with dental treatments, residual risks and (if applicable) short clinical safety and performance reports (SSCPs) can be found at https://www.vita-zahnfabrik.com/product_safety

Safety data sheets can be downloaded at www.vita-zahnfabrik.com or requested by fax at (+49) 7761-562-233.

Disclaimer


Please note: Our products must be used in accordance with the instructions for use. We do not accept any liability for damage resulting from incorrect handling or usage. The user is furthermore obligated to check the product before use with regard to its suitability for the intended area of applications.

We cannot accept any liability if the product is used in conjunction with materials or equipment from other manufacturers that are not compatible or not authorized for use with our product and this results in damage. The VITA Modulbox is not necessarily a component of the product.

Date of issue of these Instructions for Use: 2022-08
After the publication of these instructions for use, any previous versions become obsolete. The current version can be found at www.vita-zahnfabrik.com

Marking

VITA Zahnfabrik has been certified, and the following product bears the marks
VITA NP BOND Paste.

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VITA Zahnfabrik H. Rauter GmbH & Co.KG, Bad Säckingen (Germany)
Zweigniederlassung Basel c/o Perrig AG, Max Kämpf-Platz 1, 4058 Basel