VITA AMBRIA® Calibration of the press furnace

Setting the optimal press temperature with VITA AMBRIA ST, HT and T.

We recommend investing a wax grid together with the restoration for the first pressing. After that, you can use these pictures to evaluate your press results and adjust the press temperature accordingly.

- 20 °C

The grid has not flowed out, the restoration looks opaque.
The pressing temperature is far too low.
Procedure: Raise the temperature in 10°C-steps.

- 15 °C

The grid has a smooth surface, but parts of the grid have not completely flowed out.
The pressing temperature is still too low.
Procedure: Raise the temperature in 5°C-steps.

- 5 °C

The grid has flowed out almost completely.
The pressing temperature is minimally low.
Procedure: Raise the temperature slightly (approx. 5°C).

+/- 0 °C

The grid has completely flowed out, the surface is smooth and shows no reaction layer.
The pressing temperature is optimal.

+ 5 °C

The grid has flowed out, the surface shows a light reaction layer.
The pressing temperature is a bit too high.
Procedure: Minimally lower the temperature (approx. 5°C).

+ 15 °C


The grid has flowed out, the surface is whitish and partly also porous.
The reaction layer is clearly visible.
The press temperature is too high.
Procedure: Lower the temperature in 5°C-steps.

+ 20 °C

The reaction layer is very pronounced, there are small cracks and holes in the edge area.
The restoration can appear transparent / colorless.
The press temperature is too high.
Procedure: Lower the temperature in 10°C-steps.

Change in the degree of translucency in AMBRIA MO press pellets (using a 200 g investment ring as an example):

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