# VITA CERAMICS ETCH

DE

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# VITA CERAMICS ETCH

EN

VITA CERAMICS ETCH is a hydrofluoric acid gel for etching ceramic restorations.

It is used to create a retentive etched pattern on restorations made of silicate ceramic and hybrid ceramic in order to strengthen the adhesion between the ceramic surface and the adhesive composite.

- easy and precise application due to gel consistency
- easily distinguishable by its red colour

## Composition

1 ml VITA CERAMICS ETCH contains 0.047 g hydrofluoric acid (HF)

## Area of application

Silicate ceramic restorations (glass ceramic, feldspar ceramic) and hybrid ceramic (VITA ENAMIC).

## Contraindications

**VITA CERAMICS ETCH must not be used in the mouth or for etching tooth enamel.**

Oxide ceramic (e.g. zirconium oxide), since restorations made of these ceramics are cemented using a modified adhesive technique without ceramic etching.

## Hazard information / Safety advices

For extraoral use only! Contains hydrofluoric acid. H301 Toxic if swallowed. H310 Fatal in contact with skin. H314 Causes severe skin burns and eye damage. H332 Harmful if inhaled. P260 Do not breathe dust/fume/gas/mist/vapours/spray. P405 Store locked up. P280 Wear protective gloves/protective clothing/eye protection/face protection/hearing protection. P303+P361+P353 IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water or shower. P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. P310 Immediately call a POISON CENTER/doctor. Wear protective gloves/protective clothing/safety goggles. Keep locked up. This material and its container must be disposed of as hazardous waste.

## Application

Apply VITA CERAMICS ETCH to the surfaces requiring etching with a small disposable

brush or a plastic disposable applying needle.

When using the syringe with the application tip: Screw a red application tip tightly onto the VITA CERAMICS ETCH syringe. The tip enables controlled extrusion from the syringe.

Prior to application on the ceramic restoration, first apply a small amount of VITA CERAMICS ETCH to the mixing pad to ensure even application and in order to become familiar with the flow behavior.

In the event of strong resistance to application, do not continue. Instead, attach a new tip and test again. Dispose of used application tips immediately in accordance with statutory regulations in order to prevent contact.

## Etching time

- Feldspar ceramic (e.g. VITABLOCS Mark II) 60 s
- Hybrid ceramic (VITA ENAMIC) 60 s
- Zirconia reinforced lithium silicate ceramic (VITA SUPRINITY PC) 20 s

Completely remove any acid remaining on the etched surface by intensive spraying with a lot of water or steam or by cleaning in an ultrasonic bath in distilled water.

Do not brush off, since this would lead to severe impurities on the surface!

Dry well before inserting, use silane coupling agent (e.g. VITA ADIVA C-PRIME) according to manufacturer's instructions.

## Important

**The ceramic restoration should be etched immediately before insertion or after the try-in, in order to avoid contamination of the etched ceramic surface. The etched surface must no longer be touched!**

## Storage

Store at temperatures between 2°C and 28°C. Protect from direct sunlight.

## Waste disposal

The empty packing should be handed over to hazardous waste disposers.

## Labelling

The expiry date  and the batch no. **(LOT)** are clearly marked on the containers.

## Presentation

- VITA CERAMICS ETCH, dropper bottle, 6 ml
- VITA CERAMICS ETCH, disposable syringe, 3 ml, 20 plastic application needles

## Danger

UFI: J600-604V-K007-5NU5

H301 Toxic if swallowed.

H310 Fatal in contact with skin.

H314 Causes severe skin burns and eye damage.

H332 Harmful if inhaled.

VITA CERAMICS ETCH contains <5 Vol.% hydrofluoric acid (HF)

VITA CERAMICS ETCH contains <10 Vol.% sulphuric acid (H<sub>2</sub>SO<sub>4</sub>)

The safety datasheet can be downloaded at [www.vita-zahnfabrik.com](http://www.vita-zahnfabrik.com)

The product was developed for dental use and must be processed according to the instructions for use. The information contained herein does not constitute an assurance of properties and are not binding. These instructions issued: 2024-06

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