VITA VMK Master®

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VITA VMK Master®

Please read the instructions for use before using the product.
VITA VMK Master is a veneering ceramic for metal substructures in the conventional CTE range.

Intended purpose

VITA VMK Master products are ceramic materials for dental treatments.

Patient target group

No restrictions.

Intended user

Dental professionals only: dentists and dental technicians (Rx only).

Indication

VITA VMK Master is a feldspar veneering ceramic for metal substructures in the usual CTE range of 13.8 - 15.2 10⁻⁶ K⁻¹.

Contraindication

- For substructures not within the recommended CTE range
- For parafunctions (e.g., bruxism)
- If minimum layer thicknesses of the ceramic can not be adhered to
- In cases of inadequate oral hygiene

Product application/control

Workflow:

Please adhere to the following workflow when using this product:

1. Substructure preparation
2. Oxide firing in accordance with the specifications of the alloy manufacturer
3. Wash application
4. Opaque firing
5. Application of DENTINE
6. Application of ENAMEL
7. Application of corrective material (optional)
8. Adjustments of the shape
9. Characterization / glaze firing

Firing:

Preheat temp. °C	→ min.	↗ °C/min.	approx. temp °C**	→ min.	↘ °C	Vac.
500	6.00	80	Washbake 960	1.00	--	on
			Opaque firing 950	1.00	--	on
		55	First dentine firing 930	1.00	--	on
		80	Glaze firing VITA AKZENT PLUS GLAZE Paste 920	1.00	--	--

Recommended product combinations

For recommended product combinations with corresponding VITA system components, corresponding processing recommendations and restrictions regarding the combination possibilities of the product, please refer to the instructions for use | full version.

Storage/disposal

They can be disposed of with household waste. The products labelled with a pictogram for hazardous substances are to be disposed of as hazardous waste. Recyclable waste (such as attachments, paper and plastics) must be disposed of using appropriate recycling systems. If necessary, contaminated product residues should be pretreated in accordance with regional regulations and disposed of separately.

Technical data

Physical properties	Unit of measure	Value
Solubility in acids	µg/cm ²	~ 5
3-point flexural strength	MPa	~ 90
Coefficient of thermal expansion (CTE)	10 ⁻⁶ K ⁻¹	13.2 - 13.7

Chemical composition	Wt%
SiO ₂	30 - 70
Al ₂ O ₃	7 - 17
K ₂ O	5 - 13
Na ₂ O	3 - 8
BaO	≤ 4
B ₂ O ₃	≤ 2
CaO	0 - 3
SnO ₂	≤ 2
CeO ₂	≤ 21
ZrO ₂	≤ 20
TiO ₂	≤ 14
Fe ₂ O ₃	≤ 5
Sb ₂ O ₃	< 0,5
MgO	< 1
ZnO	< 0,5
others	< 1

Components	Wt%	
	Pastes	Powders
Ceramic Powder	60 – 70	100
Butanediol	20 – 27	–
Glycerol	6 – 9	–
Additional components	3 – 4	–

Hazardous substances

VITA VM OPAQUE FLUID

Danger

H314 Causes severe skin burns and eye damage.

Precautionary statements

VITA VM OPAQUE FLUID

P260 Do not breathe dust/fume/gas/mist/vapours/spray.

P264 Wash ... thoroughly after handling.

P280 Wear protective gloves/protective clothing/eye protection/face protection.

P303+P361+P353 IF ON SKIN (or hair): Remove/Take off immediately all contaminated clothing. Rinse skin with water/shower

P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

P310 Immediately call a POISON CENTER or doctor/physician.

Safety at work and health protection

- In case of formation of dust, an extraction system or a dust protection mask must be used.

Product reliability

Information on reporting serious incidents in connection with medical devices, general risks associated with dental treatments, residual risks and (if applicable) short clinical safety and performance reports (SSCPs) can be found at https://www.vita-zahnfabrik.com/product_safety

The safety data sheets can be downloaded at www.vita-zahnfabrik.com or requested by fax at (+49) 7761-562-233.

Disclaimer

Please note: Our products must be used in accordance with the instructions for use. We do not accept any liability for damage resulting from incorrect handling or usage. The user is furthermore obliged to check the product before use with regard to its suitability for the intended area of applications. We cannot accept any liability if the product is used in conjunction with materials or equipment from other manufacturers that are not compatible or not authorized for use with our product and this results in damage. The VITA Modulbox is not necessarily a component of the product.

Date of issue of this information for Use: 2024-02

After the publication of this information for use, any previous versions become obsolete. The current version can be found under www.vita-zahnfabrik.com

Marking

VITA Zahnfabrik has been certified, and the following product bears the marks **VITA VMK Master®**.

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CH REP VITA Zahnfabrik H. Rauter GmbH & Co.KG, Bad Säckingen (Germany)
Zweigniederlassung Basel c/o Perrig AG, Max Kämpf-Platz 1, 4058 Basel
VITA Zahnfabrik · H. Rauter GmbH & Co. KG