# VITA YZ<sup>®</sup> SOLUTIONS

## CAD/CAM-Material

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# VITA YZ® SOLUTIONS

## CAD/CAM material

Please read the instructions for use before using the product.

VITA YZ SOLUTIONS are zirconia blanks with matched system components for reliable shade reproduction. They can be used for the production of fully/partially veneered reconstructions and monolithic bridge restorations in the anterior and posterior tooth regions.

### Intended purpose

VITA YZ SOLUTIONS products are ceramic materials for dental treatments.

### Patient target group

No restrictions.

### Intended user

Dental professionals only: dentists and dental technicians (Rx only).

### Indications for use

VITA YZ T and VITA YZ HT are approved for:

- fully-anatomical crowns and up to 14-unit bridges\* in the anterior and posterior tooth regions,
- fully and partially veneered single-tooth and up to 14-unit bridge substructures\* in the anterior and posterior tooth regions,
- single tooth restorations and up to 14-unit bridges\* on directly screwed implant abutments in the anterior and posterior tooth regions,
- primary telescopes,
- inlays\*\*\*, onlays\*\*\*, veneers\*\*\*, partial crowns\*\*\*, occlusal veneers (table tops)\*\*\*.

VITA YZ ST is approved for:

- fully-anatomical crowns and up to 14-unit\*\* bridges\* in the anterior and posterior tooth regions,
- fully and partially veneered single-tooth and up to 14-unit\*\* bridge substructures\* in the anterior and posterior tooth regions,
- single tooth restorations and up to 14-unit\*\* bridges\* on directly screwed implant abutments in the anterior and posterior tooth regions,
- inlays\*\*\*, onlays\*\*\*, veneers\*\*\*, partial crowns\*\*\*, occlusal veneers (table tops)\*\*\*.

VITA YZ XT is approved for:

- fully anatomical single tooth crowns and up to three-unit bridges,
- fully and partially veneered single-tooth crowns and up to three-unit bridge substructures in the anterior and posterior tooth regions,
- inlays\*\*\*, onlays\*\*\*, veneers\*\*\*, partial crowns\*\*\*, occlusal veneers (table tops)\*\*\*.

VITA YZ MULTI TRANSLUCENT is approved for:

- fully-anatomical crowns and up to 14-unit bridges\* in the anterior and posterior tooth regions,
- fully and partially veneered single-tooth and up to 14-unit bridge substructures\* in the anterior and posterior tooth regions,
- single tooth restorations and up to 14-unit bridges\* on directly screwed implant abutments in the anterior and posterior tooth regions,
- inlays\*\*\*, onlays\*\*\*, veneers\*\*\*, partial crowns\*\*\*, occlusal veneers (table tops)\*\*\*.

\* bridges and bridge substructures with a maximum of two adjoining bridge pontics.

\*\* VITA YZ ST is limited in Canada to bridge indications with a max. of six units with a max. of two adjoining bridge pontics.

\*\*\* for adhesive bonding only.

### Contraindication

Contraindications for VITA YZ T, VITA YZ HT, VITA YZ ST, VITA YZ XT and VITA YZ MULTI TRANSLUCENT:

- for more than two adjoining bridge pontics
- for two or more cantilever bridge pontics
- for parafunctional activity for veneered restorations, in particular for teeth grinders and clenchers
- for insufficient preparation results
- for insufficient hard tooth substance
- in patients with allergies or sensitivities to the ingredients
- temporary seating of veneered restorations
- conventional or self-adhesive incorporation of inlays, onlays, veneers, partial crowns and occlusal veneers (table top)

Also with VITA YZ XT:

- for bridge restorations with more than three units
- for cantilever bridges
- temporary seating

### Product application/control

#### Workflow:

1. Preparation, shade determination
2. Intraoral scan or taking the impression / fabricating the model
3. CAD-Process (design)
4. CAM process (milling/grinding)
5. Reworking (removal/grinding of the restoration)
6. Manual shading technique for non-colored restorations: immersion technique (optional)
7. Sintering
8. Individualization/characterization/glazing/polishing
9. Bonding

### Design parameters:

Minimum wall thickness in mm					
All information refers to sintered restorations		VITA YZ T/HT	VITA YZ ST**	VITA YZ XT	VITA YZ MULTI TRANSLUCENT
Inlay, onlay, veneer	incisal occlusal circumferential	0.5 0.5 0.4	0.6 0.6 0.5	0.8 0.8 0.7	0.6 0.6 0.5
Anterior and posterior crowns (fully anatomical or substructure)	incisal occlusal circumferential	0.5 0.5 0.4	0.6 0.6 0.5	0.8 0.8 0.7	0.6 0.6 0.5
Fully anatomical anterior bridges and substructures with one pontic	incisal circumferential	0.5 0.5	0.6 0.6	1.0 0.8	0.6 0.6
Fully anatomical posterior bridges and substructures with one pontic	occlusal circumferential	0.6 0.5	0.7 0.6	1.2 1.0	0.7 0.6
Fully anatomical multi-unit anterior bridges and substructures with two pontics	incisal circumferential	0.6 0.5	0.8 0.6	-	0.8 0.6
Fully anatomical multi-unit posterior bridges and substructures with two pontics	occlusal circumferential	0.7 0.6	0.8 0.6	-	0.8 0.6
Cantilever bridges with an extension	incisal occlusal circumferential	0.7 0.7 0.5	0.8 0.8 0.6	-	0.8 0.8 0.6

### Connector cross-sections in mm<sup>2</sup>

All information refers to sintered restorations	VITA YZ T/HT	VITA YZ ST**	VITA YZ XT	VITA YZ MULTI TRANSLUCENT
Inlay, onlay, veneer	-	-	-	-
Anterior and posterior crowns (fully anatomical or substructure)	-	-	-	-
Fully anatomical anterior bridges and substructures with one pontic	7	9	9	9
Fully anatomical posterior bridges and substructures with one pontic	9	12	12	12
Fully anatomical multi-unit anterior bridges and substructures with two pontics	9	12	-	12
Fully anatomical multi-unit posterior bridges and substructures with two pontics	12	15	-	15
Cantilever bridges with an extension	12	15	-	15

The minimum wall thicknesses refer to fully sintered restorations.

The cantilever bridge unit should be modeled approximately one third narrower in its mesiodistal dimension.

### Sintering parameters

Programs	%	T0 °C	↑ min.	↑ °C/ min.	T1 °C	→ min.	↑ °C/ min.	T2 °C	→ min.	↓ °C	%
YZ One for all	100	25	51:28	17	900	137:30	4	1450	120:00	200	100

For each VITA YZ material, you will also find corresponding (speed) sintering programs in your firing device.

The sintering results are influenced by various parameters such as model, performance and age of the sintering unit used. Detailed information on processing and sintering procedure can be found in the document Instructions for Use I full version.

**Recommended product combinations**

For recommended product combinations with the VITA system components and corresponding processing recommendations and restrictions regarding the combination possibilities of the VITA YZ SOLUTIONS products, please refer to the Instructions for Use | full version.

**Storage/disposal**

Store the VITA YZ discs in the original packaging and in a dry place. The products labelled with a pictogram for hazardous substances are to be disposed of as hazardous waste. Recyclable waste (such as attachments, paper and plastics) must be disposed of using appropriate recycling systems. If necessary, contaminated product residues should be pretreated in accordance with regional regulations and disposed of separately.

**Technical data**

Property	Unit	VITA YZ T	VITA YZ HT	VITA YZ ST	VITA YZ XT	VITA YZ MULTI TRANSLUCENT
Coefficient of thermal expansion (20-500 °C)	$10^{-6} \cdot K^{-1}$	approx. 10.5	approx. 10.5	approx. 10.3	approx. 10.0	approx. 10.2
Chemical solubility (ISO 6872)	$\mu g/cm^2$	< 20	< 20	< 20	< 20	< 20
Density after sinter firing	$g/cm^3$	approx. 6.05	approx. 6.08	approx. 6.05	approx. 6.03	approx. 6.05
3-point flexural strength (ISO 6872)	MPa	1350	1350	1200	850	approx. 850 MPa (incisal) - approx. 1200 MPa (cervical)

**Chemical composition**

Component	VITA YZ T/ HT	VITA YZ ST	VITA YZ XT	VITA YZ MULTI TRANSLUCENT
ZrO <sub>2</sub>	90 - 95	88 - 93	86 - 91	86 - 93
Y <sub>2</sub> O <sub>3</sub>	4 - 6	6 - 8	8 - 10	6 - 10
HfO <sub>2</sub>	1 - 3	1 - 3	1 - 3	1 - 3
Al <sub>2</sub> O <sub>3</sub>	0 - 1	0 - 1	0 - 1	0 - 1
Pigments	0 - 1	0 - 1	0 - 1	0 - 2

**Safety note**

- While work is in progress, wear suitable safety goggles/face protection
- Only perform work under an extraction unit
- Wear protective gloves

Information on bonding of VITA YZ restorations in the patient's mouth can be found at [www.vita-zahnfabrik.com/adiva](http://www.vita-zahnfabrik.com/adiva).

**Product reliability**

Information on reporting serious incidents in connection with medical devices, general risks associated with dental treatments, residual risks and (if applicable) short clinical safety and performance reports (SSCPs) can be found at [www.vita-zahnfabrik.com/product\\_safety](http://www.vita-zahnfabrik.com/product_safety).

The safety data sheets can be downloaded at [www.vita-zahnfabrik.com](http://www.vita-zahnfabrik.com) or requested by fax at (+49) 7761-562-233.

**Disclaimer**

Please note: Our products must be used in accordance with the instructions for use. We accept no liability for any damage resulting from incorrect handling or usage. The user is furthermore obliged to check the product before use with regard to its suitability for the intended area of application. We cannot accept any liability if the product is used in conjunction with materials and equipment from other manufacturers that are not compatible or not authorized for use with our product and this results in damage. The VITA Modulbox is not necessarily a component of the product. Date of issue of these instructions for use: 2025-01

After the publication of this information for use, any previous versions become obsolete. The current version can be found at [www.vita-zahnfabrik.com](http://www.vita-zahnfabrik.com).

VITA Zahnfabrik has been certified and the following products bear the mark:

**VITA YZ<sup>®</sup> T, VITA YZ<sup>®</sup> HT, VITA YZ<sup>®</sup> ST, VITA YZ<sup>®</sup> XT, VITA YZ<sup>®</sup> MULTI TRANSLUCENT**

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CH REP VITA Zahnfabrik H. Rauter GmbH & Co.KG, Bad Säckingen (Germany)
Zweigniederlassung Basel c/o Perrig AG, Max Kämpf-Platz 1, 4058 Basel
VITA Zahnfabrik · H. Rauter GmbH & Co. KG