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Instructions for use EN

VITA VIONIC® BASE DISC HI Denture Base Disc

Please read the instructions for use before using the product.

VITA VIONIC BASE DISC HI is a millable blank made of a high-quality, industrially polymerized acrylic polymer that features shade stability for the CAD/CAM fabrication of removable dentures.

Indication

VITA VIONIC BASE DISC HI are pre-colored dental milling blanks made of impact-resistant PMMA for the manufacture of denture bases for removable dentures. The product is suitable for long-term use in the oral cavity for up to 10 years.

Contraindication

- Intolerance to the contained components
- Crowns, bridges, model casting with clasps, bars, narrow sublin-gual bars and transversal connectors

Product application/control

Please adhere to the following workflow when using VITA VIONIC BASE DISC HI:

Designing the restoration:

- Position the denture base in the center of the milling blank.
- Do not go below a minimum wall thickness of 1 mm.

Milling and further processing:

- Clamp the PMMA disc in the holder system of the milling machine suitable for PMMA.
- Mill the final denture base. Milling must be carried out with a milling strategy adapted to PMMA. Only use one-edged milling tools in the milling machine! Wet milling is not recommended due to the alcoholic solvents contained in the coolants.
- Separate the denture base from the disc and grind the excess (protrusions) flush with the base. Carbide cutters with FSQ toothing and black-ring cutters for titanium are suitable for roughing. The separation of the denture bases from the milling blank and the fine finishing can be performed with simple, cross-toothed carbide millers. The connectors are cut out carefully and without high pressure.
- Carefully sandblast the milled cavities in the denture bases with Al₂O₃ (50 to 125 µm, 1-2 bar) and remove the mill chippings with compressed air (with water separator). Observe the operating instructions for the blasting unit used.
- The polishing is carried out in the usual manner with polishing agents customary for denture base materials. When polishing the material, ensure that the contact pressure is kept as low as possible, otherwise the heat generated can lead to smearing / melting.

Relining and repair:

- Denture bases made of VITA VIONIC BASE DISC HI can be relined and repaired with a commercially available PMMA cold-curing resin.

Cleaning:

- They can be cleaned in an ultrasound bath (max. 40°C) with mild cleaners.
Do **not** use:
 - Steam jet air ejectors
 - Autoclaves (e.g., sterilization and disinfection)
 - Organic solvents and acids (e.g., acetone, hydrochloric acid)
 - Highly concentrated rubbing alcohols
 - Chemical denture cleaners
 - Abrasive cleaning agents
 - Denture finishing spray

Detailed information on working with this product can be found in the document “Workflow VITA VIONIC SOLUTIONS“ at <https://www.vita-zahnfabrik.com/VionicBaseHi>.

Recommended product combinations

The following product combinations with the following VITA system components are recommended: VITA VIONIC VIGO®, VITA VIONIC® DENT DISC multiColor, VITA VIONIC® BOND, VITA AKZENT® LC, VITAVM=LC flow. For an updated overview of compatible materials and systems, see www.vita-zahnfabrik.com/vionic_compatibility.

Storage/disposal

Carefully read the separate instructions for use before starting work with the polymer. The medical device may only be processed by trained personnel. Only use clean, defect-free material from undamaged packaging and ensure that the blanks are stored at room temperature, in their original cartons, in a dry location and protected from light. Avoid all contact with fluids.

The rest of the blank material can be placed in the recycling bin.

Technical data

Chemical composition:

Polymethyl methacrylate (PMMA) CAS no. 9011-14-7	> 99 %
Color pigments	< 1 %

Physical properties:

Flexural strength [ISO20795-1]	65 MPa
Modulus of elasticity [ISO20795-1]	2030 MPa
Fracture toughness [ISO20795-1]	>2.6 MPam½
Fracture work [ISO20795-1]	>2200 J/m²
Charpy impact strength +23°C [ISO179/1eU]	45 kJ/m²
Density [ISO 1183]	1.17 g/cm³
Water absorption [ISO20795-1]	< 24 µg/mm³
Solubility in water [ISO20795-1]	<0.3 µg/mm³
Residual monomer content [ISO20795-1]	< 0.5 %
Service life	10 years

The technical data are averages and may fluctuate slightly between delivered batches.

Safety information

Please take note of the information in the latest version of the safety data sheet.

Avoid the inhalation of milling dust during processing. In case of inhalation: Move subject to fresh air. Give oxygen or artificial respiration as required. Avoid all contact with fluids. Wear gloves, protective goggles and a face mask to avoid skin irritation. In case of eye contact: Flush eyes promptly with copious amounts of water for 15 minutes and consult a physician. Wash skin with soap and water. In case of ingestion: Contact your regional Poison Control Center immediately.

The safety data sheets can be downloaded at www.vita-zahnfabrik.com/VionicBaseHi or requested by fax at +49 7761 562-233.

Any serious incident that has occurred in relation to the devices should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Disclaimer

Please note: Our products must be used in accordance with the instructions for use. We do not accept liability for any damage re-sulting from incorrect handling or usage. The user is furthermore obliged to check the product before use with regard to its suitability for the intended area of applications. We cannot accept any liability if the product is used in conjunction with materials and equipment from other manufacturers that are not compatible or not authorized for use with our product and this results in damage. Date of issue of these instructions for use: 2025-03

After the publication of this information for use, any previous versions become obsolete. The current version can be found at www.vita-zahnfabrik.com.

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VITA Zahnfabrik has been certified, and the following products bear the mark:

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