VITA

VITA VIONIC® TRY-IN RESIN

DE	Gebrauchsanweisung Kurzversion.....	4
EN	Instructions for use Short version see page 119 for US version.....	8
FR	Mode d'emploi Version courte.....	12
ES	Instrucciones de uso Versión abreviada.....	16
IT	Istruzioni d'uso Versione breve.....	20
BG	Инструкция за употреба Съкратена версия.....	24
CS	Návod k použití Stručná verze.....	28
DA	Brugsanvisning kort version.....	32
EL	Οδηγίες χρήσης Συνοπτικότερη έκδοση.....	36
ET	Kasutusjuhend lühiversioon.....	40
FI	Käyttöohje Lyhyt versio.....	44
HR	Upute za uporabu kratka verzija.....	48
HU	Használati utasítás Rövid változat.....	52
JA	使用説明書 ショートバージョン.....	56
LT	Naudojimo instrukcija Trumpa versija.....	60

LV	Lietošanas instrukcija Īsā versija	64
NO	Bruksanvisning Kortversjon	68
NL	Gebruiksaanwijzing Korte versie	72
PL	Instrukcja używania Skrócona wersja	76
pt-BR	Instruções de uso Versão resumida	80
pt-PT	Instruções de utilização Versão resumida	84
RO	Instrucțiuni de utilizare Versiunea pe scurt	88
RU	Инструкция по применению Краткая версия	92
SK	Návod na použitie skrátená verzia	96
SL	Navodila za uporabo Kratka različica	100
SV	Bruksanvisning Kort version	104
TR	Kullanım kılavuzu Kısa versiyon	108
UA	Інструкція із застосування Коротка версія	112
ZH	使用說明書 簡短版	116
EN-US	Instructions for use - US Version Short version	119

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VITA VIONIC TRY-IN RESIN is a light-cured resin for the fabrication of try-in dentures using the 3D printing process.

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Contraindication

Some ingredients of VITA VIONIC TRY-IN RESIN may cause allergic reactions in predisposed persons. In such cases refrain from using the product. Insert the product only intraorally in completely polymerised state.

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General guidelines

- Follow the instructions of the hardware manufacturer with regard to parameter settings / printing and post-curing recommendations.
- A layer thickness of 100 µm is suitable for the printing process.
- Ideal processing temperature 23°C ± 2°C / 70-77°F.
- To avoid adverse effects on the material quality, do not expose the liquid material to light or UV radiation under any circumstances.
- Deviating from the described fabrication process or storage conditions could result in deviations in the mechanical and optical properties of the material.
- Always keep the container tightly closed. Carefully close the container after each use.

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An up-to-date list of **validated 3D printers, washing units and post-curing units**, as well as suitable printer templates and post-curing parameters, can be found at www.vita-zahnfabrik.com/VionicTryInResin.

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Please read the instructions for use before using the product.

VITA VIONIC TRY-IN RESIN is a light-cured resin for the fabrication of try-in dentures to control occlusion, bite registration and esthetics. It can be used in combination with a stereolithographic 3D printer using a 385 nm light source. The product processing requires a computer aided and manufacturing system (CAD/CAM) that includes the following: scanner, design software, additive printer and post-cure unit. Those devices and software is not part of the product.

Intended Use

Resin for dental 3D printing.

Patient Target Group

Persons being treated in the context of a dental procedure.

Intended Users:

Dentist, dental technician (Rx only).

Indication

Individual functional try-ins of complete dentures.

Contraindication

Contains(meth)acrylates and phosphine oxides. Some ingredients of VITA VIONIC TRY-IN RESIN may cause allergic reactions in predisposed persons. In such cases refrain from using the product. Insert the product only intraorally in completely polymerised state.

Product application

Please adhere to the following workflow when using this product:

Design process:

1. Prepare the data:
 - CAD design.
 - Generation of an STL file
 - Convert the data to the CAM software and create supporting structures. Edit the supporting structures suggested by the software. Avoid supporting structures in cavities.
 - Conversion to a printer-compatible format.

2. Select the process and material parameters in the printer software: Use the published printing template for VITA VIONIC TRY-IN RESIN from your printer manufacturer.
3. Transfer the prepared data to the 3D printer.
4. Preparing for 3D printing: Before use, the material in the bottle must be shaken well and homogenized with a bottle roller to avoid any color deviations or misprints.
5. Fill the resin tank of the 3D printer in accordance with the instructions of the printer manufacturer.
6. Start the printing process.

Reworking process:

7. Allow the printed objects to drain for approx. 10 min. in the printer. Reworking should take place as soon as possible after the design process.
8. Remove the platform from the printer and remove the printed objects from the platform. In the case of multiple objects, separate from each other.
9. Remove the outer supporting structure and remove excess resin with compressed air.

Cleaning:

10. Clean the finished printed objects twice for three minutes each time in an ultrasonic bath with fresh isopropanol $\geq 98\%$, according to the device manufacturer's instructions.
11. Dry the cleaned printed objects until there are no longer any traces of isopropyl alcohol. Check openings, hollow spaces and gaps for residue and blow out with compressed air to speed up the drying process. The surface may no longer be sticky and shiny after cleaning.

Post-exposure:

12. Perform post-curing with a suitable light-curing device according to the manufacturer's instructions. The end product's properties and shade effects depend on several factors, such as the reworking process. Correct post-exposure is important for biocompatibility. For this reason, it must be ensured that the exposure device is in proper condition and that the moulds are completely hardened.

Surface treatment:

13. The denture base is now ready for subsequent processing or finishing. Process the surface mechanically (smooth with a milling machine) and polish. Pre-polishing is done using brushes and pre-polishing paste/pumice stone, and high-gloss polishing is done using buffing wheels and high-gloss polishing agents for plastics.

Final cleaning:

14. The cleaning process is carried out by dental laboratories/dental practices in an ultrasonic bath. Do not use heat-based disinfection or sterilization methods. This could cause the denture to lose its shape.

General guidelines

- Follow the instructions of the hardware manufacturer with regard to parameter settings / printing and post-curing recommendations.
- A layer thickness of 100 µm is suitable for the printing process.
- Ideal processing temperature 23°C ±2°C / 70-77°F.
- To avoid adverse effects on the material quality, do not expose the liquid material to light or UV radiation under any circumstances.
- Deviating from the described fabrication process or storage conditions could result in deviations in the mechanical and optical properties of the material.
- Always keep the container tightly closed. Carefully close the container after each use.
- The 3D printing resin can remain in the vat for a maximum of two days (maximum 10 print cycles).
- It is not allowed to pour 3D printing resin from the reservoir of the printer back into the bottle.

Recommended product combinations

For recommended product combinations with the VITA system components, corresponding processing recommendations and restrictions regarding the combination possibilities of the product, please refer to the Instructions for Use | full version.

Use of this product is only authorized in conjunction with compatible devices.

The devices should be set up, validated and maintained in accordance with their labeling and instructions for use. Any unauthorized changes to the process equipment, parameters, or software may result in a product that is out of specification.

An up-to-date list of **validated 3D printers**, **washing units** and **post-curing units**, as well as suitable printer templates and post-curing parameters, can be found at www.vita-zahnfabrik.com/VionicTryInResin.

Storage/disposal

Storage at temperature, 15°C (59°F) to 28°C (82°F). Protect from direct sunlight. Even a small amount of light can trigger polymerization. To protect against contamination, cover the material in the tray with a lid or a glass plate. The products labelled with a pictogram for hazardous substances are to be disposed of as hazardous waste. Recyclable waste (such as attachments, paper and plastics) must be disposed of using appropriate recycling systems. If necessary, contaminated product residues should be pretreated in accordance with regional regulations and disposed of separately.

Adverse effects

Product may cause allergic reactions.

Safety information

- Avoid direct contact with the liquid material and the components before post-curing, especially in pregnant/breastfeeding women. Irritating to eyes and skin (sensitization possible).
- When handling uncovered, use equipment with local exhaust ventilation.
- Wear suitable personal protective equipment (when handling the uncured material: protective gloves and goggles, when finishing the cured material additionally face mask).
- After contact with skin wash immediately with water and soap.
- If resin gets into your eyes immediately rinse and consult a doctor.
- Biocompatibility is only guaranteed with complete polymerisation.
- Please take note of the information in the latest version of the safety data sheet.
- Any serious incident that has occurred in relation to this product must be reported immediately to info@vita-zahnfabrik.com and to the competent authority of the Member State in which the user and/or patient is established.

MRI Safety Information: VITA VIONIC TRY-IN RESIN is MR Safe as it is composed of materials that are electrically non conductive, nonmetallic and nonmagnetic.

Technical Data

Chemical composition	Wt%
Polyfunktional Acrylates	46-65%
Urethane bismethacrylate polymer	30-35 %
Methacrylate monomers	2-10 %
Initiators, pigments, stabilizers	2-6 %
MMA and phthalat free	0%


Physical properties*	Value
Flexural strength <small>DIN EN ISO 20795-1*</small>	> 100 MPa
Flexural modulus <small>DIN EN ISO 20795-1*</small>	> 2200 MPa
Hardness	> 75 ShoreD

*Applies to cured objects | according to norm at room temperature

Disclaimer

Please note: Our products must be used in accordance with the instructions for use. We do not accept liability for any damage resulting from incorrect handling or usage. We do not accept liability for any damage resulting from improper use of the impression material. The user is furthermore obliged to check the product before use with regard to its suitability for the intended area of applications. We cannot accept any liability if the product is used in conjunction with materials and equipment from other manufacturers that are not compatible or not authorized for use with our product and this results in damage. Date of issue of these instructions for use: 2025-01
After the publication of this information for use, any previous versions become obsolete. The current version can be found at www.vita-zahnfabrik.com.

Symbol Explanation

	Manufacturer
	Use-by date
	Batch code
	Distributor

	Keep away from sunlight
	Temperature limit
	Consult instructions for use
	Medical Device

	Prescription only
	Magnetic resonance
	Caution
	This way up

VITA Zahnfabrik has been certified, and the product **VITA VIONIC® TRY-IN RESIN** bears the mark:

CE 0483

MD ⓘ R Only

Made in Germany

VITA Zahnfabrik H. Rauter GmbH & Co. KG

Spitalgasse 3
79713 Bad Säckingen
Germany

Phone: +49 7761 562-0
Fax: +49 7761 562-299
E-Mail: info@vita-zahnfabrik.com

detax GmbH

Carl-Zeiss-Str. 4
76275 Ettlingen
Germany

Phone: +49 7243 510-0
Fax: +49 7243 510-100
E-Mail: post@detax.com

922-13051_0125_V01

VITA