# VITA

# VITA VIONIC® BASE RESIN IMPACT

DE	Gebrauchsanweisung   Kurzversion.....	4
EN	Instructions for use   Short version   see page 119 for US version.....	8
FR	Mode d'emploi   Version courte.....	12
ES	Instrucciones de uso   Versión abreviada.....	16
IT	Istruzioni d'uso   Versione breve.....	20
BG	Инструкция за употреба   Съкратена версия.....	24
CS	Návod k použití   Stručná verze.....	28
DA	Brugsanvisning   kort version.....	32
EL	Οδηγίες χρήσης   Συνομημεμένη έκδοση.....	36
ET	Kasutusjuhend   lühiversioon.....	40
FI	Käyttöohje   Lyhyt versio.....	44
HR	Upute za uporabu   kratka verzija.....	48
HU	Használati utasítás   Rövid változat.....	52
JA	使用説明書   ショートバージョン.....	56
LT	Naudojimo instrukcija   Trumpa versija.....	60

LV	Lietošanas instrukcija   Īsā versija .....	64
NO	Bruksanvisning   Kortversjon .....	68
NL	Gebruiksaanwijzing   Korte versie .....	72
PL	Instrukcja używania   Skrócona wersja .....	76
pt-BR	Instruções de uso   Versão resumida .....	80
pt-PT	Instruções de utilização   Versão resumida .....	84
RO	Instrucțiuni de utilizare   Versiunea pe scurt .....	88
RU	Инструкция по применению   Краткая версия .....	92
SK	Návod na použitie   skrátená verzia .....	96
SL	Navodila za uporabo   Kratka različica .....	100
SV	Bruksanvisning   Kort version .....	104
TR	Kullanım kılavuzu   Kısa versiyon .....	108
UA	Інструкція із застосування   Коротка версія .....	112
ZH	使用說明書   簡短版 .....	116
EN-US	Instructions for use   Short version .....	119

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DE

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VITA VIONIC BASE RESIN IMPACT is a light-cured resin for the fabrication of removable denture bases using the 3D printing process.

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Final cleaning:

15. The cleaning process is carried out by dental laboratories/dental practices in an ultrasonic bath. Do not use heat-based disinfection or sterilization methods. This could cause the denture to lose its shape.
16. The patient should use lukewarm water, a soft toothbrush, liquid soap or detergent to clean the dentures daily. Rinse well afterwards. Denture cleaning tablets can be used according to the manufacturer's instructions.  
Note: Do not use toothpaste; the abrasive particles are coarse and scratch the denture.  
Do not use any acids or solvents.

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**Recommended product combinations**

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### Disclaimer

Please note: Our products must be used in accordance with the instructions for use. We do not accept liability for any damage resulting from incorrect handling or usage. We do not accept liability for any damage resulting from improper use of the impression material. The user is furthermore obliged to check the product before use with regard to its suitability for the intended area of applications. We cannot accept any liability if the product is used in conjunction with materials and equipment from other manufacturers that are not compatible or not authorized for use with our product and this results in damage. Date of issue of these instructions for use: 2025-01. After the publication of this information for use, any previous versions become obsolete. The current version can be found at [www.vita-zahnfabrik.com](http://www.vita-zahnfabrik.com).

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Please read the instructions for use before using the product.

VITA VIONIC BASE RESIN IMPACT is a light-cured resin for the fabrication of removable denture bases. It can be used in combination with a stereolithographic 3D printer using a 385 nm light source. The product processing requires a computer aided and manufacturing system (CAD/CAM) that includes the following: scanner, design software, additive printer and post-cure unit. Those devices and software is not part of the product.

CAUTION STATEMENT: Caution: Federal Law restricts this device to sale by or on the order of a dentist.

### **Intended Use**

Resin for dental 3D printing.

### **Patient Target Group**

Persons being treated in the context of a dental procedure.

### **Intended Users**

Dentist, dental technician (Rx only).

### **Indication**

Denture bases.

### **Contraindication**

Contains (meth)acrylate and phosphine oxide.

Some ingredients of VITA VIONIC BASE RESIN IMPACT may cause allergic reactions in predisposed persons. In such cases refrain from using the product. Insert the product only intraorally in completely polymerised state.

### **Product application**

Please adhere to the following workflow when using this product:

#### Design process:

1. Prepare the data:

- CAD design.

Recommended **bonding gap settings:**

- VITA VIONIC BASE RESIN IMPACT in combination with VITA VIONIC DENT RESIN or VITA VIONIC DENT DISC multiColor: **150µm**
- VITA VIONIC BASE RESIN IMPACT in combination with VITA VIONIC VIGO: **100µm**

- Generate an STL file.
  - Convert the data to the CAM software and create supporting structures. Edit the supporting structures suggested by the software. Avoid supporting structures in cavities.
  - Convert to a printer-compatible format.
2. Select the process and material parameters in the printer software: Use the published printing template for VITA VIONIC BASE RESIN IMPACT from your printer manufacturer.
  3. Transfer the prepared data to the 3D printer.
  4. Preparing for 3D printing: Before use, the material in the bottle must be shaken well and homogenized with a bottle roller to avoid any color deviations or misprints.
  5. Fill the resin tank of the 3D printer in accordance with the instructions of the printer manufacturer.
  6. Start the printing process.

#### Reworking process:

7. Allow the printed objects to drain for approx. 10 min. in the printer. Reworking should take place as soon as possible after the design process.
8. Remove the platform from the printer and remove the printed objects from the platform. In the case of multiple objects, separate from each other.
9. Remove the outer supporting structure.
10. Remove excess resin with compressed air.

#### Cleaning:

11. Clean the finished printed objects twice for three minutes each time in an ultrasonic bath with fresh isopropanol  $\geq 98\%$ , according to the device manufacturer's instructions.
12. Dry the cleaned printed objects until there are no longer any traces of isopropyl alcohol. Check openings, hollow spaces and gaps for residue and blow out with compressed air to speed up the drying process. The surface may no longer be sticky and shiny after cleaning.

#### Post-exposure:

13. Perform post-curing with a suitable light-curing device according to the manufacturer's instructions. The end product's properties and shade effects depend on several factors, such as the reworking process. Correct post-exposure is important for biocompatibility. For this reason, it must be ensured that the exposure device is in proper condition and that the moulds are completely hardened.

#### Surface treatment:

14. The denture base is now ready for subsequent processing or finishing. Process the surface mechanically (smooth with a milling machine) and polish. Pre-polishing is done using brushes and pre-polishing paste/pumice stone, and high-gloss polishing is done using buffing wheels and high-gloss polishing agents for plastics.

#### Final cleaning:

15. The cleaning process is carried out by dental laboratories/dental practices in an ultrasonic bath. Do not use heat-based disinfection or sterilization methods. This could cause the denture to lose its shape.

16. The patient should use lukewarm water, a soft toothbrush, liquid soap or detergent to clean the dentures daily. Rinse well afterwards. Denture cleaning tablets can be used according to the manufacturer's instructions.

Note: Do not use toothpaste; the abrasive particles are coarse and scratch the denture.  
Do not use any acids or solvents.

### General guidelines

- Follow the instructions of the hardware manufacturer with regard to parameter settings, printing and post-curing recommendations.
- A layer thickness of 100 µm is suitable for the printing process.
- Ideal processing temperature 23°C ±2°C / 70-77°F.
- To avoid adverse effects on the material quality, do not expose the liquid material to light or UV radiation under any circumstances.
- Deviating from the described fabrication process or storage conditions could result in deviations in the mechanical and optical properties of the material.
- Always keep the container tightly closed. Carefully close the container after each use.
- It is not allowed to pour 3D printing resin from the reservoir of the printer back into the bottle.

### Recommended product combinations

For recommended product combinations with the VITA system components, corresponding processing recommendations and restrictions regarding the combination possibilities of the product, please refer to the Instructions for Use | full version.

Use of this product is only authorized in conjunction with compatible devices.

The devices should be set up, validated and maintained in accordance with their labeling and instructions for use. An up-to-date list of **validated 3D printers**, **washing units** and **post-curing units**, as well as suitable printer templates and post-curing parameters, can be found at [www.vita-zahnfabrik.com/VionicBaseResinImpact](http://www.vita-zahnfabrik.com/VionicBaseResinImpact).

### Storage/disposal

Storage at temperature, 15°C (59°F) to 28°C (82°F). Protect from direct sunlight. Even a small amount of light can trigger polymerization. To protect against contamination, cover the material in the tray with a lid or a glass plate. The products labelled with a pictogram for hazardous substances are to be disposed of as hazardous waste. Recyclable waste (such as attachments, paper and plastics) must be disposed of using appropriate recycling systems. If necessary, contaminated product residues should be pretreated in accordance with regional/national and international regulations and disposed of separately.

### Adverse effects

Product may cause allergic reactions.

## Safety information

- Avoid direct contact with the liquid material and the components before post-curing, especially in pregnant/breastfeeding women. Irritating to eyes and skin (sensitization possible).
- When handling uncovered, use equipment with local exhaust ventilation.
- Wear suitable personal protective equipment (when handling the uncured material: protective gloves and goggles, when finishing the cured material additionally face mask).
- After contact with skin wash immediately with water and soap.
- If resin gets into your eyes immediately rinse and consult a doctor.
- Biocompatibility is only guaranteed with complete polymerisation.
- Please take note of the information in the latest version of the safety data sheet.
- Any serious incident that has occurred in relation to this product must be reported immediately to [info@vita-zahnfabrik.com](mailto:info@vita-zahnfabrik.com) and to the competent authority of the Member State in which the user and/or patient is established.

MRI Safety Information: VITA VIONIC BASE RESIN IMPACT is MR Safe as it is composed of materials that are electrically non conductive, nonmetallic and nonmagnetic.

## Technical Data

Chemical composition	Wt%
Polyfunktional Acrylates	40-60%
Urethane bismethacrylate polymer	20-40 %
Methacrylate monomers	5-20 %
Fillers, pigments, stabilizers	< 5%
Photoinitiators	< 0.2-10%
MMA and phthalat free	0%

Physical properties*	Value
Flexural strength <small>DIN EN ISO 20795-1*</small>	> 65 MPa
Flexural modulus <small>DIN EN ISO 20795-1*</small>	> 2000 MPa
Fracture work <small>DIN EN ISO 179-1 unnotched</small>	≥ 25 kJ/m <sup>2</sup>
Water sorption <small>DIN EN ISO 20795-1*</small>	< 32 µg/mm <sup>3</sup>
Water solubility <small>DIN EN ISO 20795-1*</small>	< 1.6 µg/mm <sup>3</sup>
Viscosity**	400-1000 mPas


\*Applies to cured objects | according to norm at room temperature \*\*Applies to liquid resin

## Disclaimer

Please note: Our products must be used in accordance with the instructions for use. We do not accept liability for any damage resulting from incorrect handling or usage. We do not accept liability for any damage resulting from improper use of the impression material. The user is furthermore obliged to check the product before use with regard to its suitability for the intended area of applications. We cannot accept any liability if the product is used in conjunction with materials and equipment from other manufacturers that are not compatible or not authorized for use with our product and this results in damage. Date of issue of these instructions for use: 2025-01

After the publication of this information for use, any previous versions become obsolete. The current version can be found at [www.vita-zahnfabrik.com](http://www.vita-zahnfabrik.com).

## Symbol Explanation

Graphic	Title	Description	Reference
	Manufacturer	Indicates the medical device manufacturer.	ISO 15223-1:2021 Medical devices - Symbols to be used with information to be supplied by the manufacturer - Part 1: General requirements
	Use-by date	Indicates the date after which the medical device is not to be used.	ISO 15223-1:2021 Medical devices - Symbols to be used with information to be supplied by the manufacturer - Part 1: General requirements
	Batch code	Indicates the manufacturer's batch code so that the batch or lot can be identified.	ISO 15223-1:2021 Medical devices - Symbols to be used with information to be supplied by the manufacturer - Part 1: General requirements
	Medical device	Indicates that the product is a medical device.	ISO 15223-1:2021 Medical devices - Symbols to be used with information to be supplied by the manufacturer - Part 1: General requirements
	Keep away from sunlight	Indicates a medical device that needs protection from light sources.	ISO 15223-1:2021 Medical devices - Symbols to be used with information to be supplied by the manufacturer - Part 1: General requirements
	Temperature limit	Indicates the temperature limits to which the medical device can be safely exposed.	ISO 15223-1:2021 Medical devices - Symbols to be used with information to be supplied by the manufacturer - Part 1: General requirements
	Consult instructions for use	Indicates the need for the user to consult the instructions for use.	ISO 15223-1:2021 Medical devices - Symbols to be used with information to be supplied by the manufacturer - Part 1: General requirements
	Prescription only	Requires a prescription in the United States of America. Caution: In the United States of America, federal law restricts this device to sale or use by, or on the order of, a physician.	USA Code of Federal Regulations 21CFR Part 801 § 801.109 (b)(1)
	Magnetic resonance	The freeprint crown is MR Safe.	FDA Guidance document: Testing and Labeling Medical Devices for Safety in the Magnetic Resonance (MR) Environment dated May 20, 2021.
	Caution	The uncured material is considered as hazardous substance.	29 CFR 1910.1200 Occupational Safety and Health Standards
	This way up	This is the correct upright position of the distribution packages for transport and/or storage.	ISO 780:2015(E) Packaging-Distribution packaging- Graphical symbols for handling and storage of packages
	Distributor	Indicates the entity distributing the medical device into the locale	ISO 15223-1:2021 Medical devices - Symbols to be used with information to be supplied by the manufacturer - Part 1: General requirements

VITA Zahnfabrik has been certified, and the product **VITA VIONIC®** BASE RESIN IMPACT bears the mark:

CE 0483

MD  R Only

Made in Germany

## **VITA Zahnfabrik H. Rauter GmbH & Co. KG**

Spitalgasse 3  
79713 Bad Säckingen  
Germany

Phone: +49 7761 562-0  
Fax: +49 7761 562-299  
E-Mail: [info@vita-zahnfabrik.com](mailto:info@vita-zahnfabrik.com)

## **detax GmbH**

Carl-Zeiss-Str. 4  
76275 Ettlingen  
Germany

Phone: +49 7243 510-0  
Fax: +49 7243 510-100  
E-Mail: [post@detax.com](mailto:post@detax.com)

922-13052\_0125\_V01

# VITA