# VITA

# VITA VIONIC® DENT RESIN

|    |  |    |
|----|--|----|
| DE | Gebrauchsanweisung   Kurzversion.....                                    | 4  |
| EN | Instructions for use   Short version   see page 120 for US version ..... | 8  |
| FR | Mode d'emploi   Version courte.....                                      | 12 |
| ES | Instrucciones de uso   Versión abreviada.....                            | 16 |
| IT | Istruzioni d'uso   Versione breve.....                                   | 20 |
| BG | Инструкция за употреба   Съкратена версия.....                           | 24 |
| CS | Návod k použití   Stručná verze .....                                    | 28 |
| DA | Brugsanvisning   kort version .....                                      | 32 |
| EL | Οδηγίες χρήσης   Συνοπτικότερη έκδοση .....                              | 36 |
| ET | Kasutusjuhend   lühiversioon.....  | 40 |
| FI | Käyttöohje   Lyhyt versio.....   | 44 |
| HR | Upute za uporabu   kratka verzija.....                                   | 48 |
| HU | Használati utasítás   Rövid változat .....                               | 52 |
| JA | 使用説明書   ショートバージョン .....  | 56 |
| LT | Naudojimo instrukcija   Trumpa versija.....                              | 60 |

|       |   |     |
|-------|---|-----|
| LV    | Lietošanas instrukcija   Īsā versija .....            | 64  |
| NO    | Bruksanvisning   Kortversjon .....                    | 68  |
| NL    | Gebruiksaanwijzing   Korte versie .....               | 72  |
| PL    | Instrukcja używania   Skrócona wersja .....           | 76  |
| pt-BR | Instruções de uso   Versão resumida .....             | 80  |
| pt-PT | Instruções de utilização   Versão resumida .....      | 84  |
| RO    | Instrucțiuni de utilizare   Versiunea pe scurt .....  | 88  |
| RU    | Инструкция по применению   Краткая версия .....       | 92  |
| SK    | Návod na použitie   skrátená verzia .....             | 96  |
| SL    | Navodila za uporabo   Kratka različica .....          | 100 |
| SV    | Bruksanvisning   Kort version .....                   | 104 |
| TR    | Kullanım kılavuzu   Kısa versiyon .....               | 108 |
| UA    | Інструкція із застосування   Коротка версія .....     | 112 |
| ZH    | 使用說明書   簡短版 .....                                     | 116 |
| EN-US | Instructions for use US version   Short version ..... | 120 |

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VITA VIONIC DENT RESIN is a light-cured resin for the fabrication of removable denture teeth using the 3D printing process.

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Dentist, dental technician.

### Indication

Permanent denture teeth and long-term provisional crowns and bridges for a wearing time up to 12 months.

### Contraindication

Contains (meth)acrylate and phosphine oxides. Some ingredients of VITA VIONIC DENT RESIN may cause allergic reactions in predisposed persons. In such cases refrain from using the product. Only insert the product intraorally in completely polymerised state.

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#### Design process:

1. Prepare the data:
  - CAD design.  
Recommended ***bonding gap settings*** in combination with VITA VIONIC BASE RESIN IMPACT: ***150µm***
  - Generation of an STL file
  - Convert the data to the CAM software and create supporting structures. Edit the supporting structures suggested by the software.
  - Conversion to a printer-compatible format.
2. Select the process and material parameters in the printer software: Use the published printing template for VITA VIONIC DENT RESIN from your printer manufacturer.
3. Transfer the prepared data to the 3D printer.

4. Preparing for 3D printing: Before use, the material in the bottle must be shaken well and homogenized with a bottle roller to avoid any color deviations or misprints.
5. Fill the resin tank of the 3D printer in accordance with the instructions of the printer manufacturer.
6. Start the printing process.

#### Post-processing:

7. Allow the printed objects to drain for approx. 10 min. in the printer. Reworking should take place as soon as possible after the design process.
8. Remove the platform from the printer and remove the printed objects from the platform. In the case of multiple objects, separate from each other.
9. Remove the outer supporting structure.
10. Remove excess resin with compressed air.

#### Cleaning:

11. Clean the finished printed objects twice for 1 Min. minutes each time in an ultrasonic bath with fresh isopropanol  $\geq 98\%$ , according to the device manufacturer's instructions.
12. Dry the cleaned printed objects until there are no longer any traces of isopropyl alcohol. Check openings, hollow spaces and gaps for residue and blow out with compressed air to speed up the drying process. The surface may no longer be sticky and shiny after cleaning. Store components protected from light at room temperature for 30 minutes before post-exposure.

#### Bonding:

13. Leave the printed denture base and teeth in their green body state (not post-exposed). Place the printed denture teeth in the corresponding tooth receptacles of the denture base, check the correct seating of the teeth and remove them again.
14. Apply liquid VITA VIONIC BASE RESIN IMPACT to the cavities with a small brush or a syringe.
15. Connect the repositioned denture teeth by exposure in a Xenon flash unit (Otoflash G171) with 400 flashes under inert gas conditions (nitrogen). Finally post-cure the complete denture in a xenon flash unit (Otoflash G171) with  $2 \times 1800$  flashes under inert gas conditions (nitrogen).

#### Post-exposure:

16. Perform post-curing with a suitable light-curing device according to the manufacturer's instructions. The properties and shade effects of the end product depend on a number of factors, such as the reworking process. Correct post-exposure is important for biocompatibility. For this reason, it must be ensured that the exposure device is in proper condition and that the moulds are completely hardened.

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**General guidelines**

- Follow the instructions of the hardware manufacturer with regard to parameter settings / printing and post-curing recommendations.
- A layer thickness of 50 µm is suitable for the printing process.
- Ideal processing temperature 23°C ± 2°C / 70-77°F.
- To avoid adverse effects on the material quality, do not expose the liquid material to light or UV radiation under any circumstances.
- Deviating from the described fabrication process or storage conditions could result in deviations in the mechanical and optical properties of the material.
- Always keep the container tightly closed. Carefully close the container after each use.
- For long-term provisional crowns and bridges up to 3-units, the minimum occlusal thickness is 1.5 mm, circular 1.5 mm and cervical 1.0 mm. Failure to observe this information may lead to an undesirable result. The temporary crown and bridge restorations can be used with the commonly used permanent luting cements, e.g. Variolink Esthetic DC by Ivoclar or RelyX Unicem2 by 3M Espe. Observe the manufacturer's instructions.

**Recommended product combinations**

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### **Disclaimer**

Please note: Our products must be used in accordance with the instructions for use. We do not accept liability for any damage resulting from incorrect handling or usage. We do not accept liability for any damage resulting from improper use of the impression material. The user is furthermore obliged to check the product before use with regard to its suitability for the intended area of applications. We cannot accept any liability if the product is used in conjunction with materials and equipment from other manufacturers that are not compatible or not authorized for use with our product and this results in damage. Date of issue of these instructions for use: 2025-01  
After the publication of this information for use, any previous versions become obsolete. The current version can be found at [www.vita-zahnfabrik.com](http://www.vita-zahnfabrik.com).

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Please read the instructions for use before using the product.

VITA VIONIC DENT RESIN is a light-cured resin for the fabrication of removable denture teeth and temporary restorations such as crowns and bridges. It can be used in combination with a stereolithographic 3D printer using a 385nm light source. The product processing requires a computer aided and manufacturing system (CAD/CAM) that includes the following: scanner, design software, additive printer and post-cure unit. Those devices and software are not part of the product.

**CAUTION STATEMENT** Caution: Federal Law restricts this device to sale by or on the order of a dentist.

### **Intended Use**

Resin for dental 3D printing.

### **Patient Target Group**

Persons being treated in the context of a dental procedure.

### **Intended Users:**

Dentist, dental technician (Rx only).

### **Indication**

Permanent denture teeth and long-term provisional crowns and bridges.

### **Contraindication**

Contains (meth)acrylate, phosphine oxides and silanized dental glass. Overall share of inorganic fillers (particle size 1.5 µm) totals 20–40 % by mass. Some ingredients of VITA VIONIC DENT RESIN may cause allergic reactions in predisposed persons. In such cases refrain from using the product. Only insert the product intraorally in completely polymerised state.

### **Product application**

Please adhere to the following workflow when using this product:

#### Design process:

1. Prepare the data:
  - CAD design.  
Recommended ***bonding gap settings*** in combination with VITA VIONIC BASE RESIN IMPACT: **150µm**
  - Generation of an STL file
  - Convert the data to the CAM software and create supporting structures. Edit the supporting structures

suggested by the software.

- Conversion to a printer-compatible format.
- 2. Select the process and material parameters in the printer software: Use the published printing template for VITA VIONIC DENT RESIN from your printer manufacturer. Transfer the prepared data to the 3D printer.
- 3. Transfer the prepared data to the 3D printer.
- 4. Preparing for 3D printing: Before use, the material in the bottle must be shaken well and homogenized with a bottle roller to avoid any color deviations or misprints.
- 5. Fill the resin tank of the 3D printer in accordance with the instructions of the printer manufacturer.
- 6. Start the printing process.

#### Post-processing:

- 7. Allow the printed objects to drain for approx. 10 min. in the printer. Reworking should take place as soon as possible after the design process.
- 8. Remove the platform from the printer and remove the printed objects from the platform. In the case of multiple objects, separate from each other.
- 9. Remove the outer supporting structure.
- 10. Remove excess resin with compressed air.

#### Cleaning:

- 11. Clean the finished printed objects twice for 1 Min. minutes each time in an ultrasonic bath with fresh isopropanol  $\geq 98\%$ , according to the device manufacturer's instructions. Do not use heat-based disinfection or sterilization methods. This could cause the denture to lose its shape.
- 12. Dry the cleaned printed objects until there are no longer any traces of isopropyl alcohol. Check openings, hollow spaces and gaps for residue and blow out with compressed air to speed up the drying process. The surface may no longer be sticky and shiny after cleaning. Store components protected from light at room temperature for 30 minutes before post-exposure.

#### Bonding:

- 13. Leave the printed denture base and teeth in their green body state (not post-exposed). Place the printed denture teeth in the corresponding tooth receptacles of the denture base, check the correct seating of the teeth and remove them again.
- 14. Apply liquid VITA VIONIC BASE RESIN IMPACT to the cavities with a small brush or a syringe.
- 15. Connect the repositioned denture teeth by exposure in a Xenon flash unit (Otoflash G171) with 400 flashes under inert gas conditions (nitrogen). Finally post-cure the complete denture in a xenon flash unit (Otoflash G171) with  $2 \times 1800$  flashes under inert gas conditions (nitrogen).

#### Post-exposure:

- 16. Perform post-curing with a suitable light-curing device according to the manufacturer's instructions. The properties and shade effects of the end product depend on a number of factors, such as the reworking process. Correct post-exposure is important for biocompatibility. For this reason, it must be ensured that the exposure device is in proper condition and that the moulds are completely hardened.

#### Surface treatment:

17. The denture is now ready for subsequent processing or finishing. Process the surface mechanically (smooth with a milling machine) and polish. Pre-polishing is done using brushes and pre-polishing paste/pumice stone, and high-gloss polishing is done using buffing wheels and high-gloss polishing agents for plastics.

#### Final cleaning:

18. The cleaning process is carried out by dental laboratories/dental practices in an ultrasonic bath. Do not use heat-based disinfection or sterilization methods. This could cause the denture to lose its shape.
19. The patient should use lukewarm water, a soft toothbrush, liquid soap or detergent to clean the dentures daily. Rinse well afterwards. Denture cleaning tablets can be used according to the manufacturer's instructions.

Note: Do not use toothpaste; the abrasive particles are coarse and scratch the denture.  
Do not use any acids or solvents.

### **General guidelines**

- Follow the instructions of the hardware manufacturer with regard to parameter settings / printing and post-curing recommendations.
- A layer thickness of 50 µm is suitable for the printing process.
- Ideal processing temperature 23°C ± 2°C / 70-77°F.
- To avoid adverse effects on the material quality, do not expose the liquid material to light or UV radiation under any circumstances.
- Deviating from the described fabrication process or storage conditions could result in deviations in the mechanical and optical properties of the material.
- Always keep the container tightly closed. Carefully close the container after each use.
- The 3D printing resin can remain in the vat for a maximum of 14 days (maximum 10 print cycles).
- It is not allowed to pour 3D printing resin from the reservoir of the printer back into the bottle.
- For long-term provisional crowns and bridges up to 3-units, the minimum occlusal thickness is 1.5 mm, circular 1.5 mm and cervical 1.0 mm. Failure to observe this information may lead to an undesirable result. The temporary crown and bridge restorations can be used with the commonly used permanent luting cements, e.g. Variolink Esthetic DC by Ivoclar or RelyX Unicem2 by 3M Espe. Observe the manufacturer's instructions.

## Recommended product combinations

For recommended product combinations with the VITA system components, corresponding processing recommendations and restrictions regarding the combination possibilities of the product, please refer to the Instructions for Use | full version.

Use of this product is only authorized in conjunction with compatible devices. The devices should be set up, validated and maintained in accordance with their labeling and instructions for use. Any unauthorized changes to the process equipment, parameters, or software may result in a product that is out of specification.

An up-to-date list of **validated 3D-printers, washing units and post curing units**, as well as suitable printer templates and post-curing parameters, can be found at [www.vita-zahnfabrik.com/VionicDentResin](http://www.vita-zahnfabrik.com/VionicDentResin).

## Storage/disposal

Storage at temperature, 15°C (59°F) to 28°C (82°F). Protect from direct sunlight. Even a small amount of light can trigger polymerization. To protect against contamination, cover the material in the tray with a lid or a glass plate. The products labelled with a pictogram for hazardous substances are to be disposed of as hazardous waste. Recyclable waste (such as attachments, paper and plastics) must be disposed of using appropriate recycling systems. If necessary, contaminated product residues should be pretreated in accordance with regional regulations and disposed of separately.

## Adverse effects

Product may cause allergic reactions.

## Safety information

- Avoid direct contact with the liquid material and the components before post-curing, especially in pregnant/breastfeeding women. Irritating to eyes and skin (sensitization possible).
- When handling uncovered, use equipment with local exhaust ventilation.
- Wear suitable personal protective equipment (when handling the uncured material: protective gloves and goggles, when finishing the cured material additionally face mask).
- After contact with skin wash immediately with water and soap.
- If resin gets into your eyes immediately rinse and consult a doctor.
- Biocompatibility is only guaranteed with complete polymerisation.
- Please take note of the information in the latest version of the safety data sheet.
- Any serious incident that has occurred in relation to this product must be reported immediately to [info@vita-zahnfabrik.com](mailto:info@vita-zahnfabrik.com) and to the competent authority of the Member State in which the user and/or patient is established.

MRI Safety Information: VITA VIONIC DENT RESIN is MR Safe as it is composed of materials that are electrically non conductive, nonmetallic and nonmagnetic.

**Technical Data**

| Chemical composition                   | Wt%       | Physical properties*                | Value                    |
|--|-----------|-------------------------------------|--------------------------|
| Polyfunctional acrylates               | 45-80%    | Flexural strength DIN EN ISO 10477* | > 100 MPa                |
| Urethane bismethacrylate polymer       | 5-20 %    | Flexural modulus DIN EN ISO 10477*  | > 2800 MPa               |
| Methacrylate monomer                   | 0.1-5 %   | Water sorption DIN EN ISO 10477*    | < 40 µg/mm <sup>3</sup>  |
| Fillers (particle size 1.5 µm)         | < 50%     | Water solubility DIN EN ISO 10477** | < 7.5 µg/mm <sup>3</sup> |
| Photoinitiators, pigments, stabilizers | < 0.2-10% | Viscosity**                         | 1400-2200 mPas           |

\*Applies to cured objects | according to norm at room temperature \*\*Applies to liquid resin

**Disclaimer**

Please note: Our products must be used in accordance with the instructions for use. We do not accept liability for any damage resulting from incorrect handling or usage. We do not accept liability for any damage resulting from improper use of the impression material. The user is furthermore obliged to check the product before use with regard to its suitability for the intended area of applications. We cannot accept any liability if the product is used in conjunction with materials and equipment from other manufacturers that are not compatible or not authorized for use with our product and this results in damage. Date of issue of these instructions for use: 2025-01

After the publication of this information for use, any previous versions become obsolete. The current version can be found at [www.vita-zahnfabrik.com](http://www.vita-zahnfabrik.com).

**Symbol Explanation**

| Graphic  | Title        | Description   | Reference   |
|--|--------------|---|---|
|  | Manufacturer | Indicates the medical device manufacturer.  | ISO 15223-1:2021<br>Medical devices - Symbols to be used with information to be supplied by the manufacturer - Part 1: General requirements |
|  | Use-by date  | Indicates the date after which the medical device is not to be used.                | ISO 15223-1:2021<br>Medical devices - Symbols to be used with information to be supplied by the manufacturer - Part 1: General requirements |
|  | Batch code   | Indicates the manufacturer's batch code so that the batch or lot can be identified. | ISO 15223-1:2021<br>Medical devices - Symbols to be used with information to be supplied by the manufacturer - Part 1: General requirements |

| Graphic  | Title                         | Description   | Reference   |
|--|-------------------------------|---|---|
|  | Distributor                   | Indicates the medical device manufacturer.  | ISO 15223-1:2021<br>Medical devices - Symbols to be used with information to be supplied by the manufacturer - Part 1: General requirements |
|  | Contains hazardous substances | Indicates a medical device that contains substances that can be carcinogenic, mutagenic, reprotoxic (CMR), or substances with endocrine-disrupting properties.                          | ISO 15223-1:2021<br>Medical devices - Symbols to be used with information to be supplied by the manufacturer - Part 1: General requirements |
|  | Medical device                | Indicates that the product is a medical device.   | ISO 15223-1:2021<br>Medical devices - Symbols to be used with information to be supplied by the manufacturer - Part 1: General requirements |
|  | Keep away from sunlight       | Indicates a medical device that needs protection from light sources.  | ISO 15223-1:2021<br>Medical devices - Symbols to be used with information to be supplied by the manufacturer - Part 1: General requirements |
|  | Temperature limit             | Indicates the temperature limits to which the medical device can be safely exposed.   | ISO 15223-1:2021<br>Medical devices - Symbols to be used with information to be supplied by the manufacturer - Part 1: General requirements |
|  | Consult instructions for use  | Indicates the need for the user to consult the instructions for use.  | ISO 15223-1:2021<br>Medical devices - Symbols to be used with information to be supplied by the manufacturer - Part 1: General requirements |
|  | Prescription only             | Requires a prescription in the United States of America Caution: In the United States of America, federal law restricts this device to sale or use by, or on the order of, a physician. | USA Code of Federal Regulations<br>21CFR Part 801<br>§ 801.109 (b)(1)   |
|  | Magnetic resonance            | The freeprint crown is MR Safe.   | FDA Guidance document: Testing and Labeling Medical Devices for Safety in the Magnetic Resonance (MR) Environment dated May 20, 2021.       |
|  | Caution                       | The uncured material is considered as hazardous substance.  | 29 CFR 1910.1200<br>Occupational Safety and Health Standards  |
|  | This way up                   | This is the correct upright position of the distribution packages for transport and/or storage.   | ISO 780:2015(E)<br>Packaging-Distribution packaging- Graphical symbols for handling and storage of packages                                 |

VITA Zahnfabrik has been certified, and the product **VITA VIONIC® DENT RESIN** bears the mark:

CE 0483

MD  R Only

Made in Germany

## **VITA Zahnfabrik H. Rauter GmbH & Co. KG**

Spitalgasse 3  
79713 Bad Säckingen  
Germany

Phone: +49 7761 562-0  
Fax: +49 7761 562-299  
E-Mail: [info@vita-zahnfabrik.com](mailto:info@vita-zahnfabrik.com)

## **detax GmbH**

Carl-Zeiss-Str. 4  
76275 Ettlingen  
Germany

Phone: +49 7243 510-0  
Fax: +49 7243 510-100  
E-Mail: [post@detax.com](mailto:post@detax.com)

922-13053\_0125\_V01

# VITA