n!ce^{®1} Zirconia GT

CAD/CAM-Material

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n!ce®¹ Zirconia GT

CAD/CAM material

Please read the instructions for use before using the product.

n!ce® Zirconia GT are zirconia blanks with matched system components for reliable shade reproduction. They can be used for the production of fully/partially veneered reconstructions and monolithic bridge restorations in the anterior and posterior tooth regions.

Intended purpose

n!ce® Zirconia GT products are ceramic materials for dental treatments.

Patient target group

No restrictions.

Intended user

Dental professionals only: dentists and dental technicians (Rx only).

Indications for use

n!ce® Zirconia GT is approved for:

- fully-anatomical crowns and up to 14-unit** bridges* in the anterior and posterior tooth regions,
 - fully and partially veneered single-tooth and up to 14-unit** bridge substructures* in the anterior and posterior tooth regions,
 - single tooth restorations and up to 14-unit** bridges* on directly screwed implant abutments in the anterior and posterior tooth regions,
 - inlays***, onlays***, veneers***, partial crowns***, occlusal veneers (table tops)***.
- * bridges and bridge substructures with a maximum of two adjoining bridge pontics.

** n!ce® Zirconia GT is limited in Canada to bridge indications with a max. of six units with a max. of two adjoining bridge pontics.

*** for adhesive bonding only.

Contraindication

Contraindications for n!ce® Zirconia GT:

- for more than two adjoining bridge pontics
- for two or more cantilever bridge pontics
- for parafunctional activity for veneered restorations, in particular for teeth grinders and clenchers
- for insufficient preparation results
- for insufficient hard tooth substance
- in patients with allergies or sensitivities to the ingredients
- temporary seating of veneered restorations
- conventional or self-adhesive incorporation of inlays, onlays, veneers, partial crowns and occlusal veneers (table top)

Product application/control

Workflow:

1. Preparation, shade determination
2. Intraoral scan or taking the impression / fabricating the model
3. CAD process (design)
4. CAM process (milling/grinding)
5. Reworking (removal/grinding of the restoration)
6. Manual shading technique for non-colored restorations: immersion technique (optional)
7. Sintering
8. Individualization/characterization/glazing/polishing
9. Bonding

Sintering parameters

The sintering results are influenced by various parameters such as model, performance and age of the sintering unit used. Detailed information on processing and sintering procedure can be found in the document Instructions for Use I full version.

n!ce® Zirconia GT Universal for VITA ZYRCOMAT 6000/6100MS (incl. esthetic cooling)*															
T ₀ [°C]	R ₁ ↗ [°C/min]	T ₁ [°C]	H ₁ → [min]	R _{C1} ↘ [°C/min]	T _{C1} [°C]	H _{C1} → [min]	Lift _{C1} [%]	R _{C2} ↘ [°C/min]	T _{C2} [°C]	H _{C2} → [min]	Lift _{C2} [%]	R _{C3} ↘ [°C/min]	T _{C3} [°C]	H _{C3} → [min]	Lift _{C3} %
≤ 200	17	1450	30	-	1350	15	100	-	1050	15	100	-	200	-	100

* Light-colored and delicate restorations do not necessarily have to be cooled slowly.

n!ce® Zirconia GT Speed for VITA ZYRCOMAT 6000/6100MS**																		
T ₀ [°C]	R ₁ ↗ [°C/min]	T ₁ [°C]	H ₁ → [min]	R ₂ ↗ [°C/min]	T ₂ [°C]	H ₂ → [min]	R _{C1} ↘ [°C/min]	T _{C1} [°C]	H _{C1} → [min]	Lift _{C1} [%]	R _{C2} ↘ [°C/min]	T _{C2} [°C]	H _{C2} → [min]	Lift _{C2} [%]	R _{C3} ↘ [°C/min]	T _{C3} [°C]	H _{C3} → [min]	Lift _{C3} [%]
200	90	800	5	50	1450	30	-	750	-	100	-	600	-	70	-	400	-	40

** In the case of extensive and/or chromatic restorations, slow cooling is recommended. Slow cooling generally results in a more intense color.

n!ce® Zirconia GT Universal for other dental furnaces***											
T ₀ [°C]	R ₁ ↗ [°C/min]	T ₁ [°C]	H ₁ → [min]	R _{C1} ↘ [°C/min]	T _{C1} [°C]	H _{C1} → [min]	Lift _{C1} [%]	R _{C2} ↘ [°C/min]	T _{C2} [°C]	H _{C2} → [min]	Lift _{C2} [%]
≤ 200	17	1520	30	-5	1000	-	100	-	200	-	100

*** Light-colored and delicate restorations do not necessarily have to be cooled slowly.

n!ce® Zirconia GT for industrial furnaces																		
T ₀ [°C]	R ₁ ↗ [°C/min]	T ₁ [°C]	H ₁ → [min]	R ₂ ↗ [°C/min]	T ₂ [°C]	H ₂ → [min]	R _{C1} ↘ [°C/min]	T _{C1} [°C]	H _{C1} → [min]	Lift _{C1} [%]	R _{C2} ↘ [°C/min]	T _{C2} [°C]	H _{C2} → [min]	Lift _{C2} [%]	R _{C3} ↘ [°C/min]	T _{C3} [°C]	H _{C3} → [min]	Lift _{C3} [%]
≤ 200	17	1150	30	8	1520	30-60****	-	1350	15	100	-	1050	15	100	-	200	-	100

**** Depending on the total load of the furnace, a longer holding time is recommended

n!ce® Zirconia GT for industrial furnaces alternatively (if regulated cooling is possible)			
R _{C1} ↘ [°C/min]	T _{C1} [°C]	H _{C1} → [min]	Lift _{C1} [%]
-5	1000	-	100

¹ n!ce® is a registered trademark of Straumann Holding AG.

Design parameters:

Minimum wall thickness in mm		
All information refers to sintered restorations		n!ce® Zirconia GT
Inlay, onlay, veneer	incisal occlusal circumferential	0.6 0.6 0.5
Anterior and posterior crowns (fully anatomical or substructure)	incisal occlusal circumferential	0.6 0.6 0.5
Fully anatomical anterior bridges and substructures with one pontic	incisal circumferential	0.6 0.6
Fully anatomical posterior bridges and substructures with one pontic	occlusal circumferential	0.7 0.6
Fully anatomical multi-unit anterior bridges and substructures with two pontics	incisal circumferential	0.8 0.6
Fully anatomical multi-unit posterior bridges and substructures with two pontics	occlusal circumferential	0.8 0.6
Cantilever bridges with an extension	incisal occlusal circumferential	0.8 0.8 0.6

Connector cross-sections in mm ²	
All information refers to sintered restorations	n!ce® Zirconia GT
Inlay, onlay, veneer	-
Anterior and posterior crowns (fully anatomical or substructure)	-
Fully anatomical anterior bridges and substructures with one pontic	9
Fully anatomical posterior bridges and substructures with one pontic	12
Fully anatomical multi-unit anterior bridges and substructures with two pontics	12
Fully anatomical multi-unit posterior bridges and substructures with two pontics	15
Cantilever bridges with an extension	15

The minimum wall thicknesses refer to fully sintered restorations.

The cantilever bridge unit should be modeled approximately one third narrower in its mesiodistal dimension.

Recommended product combinations

For recommended product combinations with the VITA system components and corresponding processing recommendations and restrictions regarding the combination possibilities of the n!ce® Zirconia GT products, please refer to the Instructions for Use | full version.

Storage/disposal

Store the n!ce® Zirconia GT discs in the original packaging and in a dry place.

Technical data

Property	Unit	n!ce® Zirconia GT
Coefficient of thermal expansion (20-500 °C)	10 ⁻⁶ · K ⁻¹	approx. 10.2
Chemical solubility (ISO 6872)	µg/cm ²	< 20
Density after sinter firing	g/cm ³	approx. 6.05
3-point flexural strength (ISO 6872)	MPa	approx. 850 MPa (incisal) - approx. 1200 MPa (cervical)

Chemical composition

Component	n!ce® Zirconia GT
ZrO ₂	86 - 93
Y ₂ O ₃	6 - 10
HfO ₂	1 - 3
Al ₂ O ₃	0 - 1
Pigments	0 - 2

Safety note

- While work is in progress, wear suitable safety goggles/face protection
- Only perform work under an extraction unit
- Wear protective gloves

Information on bonding of n!ce® Zirconia GT restorations in the patient's mouth can be found at www.vita-zahnfabrik.com/adiva.

Product reliability

Information on reporting serious incidents in connection with medical devices, general risks associated with dental treatments, residual risks and (if applicable) short clinical safety and performance reports (SSCPs) can be found at www.vita-zahnfabrik.com/product_safety.

The safety data sheets can be downloaded at www.vita-zahnfabrik.com or requested by fax at (+49) 7761-562-233.

Disclaimer

Please note: Our products must be used in accordance with the instructions for use. We accept no liability for any damage resulting from incorrect handling or usage. The user is furthermore obliged to check the product before use with regard to its suitability for the intended area of application. We cannot accept any liability if the product is used in conjunction with materials and equipment from other manufacturers that are not compatible or not authorized for use with our product and this results in damage. The VITA Modulbox is not necessarily a component of the product. Date of issue of these instructions for use: 2025-11

After the publication of this information for use, any previous versions become obsolete.

The current version can be found at www.vita-zahnfabrik.com.

VITA Zahnfabrik has been certified and the following products bear the mark:

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CH REP VITA Zahnfabrik H. Rauter GmbH & Co.KG, Bad Säckingen (Germany)
Zweigniederlassung Basel c/o Perrig AG, Max Kämpf-Platz 1, 4058 Basel
VITA Zahnfabrik · H. Rauter GmbH & Co. KG
Institut Straumann AG