VITA ENAMIC® GLAZE

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Produktbeschreibung

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Working instructions

Product description

VITA ENAMIC GLAZE is a clear transparent light-curing varnish for reliable surface sealing of restorations that have been customized using VITA ENAMIC STAINS shades and made from VITA ENAMIC, the hybrid ceramic with a dual-network structure. Following polymerization **within a spectral range of 350 – 500 nm**, this extremely thin liquid varnish provides a homogenous and highly abrasion-resistant surface, which, thanks to a special surface energy, significantly resists adhesion of biofilm.

Indication

Varnish for surface sealing of VITA ENAMIC hybrid ceramic restorations on which shade characterization has been performed using VITA ENAMIC STAINS.

Contraindication

Due to the thin layers it creates, VITA ENAMIC GLAZE is not suitable as a substitute for polishing.

Contraindications

In the case of persons who are allergic to VITA ENAMIC GLAZE product components.

Technical data

Chemical composition: Methyl methacrylate, acrylic resin, pentaerythritol tetraacrylate, photoinitiator

Processing

For optimum coverage and retentive adhesion, the surfaces to be sealed should be roughened, oil-free and dry. There should not be any residual milling and lubricating fluid on the surface. So that the varnish can be applied evenly and homogenously, the packaging includes 20 disposable Microbrush applicators. These ensure optimum distribution of even very small amounts of material on the surface of the VITA ENAMIC restorations on which shade characterization has been performed.

Do not use on polished surfaces!

Conditioning the surface

The surface must be conditioned as follows:

- Sandblasting using Al₂O₃, at max. 50 µm and a blasting pressure of max. 1 bar or
 - etching with a 5% hydrofluoric acid gel such as VITA CERAMICS ETCH as follows: Using a small disposable brush, apply VITA CERAMICS ETCH to the surfaces to be etched. Etching interval: 60 seconds Once the application time has come to an end, completely remove any residual acid from the etched surface by spraying down with copious amounts of water, by cleaning thoroughly using a steam jet device, or by cleaning in an oil-free ultrasonic bath using distilled water.
- Do not brush off, as this would lead to significant surface contamination.

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Surfaces sandblasted using Al₂O₃ must also be cleaned thoroughly. The roughened surface must then be silanized, e.g. using VITASIL.

The conditioned surface should no longer be touched.

Place a small drop of VITA ENAMIC GLAZE on the ceramic mixing palette and close the bottle again immediately.

Important:

The sensitive photoinitiator contained in VITA ENAMIC GLAZE may already be activated by lighting at the work surface or by OR lighting. Using the disposable MicroBrush applicator, quickly apply a single coat of VITA ENAMIC GLAZE to the entire restoration to create a thin and even layer. Do not allow areas of excess glaze to form.

Polymerization*

The varnish can be set using any conventional light-curing device intended for dental applications **with a spectral range of 350 – 500 nm**. All coated surfaces must be fully polymerized.

- Using polymerization lamps for dental clinic

Manufacturer	Device	Programs	Duration ¹⁾
Ivoclar Vivadent	Bluephase C8	High (800 mW/cm²)	60 seconds
Ivoclar Vivadent	Astralis 10	HIP (1200 mW/cm²)	60 seconds
Ultradent	Valo LED	Xtra power (3200 mW/cm²)	30 seconds

- Using polymerization devices for dental laboratories

Manufacturer	Device	Duration ¹⁾
Dreve	Polylux PT	10 - 15 minutes
Hager & Werken	Speed Labolight	10 minutes
Heraeus Kulzer	Heraflash	3 x 180 seconds

* Due to statutory provisions regarding product liability, details cannot be provided on light polymerization devices from other manufacturers. The quality of polymerization depends on factors that include the luminous intensity and wavelength of the radiated light.

Important:

¹⁾Duration times indicated are minimum times. To prevent irritation of the mucosa, correct procedure must be observed during polymerization of VITA ENAMIC GLAZE. Once set, the surface should no longer be sticky. Please make sure that the luminous intensity and wavelength of the polymerization device used enables complete hardening.

Enthält: Methylmethacrylat, CAS: 80-62-6 Pentaerythrittetraacrylat, CAS: 4986-89-4 Pentaerythrittriacrylat, CAS: 3524-68-3 Acrylharz Photoinitiator

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Highly flammable liquid and vapor. Health hazard Suspected of damaging fertility or the unborn child. Causes skin irritation. Causes serious eye irritation. May cause an allergic skin reaction. May cause respiratory irritation.

The corresponding safety data sheet can be downloaded from www.vita-zahnfabrik.com or requested via fax from (+49) 7761-562-233.

Storage

Storage: At room temperature (15°C – 28°C). Keep the bottle tightly closed. Protect from light.

Shelf life:

When stored correctly, until the expiry date (24 months).

The material should not be used after the expiry date.

Labeling

The expiry date and the lot number are shown on the container. Form of delivery: 1 20 ml bottle of VITA ENAMIC GLAZE with a dropper mechanism and 1 pack of 20 disposable Microbrush applicators, prod. no. EENSTG20

This product has been developed for use in dentistry and must be used in accordance with the instructions for use. The information/illustrations provided do not constitute a guarantee of features and are not binding. Information updated: 01.2013 Rx only

Manufactured by

Innovation MediTech GmbH, Max-Planck-Str. 31, 59423 Unna / Germany
Tel.: +49/(0)2303 8807-0 Fax: +49/(0)2303 8290-9
www.innovation-meditech.com
info@innovation-meditech.com

MediTech GmbH is certified according to the Medical Device Directive. The following product bears the mark **CE** 0044

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VITA Zahnfabrik H. Rauter GmbH & Co. KG