**Autor / Author**  
 Anwender / User  
**Status / Status**  
 Aktuell / Newsworthy  
**Kategorie / Category**  
 Anwenderbericht /  
 Case report

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### VITABLOCS RealLife

## Practicability test of a “novel” CAD/CAM material

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Fig. 1: VITABLOCS RealLife have a three-dimensional block structure.

In June 2010, the company VITA Zahnfabrik (D-Bad Säckingen) launches a new ceramic material for CAD/CAM production. The novel VITABLOCS RealLife combine tradition and innovation. The blocks are based on the fine-structure feldspar ceramic material VITABLOCS, which have proven themselves millions of times. The innovative part of the block concept is the three-dimensional block structure, which imitates the buildup of the natural tooth (Fig. 1).

VITABLOCS RealLife were specifically developed for highly esthetic anterior restorations. They are processed in the CAD/CAM systems CEREC or inLab MC XL (Sirona Dental Systems, D-Bensheim). Precondition for the use of the new blocks is the version 3.80 of the design software, which will also be available from June 2010.

In order to test the practicability of the “novel” block concept a zirconia crown as well as a crown made from VITABLOCS RealLife were produced for a real patient case (Fig. 2).

### Initial situation

In February 2010, an 18-year-old female patient presented in the practice to receive a restoration of tooth 21, primarily due to esthetic reasons (Figs. 3 and 4). After an accident in 2000, the tooth had been restored with an edge build-up and re-

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First information about VITABLOCS RealLife was published in DIGITAL\_DENTAL\_NEWS, issue April 2009, p. 10.

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
Abb. 2: Testweise angefertigte Krone.

Fig. 2: Test crown.

ceived an endodontic treatment a few months later. However, a crown was not produced before 2010, because growth of the patient was not yet completed.

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Fig. 4: Tooth 21 had received an endodontic treatment.

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## Preparation and digitization

In April 2010, a part of the existing root canal filling was removed to produce an inverted build-up following the shape of an endo crown. For this, the tooth surface was treated with the light-curing dentin / enamel bonding agent Syntac (Ivoclar Vivadent, FL-Schaan). The build-up was done with the universal composite Tetric EvoCeram (Ivoclar Vivadent). The patient decided against an internal bleaching for reasons of cost. Amongst other things due to a discolored tooth stump, the decision was made to produce a zirconia crown. According to the guidelines for all-ceramic restorations, a chamfer preparation was performed after having placed the retraction cord (Fig. 5).

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A subgingival preparation was not regarded as adequate in this patient case, since the patient has a thin and narrow tight gingiva and a shift of the bone margin in apical direction may have a negative impact on the biological width. If a new restoration is required at a later time, an even deeper preparation becomes necessary. After having informed the patient about the situation, she was content with the preparation of an equigingival finishing line and also ready to tolerate a possible difference in color at the transition between the crown and the tooth stump, since this area is covered by her upper lip.

With the polyether impression material Impregum (3M ESPE, D-Seefeld), a conventional impression was taken and a temporary restoration produced with Dentalon Plus (Heraeus Kulzer, D-Hanau). For the computer-aided production of the final restoration, a scan model was fabricated using super hard cast material CAM-base (dentona, D-Dormund) which was utilized to digitize the preparation with the scanner inEos (Sirona) (Fig. 6).

## Design and production

The design with the modeling software of the system inLab (Sirona) follows the proven work steps. To enable a comparison, two crowns were designed – one anatomically reduced for the crown with a zirconia coping and one in full contour for the crown made from VITABLOCS RealLife (Fig. 7).

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Fig. 6: Digitized model with marked finishing line.

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## Block structure

While the polychromatic VITABLOCS TriLuxe and TriLuxe forte are composed of several different layers of material, VITABLOCS RealLife have a three-dimensional block structure with dentine core and enamel coat: Due to the curved shade transitions between dentin and incisal edge, the natural play of colors of anterior teeth is imitated exactly. The software automatically places the planned restoration on the enamel-dentin margin of the virtual block. Thus the restoration is completely covered with enamel. Subsequently, the CEREC or inLab user can individually modify the position and thus precisely adjust the restoration to the residual dentition. For this, the position of the construction can be modified as desired in all three spatial dimensions (Figs. 8 to 12).

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Fig. 9: ... the construction is placed in the virtual block.

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Fig. 10: The positioning of the construction can be modified in all three dimensions ...

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Fig. 11: ... and individually adapted to the residual dentition of the patient.

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Fig. 12: When the final position is achieved the data is transferred to the milling machine.

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## Comparison

Figure 13 shows the crown produced from VITABLOCS RealLife with the milling machine in Lab MC XL on the right handside, and the crown with zirconia coping on the left handside. In figure 14, the opaque core of the zirconia restoration (l.) becomes clearly visible compared to the feldspathic crown (r.) due to transmitted light.

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Fig. 13: Crown made from VITABLOCS RealLife (r.) and crown with zirconia coping (l.).

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Abb. 14: Crown with zirconia coping (l.) and crown made from VITABLOCS RealLife (r.) in transmitted light.

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Fig. 15: Crown with zirconia coping ...

At try-in, it became apparent that the restoration with zirconia coping could be assessed as suboptimal (Figs. 15 and 16). The result could have been optimized by an adjustment of shade, shape and transparency. However, the crown made from the new three-dimensional block convinced with a great performance and thus, the original plan was dismissed and it was decided to permanently cement the "test" crown. Prior to placement, the restoration made from VITABLOCS RealLife was finalized: The distal incisal edge was lengthened and a high-gloss polishing was done (Figs. 17 to 20). To allow for an optimal material effect, finally, the crown was definitively bonded with transparent Multilink (Ivoclar Vivadent).

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Fig. 16: ... at the try-in in the patient's mouth.

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Fig. 19: Crown after lengthening of the distal incisal edge ...

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Fig. 20: ... and final polishing.

## Material characteristics

The fine-structure feldspar ceramic material VITABLOCS consist of regularly achieves excellent results in clinical studies. For example, a survival rate of 94 % was found after seven years in a survey of G. Christensen (Christensen, G.: Clinical trial revealed a 94 % success rate for VITABLOCS Mark II. In: CRA News, April 2006). Additionally, the ceramic material exhibits enamel-like abrasion properties and, due to standardized industrial production of the blocks, a highly homogeneous microstructure with consistent material quality is ensured.

VITABLOCS RealLife are available in the geometry RL-14/14. This corresponds to dimensions of 14 x 14 x 18 mm. The shade offering according to the standard VITA SYSTEM 3D MASTER includes IM2C, 2M2C as well as 3M2C. From 2011, three additional shades – 0MIC, 1MIC as well as 2MIC – will complement the portfolio according to the manufacturer.

## Conclusion

The doubts that – due to a higher translucency of feldspar ceramic material compared to zirconia – the dark tooth stump of the endodontically treated tooth might lead to an esthetic constraint was, as shown by the case example, unfounded. Dentist, dental technician and patient were pleasantly surprised by the result, particularly in direct comparison with the restoration that had been originally planned. The three-dimensional block structure offers an outstanding esthetic result (Fig. 21).

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Fig. 21: The patient was pleasantly surprised by the final result.

The correct positioning of the crown construction in the milling preview seems to be slightly complex at the first attempt, but turns out to be easily practicable after a short period of training thanks to a simple use of the software. With a view to the CAD/CAM materials for esthetically demanding areas which have been available on the market so far, VITABLOCS RealLife represent a quantum leap. In this way, the practiced CEREC or inLab MC XL user can very easily and quickly produce highly esthetic restorations.

My thanks go to my nieces Dr. Ute Werling for the dentist and MDT Cathrin Werling for the technological implementation of this patient case.

DDN online Auch online unter: [www.ddn-online.net](http://www.ddn-online.net) ☆ ▼

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- 1983-1985 training as a dental technician
- 1985-1989 studies of dentistry in Frankfurt / Main, Germany
- 1992 establishment of a joint practice with his sister Ursula Werling in Bellheim, Germany
- since 1993 national and international lecturer for implantology
- 2005 establishment of a special laboratory for CAD/CAM technology
- supervision and visiting courses related to the topics of implantology and CAD/CAM
- scientific consultant for dental companies and medical faculties

### Kontakt / Contact

[dr.werling@t-online.de](mailto:dr.werling@t-online.de)