

VITABLOCS® for Rapid Layer Technology

For the highly efficient CAD/CAM veneering of bridge frameworks with CEREC/inLab MC XL



VITA shade taking

VITA shade communication

VITA shade reproduction

VITA shade control

Date of issue: 07.14

VITA shade, VITA made.

VITA

Product description

VITABLOCS Mark II and TriLuxe forte blocks are available in the large geometry I-/TF-40/19. This block geometry was specially developed for CAD/CAM veneering (VITA Rapid Layer Technology) of multi-unit bridge sub-structures. This technology enables both the framework (oxide ceramic) and veneer (fine-structure feldspar ceramic) to be fabricated

in a digital manufacturing process. The harmoniously matched shade nuances of VITABLOCS TriLuxe forte facilitate a veneer structure with a natural colour gradient. With VITA Rapid Layer Technology, the framework and veneer component are adhesively bonded with composite.

VITABLOCS® for Rapid Layer Technology

Indications

VITABLOCS Mark II and TriLuxe forte I-/TF-40/19 are suitable for the manufacture of CAD/CAM veneers of posterior bridge substructures with up to four units.

Geometry

VITABLOCS Mark II and TriLuxe forte I-/TF-40/19 are available in the geometry I-/TF-40/19 (15.5 x 19 x 39 mm).

Features/benefits

Time savings:

Compared to the classical layering technique, this innovative procedure offers a great time savings, since framework and veneer are both manufactured in an entirely computer-aided process.

Simple to bond:

The framework and veneer are attached to one another simply and securely by means of an adhesive composite. The more complex and time-consuming ceramic sintering procedure is no longer necessary.

High degree of reliability:

The framework and veneer component are perfectly matched to enable optimum processing and clinical reliability.

Wide spectrum of indications:

The combination of VITABLOCS Mark II with TriLuxe forte I-/TF-40/19 and VITA In-Ceram YZ-40/19 enable the fabrication of fully anatomical bridge restorations with up to four units.

VITABLOCS Mark II for CEREC/inLab MC XL* – available items

Shade	Designation	Size	Content	Art. No.
1M1C	I-40/19	15.5 x 19 x 39 mm	2 pcs	EC41M1CI40192
1M2C	I-40/19	15.5 x 19 x 39 mm	2 pcs	EC41M2CI40192
2M2C	I-40/19	15.5 x 19 x 39 mm	2 pcs	EC42M2CI40192
3M2C	I-40/19	15.5 x 19 x 39 mm	2 pcs	EC43M2CI40192

VITABLOCS TriLuxe forte for CEREC/inLab MC XL* – available items

Shade	Designation	Size	Content	Art. No.
1M2C	TF-40/19	15.5 x 19 x 39 mm	2 pcs	EC41M2TF40192
2M2C	TF-40/19	15.5 x 19 x 39 mm	2 pcs	EC42M2TF40192
3M2C	TF-40/19	15.5 x 19 x 39 mm	2 pcs	EC43M2TF40192

* VITABLOCS Mark II and TriLuxe forte for the Rapid Layer Technology can be used with the Sirona CEREC/ inLab MC XL system from the inLab 3D software version V3.80 onwards. IMPORTANT: CEREC users require the inLab 3D software in order to use this technique.