## VITA YZ HT SHADE LIQUID

**D**

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## VITA YZ HT SHADE LIQUID

VITA YZ HT SHADE LIQUID is perfectly tailored to highly-translucent VITA YZ HT<sup>White</sup> for coloring using the brush technique. The best results are obtained with the coloring liquids when this combination is used. The liquids can also be used in combination with VITA YZ HT<sup>Color</sup> restorations. However, due to the pigments that are already present, deviations in terms of shade may occur prior to sintering.

### Indication:

For complete or partial coloring of VITA YZ HT restorations.

### Using the brush technique:

- Prior to final sintering, air-blow the VITA YZ HT framework structure.
- Protect the structure against damp prior to coloring.
- Always shake the bottle thoroughly (VITA YZ HT SHADE LIQUID) before use!**  
If the storage location is too warm, vapors could develop inside the bottle and pressure could build up. In this case, open the bottle carefully under a suction hood. You can then continue working as normal.
- Dip a metal-free brush (YZ HT SHADE LIQUID Brush) fully into the liquid and then wipe off on the edge of the container or dab gently with a tissue.  
Proceed in the same way each time liquid is applied, otherwise, differences in shade could occur.
- Then paint the structure using the brush in accordance with the pre-defined pattern. In order to prevent the liquid from being thinned out and to avoid contamination, always wash and dry the brush thoroughly after use and whenever a different color is used.
- SDry the structure under an appropriate infrared lamp (250 watts) for at least 45 minutes at 70 °C. In the case of larger, dense structures, a longer pre-drying time is recommended. A longer pre-drying time does not affect the shade result.  
Do not perform sintering while moist, as differences in shade could occur as well as problems during sintering!
- After coloring, always close the liquids tightly in order to prevent them drying out.

Tooth shade	Shade Liquid	Number of brush strokes*		
		Cervical	Body	Incisal
A1	A1	4	3	2
A2	A2	4	3	2
A3	A3	4	3	2
A3.5	A3.5	4	3	2
B2	B2	4	3	2
C2	C2	4	3	2
D2	D2	4	3	2
1M1	1M1	3	2	1
1M2	1M2	3	2	1
2L1.5	2L1.5	3	2	1
2M2	2M2	3	2	1
3M2	3M2	4	3	2
3M3	3M3	4	3	2
4M2	4M2	3	2	1

\* All statements are applicable for VITA YZ HT<sup>White</sup>. These are reference values that may vary depending on brush pressure, handling, and on the amount of liquid.

**Note:** Keep the brush and the liquid away from metal (e.g. shaping instruments, devices, etc.) in order to prevent contamination. If using additional shades, make sure to clean the brush thoroughly in order to prevent contamination with other shades.

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**GB**

### VITA YZ HT SHADE LIQUID Stabilizer:

After some time, the zirconia may become blotchy during coloring. This is due to evaporation of the stabilizer in the liquid. In this case, you should once again add 5% of the stabilizer to the SHADE LIQUID.

### VITA YZ HT SHADE LIQUID Indicator (adding shade pigments to the mixture):

Over time, the shade pigments that have been mixed with the product in order to achieve improved perceptibility of the dentin shade will evaporate (e.g. if stored in too bright a location). These shade pigments can be added again using SHADE LIQUID Indicator (1 cm of liquid each in the bottle corresponds to 1-2 drops of SHADE LIQUID Indicator).

### Storage and shelf life:

The storage temperature may not exceed 25°C / 77°F; short-term storage at higher temperatures up to max. 45°C / 113°F is permitted (e.g. during transport).

VITA YZ HT SHADE LIQUID must be stored in a dark and cool location, as otherwise the shade pigments lose their effect and must be added again more frequently (see the section on adding shade pigments to the mixture). Ideally, the liquid should be stored at less than 18°C / 64°F. If the storage location is too cold, the liquid may freeze and must be thawed out again at room temperature before use.

Do not use after the expiry date.

### Safety measures:

- Use gloves
- Wear protective goggles with side shields
- Ensure sufficient ventilation during use
- Ensure sufficient ventilation while drying

### Safety information:

- In the case of contact with the eyes, immediately rinse the eyes with clean running water for a few minutes with the eyelids open. We recommend consulting a physician.
- In the case of contact with the skin, rinse with plenty of soap and water. Consult a physician if the skin becomes irritated.
- Swallowing: If swallowed, do not induce vomiting or drink any water. Immediately consult a physician or attend a hospital. Present the safety data sheet.
- Inhalation: In the case of inhalation, bring the affected person outside into the open air. Ensure that fresh air is provided and bring affected persons away from the danger area. Consult a physician.

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**Please note:** Our products must be used in accordance with the instructions for use. We accept no liability for any damage resulting from incorrect handling or usage. The user is furthermore obliged to check the product before use with regard to its suitability for the intended area of application. We cannot accept any liability if the product is used in conjunction with materials and equipment from other manufacturers that are not compatible or not authorized for use with our product. Furthermore, our liability for the accuracy of this information is independent of the legal basis and, in as far as legally permissible, shall always be limited to the value as invoiced of the goods supplied, excluding value-added tax. In particular, as far as legally permissible, we do not assume any liability for loss of earnings, indirect damages, ensuing damages or for third-party claims against the purchaser. Claims for damages based on fault liability (culpa in contrahendo, breach of contract, unlawful acts, etc.) can only be made in the case of intent or gross negligence. Date of issue of these instructions for use: 03.15

Further information at: [www.vita-zahnfabrik.com/cadcam](http://www.vita-zahnfabrik.com/cadcam)

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Proveedor : VITA Zahnfabrik · H. Rauter GmbH & Co. KG · Spitalgasse 3 · 79713 Bad Säckingen, Alemania
Zirkonzahn GmbH · An der Ahr 7 · 39030 Gais · Italy
Fornitore : VITA Zahnfabrik · H. Rauter GmbH & Co. KG · Spitalgasse 3 · 79713 Bad Säckingen, Germania
VITA Zahnfabrik H. Rauter GmbH & Co.KG