VITA Zirconia T / HT / ST / XT

Recommodation for CNC machines

Machining mode : Milling

Information and tips

Tools

Machining strategy

Parameters

VITA

Version: 09.04.18

Information

The information presented here, are intended as a recommendation. Depending on the available CNC machines, CAM software, tools, etc. the information have to be adapted to your own production situation. As a result, different results may obtained.

The development of the strategies and parameter was done with following system:

- DMG MORI Ultrasonic 20 und Imesslcore CoriTech 350i
- CAM Software: Hyperdent V8.2 Beta

According to this recommendation, a coping bridge (tooth 11-17) can be finished in 01:10:00 h (14mm Disk), with a good surface and fit.

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www.emuge-franken.de

Tips for YZ

- Plunge into the material with a ramp (5-10 degree) or helically and use a reduced plunge feed (feed Z)
- YZ should always be milled dry with air cooling
- The diameter of the restoration holding pin should range from 1,5 - 2 mm

Strategy

- A two side machining and 3+2 strategies are sufficient in most cases.
- In order to maintain a good fit, even by restoration with undercuts, the last finishing of the cavity should be done with a 5 axis strategy.
- In order to maintain a good occlusal fit, the complete occlusal side should be finished with max. a $\varnothing 1.2\text{mm}$ tool (or less). A special finishing of the fissures isn't necessary.
- If chipping occurs, the feed speed and step should be lowered. Also the Oversize should be checked. Too much or too less oversize will also lead to chipping.

Recommended tools


| Diameter | Description | Manufacturer | Order-code |
|---------------------------|------------------------------------|--------------|--------------------|
| $\varnothing 2\text{ mm}$ | Ball nose end mill, Diamond coated | Zecha | 421.B2.200.100.200 |
| $\varnothing 1\text{ mm}$ | Ball nose end mill, Diamond coated | Zecha | 421.B2.100.050.100 |
| $\varnothing 2\text{ mm}$ | Ball nose end mill, Diamond coated | Franken | 2600E.200616 |
| $\varnothing 1\text{ mm}$ | Ball nose end mill, Diamond coated | Franken | 2600E.100612 |


Tool life


| Tool | Units | Restoration |
|-----------------------------------|-------|---------------------------------|
| Zecha - 2 mm Ball nose end mill | >150 | Fully anatomical crown tooth 26 |
| Zecha - 1 mm Ball nose end mill | >150 | Fully anatomical crown tooth 26 |
| Franken - 2 mm Ball nose end mill | >150 | Fully anatomical crown tooth 26 |
| Franken - 1 mm Ball nose end mill | >150 | Fully anatomical crown tooth 26 |


Order of machining


| Step | Machining side | Machining | Tool |
|------|----------------|----------------------------------|-------|
| 1 | Cavity side | Roughing | Ø 2mm |
| 2 | Occlusal side | Roughing | Ø 2mm |
| 3 | Cavity side | Finishing in cavity | Ø 2mm |
| 4 | Cavity side | Finishing outside cavity | Ø 2mm |
| 5 | Occlusal side | Finishing | Ø 2mm |
| 6 | Cavity side | Residual material in cavity | Ø 1mm |
| 7 | Cavity side | Residual material outside cavity | Ø 1mm |
| 8 | Margin line | Finishing | Ø 1mm |
| 9 | Cavity side | Finishing in cavity | Ø 1mm |
| 10 | Occlusal side | Finishing, Fissures | Ø 1mm |
| | | | |


| Step 1 | Cavity side - Roughing | | | | 3+2 axis |
|---|------------------------|-------|-------|--------|----------|
|  | Tool | Ø 2mm | | | notes: |
| | Tolerance | 0,01 | | | |
| | Spindel speed | [n] | 36000 | rpm | |
| | Feed speed XY | [Vf] | 2000 | mm/min | |
| | Feed speed Z | [Vf] | 1000 | mm/min | |
| | Width of cut XY | [ae] | 1,4 | mm | |
| | Depth of cut Z | [ap] | 1,0 | mm | |
| | Oversize | | 0,4 | mm | |


| Step 2 | Occlusal side - Roughing | | | | 3+2 axis |
|--|--------------------------|-------|-------|--------|----------|
|  | Tool | Ø 2mm | | | notes: |
| | Tolerance | 0,01 | | | |
| | Spindel speed | [n] | 36000 | rpm | |
| | Feed speed XY | [Vf] | 2000 | mm/min | |
| | Feed speed Z | [Vf] | 1000 | mm/min | |
| | Width of cut XY | [ae] | 1,4 | mm | |
| | Depth of cut Z | [ap] | 1,0 | mm | |
| | Oversize | | 0,4 | mm | |


| Step 3 | Cavity side - Finishing in cavity | | | | 3+2 axis |
|---|-----------------------------------|-------|-------|--------|----------|
|  | Tool | Ø 2mm | | | notes: |
| | Tolerance | 0,01 | | | |
| | Spindel speed | [n] | 36000 | rpm | |
| | Feed speed XY | [Vf] | 1300 | mm/min | |
| | Feed speed Z | [Vf] | 1000 | mm/min | |
| | Width of cut XY | [ae] | 0,2 | mm | |
| | Depth of cut Z | [ap] | - | mm | |
| | Oversize | | 0,1 | mm | |


| Step 4 | Cavity side - Finishing outside cavity | | | | 3+2 axis |
|---|--|-------|-------|--------|----------|
|  | Tool | Ø 2mm | | | notes: |
| | Tolerance | 0,01 | | | |
| | Spindel speed | [n] | 36000 | rpm | |
| | Feed speed XY | [Vf] | 1500 | mm/min | |
| | Feed speed Z | [Vf] | 1000 | mm/min | |
| | Width of cut XY | [ae] | 0,2 | mm | |
| | Depth of cut Z | [ap] | - | mm | |
| | Oversize | | 0,0 | mm | |


| Step 5 | Occlusal side - Finishing | | | | 3+2 axis |
|---|---------------------------|-------|-------|--------|----------|
|  | Tool | Ø 2mm | | | notes: |
| | Tolerance | 0,01 | | | |
| | Spindel speed | [n] | 36000 | rpm | |
| | Feed speed XY | [Vf] | 1300 | mm/min | |
| | Feed speed Z | [Vf] | 1000 | mm/min | |
| | Width of cut XY | [ae] | 0,2 | mm | |
| | Depth of cut Z | [ap] | - | mm | |
| | Oversize | | 0,1 | mm | |

| Step 6 | Cavity side - Residual material in cavity | | | | 5 axis |
|--|---|-------|-------|--------|--------|
|  | Tool | Ø 1mm | | | notes: |
| | Tolerance | 0,01 | | | |
| | Spindel speed | [n] | 45000 | rpm | |
| | Feed speed XY | [Vf] | 600 | mm/min | |
| | Feed speed Z | [Vf] | 300 | mm/min | |
| | Width of cut XY | [ae] | 0,2 | mm | |
| | Depth of cut Z | [ap] | 0,1 | mm | |
| | Oversize | | 0,1 | mm | |

| Step 7 | Cavity side - Residual mat. outside cavity | | | | 5 axis |
|---|--|-------|-------|--------|--------|
|  | Tool | Ø 1mm | | | notes: |
| | Tolerance | 0,01 | | | |
| | Spindel speed | [n] | 45000 | rpm | |
| | Feed speed XY | [Vf] | 600 | mm/min | |
| | Feed speed Z | [Vf] | 300 | mm/min | |
| | Width of cut XY | [ae] | 0,2 | mm | |
| | Depth of cut Z | [ap] | 0,1 | mm | |
| | Oversize | | 0,0 | mm | |

| Step 8 | Margin Line - Finishing | | | | 5 axis |
|---|-------------------------|-------|-------|--------|--------|
|  | Tool | Ø 1mm | | | notes: |
| | Tolerance | 0,01 | | | |
| | Spindel speed | [n] | 45000 | rpm | |
| | Feed speed XY | [Vf] | 1000 | mm/min | |
| | Feed speed Z | [Vf] | 1000 | mm/min | |
| | Width of cut XY | [ae] | 0,1 | mm | |
| | Depth of cut Z | [ap] | - | mm | |
| | Oversize | | 0,0 | mm | |

| Step 9 | Cavity side - Finishing in cavity | | | 5 axis | |
|---|-----------------------------------|-------|-------|--------|--------|
|  | Tool | Ø 1mm | | notes: | |
| | Tolerance | 0,01 | | | |
| | Spindel speed | [n] | 45000 | | rpm |
| | Feed speed XY | [Vf] | 1300 | | mm/min |
| | Feed speed Z | [Vf] | 1000 | | mm/min |
| | Width of cut XY | [ae] | 0,2 | | mm |
| | Depth of cut Z | [ap] | - | | mm |
| | Oversize | | 0,0 | | mm |

| Step 10 | Occlusal side - Finishing, Fissures | | | 3+2 axis | |
|--|-------------------------------------|-------|-------|----------|--------|
|  | Tool | Ø 1mm | | notes: | |
| | Tolerance | 0,01 | | | |
| | Spindel speed | [n] | 45000 | | rpm |
| | Feed speed XY | [Vf] | 1300 | | mm/min |
| | Feed speed Z | [Vf] | 1000 | | mm/min |
| | Width of cut XY | [ae] | 0,15 | | mm |
| | Depth of cut Z | [ap] | - | | mm |
| | Oversize | | 0,0 | | mm |

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ZECHA Hartmetall - Werkzeugfabrikation GmbH
FRANKEN GmbH & Co. KG - Fabrik für Präzisionswerkzeuge