VITA CEREC POWDER, VITA CEREC LIQUID & VITA CEREC PROPELLANT

D

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VITA CEREC POWDER, VITA CEREC LIQUID & VITA CEREC PROPELLANT

GB

VITA CEREC POWDER is sprayed (using the VITA CEREC spray head and VITA CEREC PROPELLANT) onto the tooth surface and adjacent teeth after they are coated and prepared with VITA CEREC LIQUID. VITA CEREC POWDER conditions the surface for the opto-electronic impression with the CEREC method to achieve a perfect picture (scan) quality for the design of the restoration.

Composition

- VITA CEREC POWDER
 - Mixture of zinc stearate, titanium dioxide and zirconium dioxide
- VITA CEREC PROPELLANT
 - Medical CFC-free propellant
- VITA CEREC LIQUID
 - Aqueous, polysorbate-containing adhesive for the CEREC POWDER

Application areas

- VITA CEREC POWDER

Intraoral use:

Contrast powder for the opto-electronic impression with the CEREC intraoral camera to achieve perfect diffuse scattering on the surface of the tooth substance (dentine and enamel).

Extraoral use:

Contrast powder for model materials such as plaster die to allow capturing using a scanner (e.g. inEos or inLab).

- VITA CEREC PROPELLANT
 - Propellant for intraoral and extraoral use for the application of the VITA CEREC POWDER.
- VITA CEREC LIQUID
 - Adhesive for the CEREC powder to be used intraorally to ensure adhesion of the powder to the tooth substance and easy removal of the powder later on. It can be completely removed using water.

Contraindication

For patients with allergic reactions to constituents of VITA CEREC POWDER, VITA CEREC LIQUID and VITA CEREC PROPELLANT.

Side effects

Side effects are not yet known.

Hazard warnings / Safety recommendations

VITA CEREC PROPELLANT

- Extremely flammable aerosol. Pressurised container:
 - May burst if heated.
- Keep away from heat/sparks/open flames/hot surfaces.
 - No smoking.
- Pressurized container: Do not pierce or burn, even after use.
- Do not spray on an open flame or other ignition source.
- Protect from sunlight. Do not expose to temperatures exceeding 50° C/122° F.

Application

• Preparation of the treatment area

The treatment area needs to be prepared to ensure perfect control (e.g. rubber dam) prior to taking the opto-electronic impression with the CEREC intraoral camera.

Optional: Place a rubber dam to allow a more clearly arranged treatment area which can be kept dry during adhesive cementation of the ceramic restoration.

• Application of VITA CEREC LIQUID

Apply the VITA CEREC LIQUID onto the preparation and the adjacent teeth using a disposable brush. Then disperse the liquid to obtain a thin film. As a result the liquid is spread precisely and enhanced adhesion of the POWDER to the tooth substance is achieved. Also, it can be removed more easily later.

• Application of VITA CEREC POWDER

Screw the VITA CEREC spray head onto the bottle containing the VITA CEREC POWDER and then attach it to the VITA CEREC PROPELLANT bottle.

Note:

The POWDER must be stored dry. Do not shake or tilt the propellant bottle prior to the use to avoid the discharge of liquid propellant and the risk of clumping of the POWDER!

Before CEREC POWDER is applied the POWDER bottle should be held vertically and tapped lightly several times to remove deposited powder from discharge hole.

First test the application of the powder outside the oral cavity. Apply a thin, uniform coat of the contrast powder using several spray blasts. This way excessive stress relief cooling of the propellant is avoided, which may be caused by clumping of the POWDER in the spray head owing to moisture deposited from the air.

Note:

Please adhere also to the respective instructions in the processing manuals of Sirona Dental Systems.

• Cleaning of the tooth surfaces after successful optical impression

After the optical impression, the contrast agent must be cleaned using the water jet and/or the air syringe. The surface of the tooth must be completely clean.

Additionally, rotary dental brushes may be used for cleaning.

Storage

Store at room temperature. VITA CEREC PROPELLANT:

protect from sunlight and temperatures above 50° C/122° F.

Disposal

Dispose of containers only when empty.

Marking

The expiry date and the lot number are printed on the containers.

Delivery forms

- VITA CEREC POWDER, Bottle, cont. 12 g. Prod. No. ECPO
- VITA CEREC LIQUID, Bottle, cont. 6 ml. Prod. No. ECL
- VITA CEREC PROPELLANT, Bottle, cont. 250 ml. Prod. No. ECPN
- VITA CEREC spray head for PROPELLANT
 - The cannula can be turned by 360°. Prod. No. ECS

VITA Zahnfabrik is certified according to the Medical Device Directive and the following VITA products bear the mark CE:

VITA CEREC® POWDER

VITA CEREC® LIQUID

VITA CEREC® PROPELLANT

CEREC®, inLab® and inEos® are registered trademarks of Sirona Dental GmbH, A-Wals.

The safety datasheet can be downloaded at www.vita-zahnfabrik.com or requested by sending a fax to (+49) 7761-562-233.

The product has been developed for the use in the dental sector and has to be processed in accordance with the instructions for use. Illustrations and written information do not include any undertaking as to qualities and are not binding. Date of issue of this information: 05.2015

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