VITA SIMULATE Preparation Material

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Produktbeschreibung

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Working Instructions

Product description

VITA SIMULATE Preparation Material is a light-curing composite for the fabrication of artificial dies to simulate the shade of the prepared tooth. The material is used especially for the fabrication of translucent all-ceramic restorations in which the shade effect is mainly influenced by the color of the die. The product allows dental technicians to reproduce the tooth shade more easily and safely since the shade of the restoration can be controlled in advance and adjusted if required.

Processing

The dentist uses the VITA SIMULATE Preparation Material Guide to determine the shade of the prepared die. The guide includes 6 VITA SYSTEM 3D-MASTER shades. "S" stands for VITA SIMULATE Preparation Material.

0M1 S,* 1M1 S, 2M3 S, 3M2 S, 4M3 S, 5M3 S

*for bleached dies.

Once the dentist has communicated the selected shade to the dental technician, the technician uses the respective shade of VITA SIMULATE Preparation Material to fabricate an artificial control die in the following way.

- Use a brush to apply VITA SIMULATE Insulation Liquid without puddle formation in a thin, homogeneous layer to the interior of the ceramic crown or veneer.
- Fill the interior of the ceramic restoration with VITA SIMULATE Preparation Material and condense the material using a modeling instrument in order to avoid cavity formation.
- Press the application stick into the uncured die material and ensure that the tip of the application stick is positioned in the center of the restoration and does not touch its sides.
- Fill the remaining cavities and adapt the material firmly around the application stick. Close syringe immediately after use!
- Then remove the excess die material from the area of the margins using a modeling instrument.
- Cure the die material with a light-curing unit or a hand-held light-curing unit (see table).

The light-curing unit must be equipped with a light source which has a wavelength of 350-500 nm and a maximum intensity of 470 nm. These values are exclusively based on correctly functioning equipment. The manufacturers' information on the unit, in particular the maintenance instructions, must be observed!

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- Remove the cured VITA SIMULATE die from the restoration and clean the restoration in the ultrasonic unit. Then the VITA SIMULATE die is used to check whether the shade of the completed restoration corresponds to the desired shade. If necessary, shade adjustments can be carried out by staining with VITA AKZENT Plus or layering with VITA VM 9 or VITA VM 7 (depending on the CTE).

Hazard information

Contains triethylene glycol dimethacrylate (CAS: 109-16-0), 2-dimethylaminoethyl methacrylate (CAS: 2867-47-2).

- Causes severe irritation of eyes and skin.
- May cause allergic skin reaction.
- Wear protective gloves / protective clothing / safety goggles.
- Wash thoroughly after use.
- In case of contact with skin: wash with copious amounts of water and soap.
- In case of contact with eyes: rinse carefully with water for some minutes. Remove contact lenses if worn and easy to do so. Continue rinsing.
- In case of persistent eye irritation: seek medical advice/consult a doctor.

The respective safety data sheets can be downloaded at www.vita-zahnfabrik.com or requested by fax at (+49) 7761-562-233.

Storage and durability

Do not store above 25°C. Do not expose to direct sunlight. It is generally recommended to store composites under cold conditions. To ensure perfect storage of the material, it should be stored in the firmly sealed container/syringe in a suitable refrigerator at 5-7°C. The material should not be used after the expiry date.

Labeling

The expiry data and the lot numbers are shown on the syringes.

Forms of delivery:

VITA SIMULATE Preparation Material KIT, Prod. No. ESPKIT

Contents:

6 syringes cont. 4 g of VITA SIMULATE Preparation Material each
 1 x 0M1 S, Prod. No. E440314
 1 x 1M1 S, Prod. No. E440344
 1 x 2M3 S, Prod. No. E440404
 1 x 3M2 S, Prod. No. E440464
 1 x 4M3 S, Prod. No. E440544
 1 x 5M3 S, Prod. No. E440594
 1 bottle of VITA SIMULATE Insulation Liquid, 16 ml, Prod. No. ESIL16
 150 pcs. VITA SIMULATE application sticks, Prod. No. E017
 1 VITA SIMULATE Preparation Material Guide (shade guide), Prod. No. B430

The product has been developed for use in the dental sector and must be processed in accordance with the instructions for use. The information/illustrations contained herein do not constitute a guarantee of features and are not binding. Date of issue: 01.2013

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VITA Zahnfabrik H. Rauter GmbH & Co.KG