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GB

Important

Both tube and bottle must be closed immediately after use – failure to observe this rule will lead to thickening of the hardener liquid.

Slight cloudiness or flocculation of the hardener liquid is not important and will not affect its use.

Although the bottle of hardener liquid in the standard assortment is only half full, it does nevertheless contain sufficient hardener liquid for the preparation of the full contents of the accompanying tube of paste.

If the paste has separated out during a longer period of storage this can be rectified by lightly kneading the flexible plastic tube.

Working Instructions

VITAFOL H should be dosed out as shown in the illustration, and then mixed together before the flasking of the fully waxed-up denture.

1. Squeeze sufficient paste out of the tube onto the mixing block or a glass slab with the rings underneath depicted here above left to fill the required circle:

outer circle I = sufficient for 28 teeth
middle circle II = sufficient for 14 teeth
inner circle III = sufficient for 7 teeth

2. Add the relevant amount of hardener liquid:

8 to 9 drops for circle I (28 teeth)
4 to 5 drops for circle II (14 teeth)
2 to 3 drops for circle III (7 teeth)

3. Now mix together and spatulate thoroughly for about 1 minute.

4. Using the brush provided, apply the VITAFOL H in an evenly thin layer either only onto the teeth that have to be insulated together with the immediately surrounding areas of wax, or else onto the entire waxed-up denture taking especial care to completely fill the embrasures.

5. Whilst this insulating film is still wet, sprinkle retention crystals evenly onto it. At normal room temperature the insulating film will take roughly 8 to 12 minutes to harden. ("Hard" being definable as being when the surface is no longer tacky.) If either the room temperature is warmer than normal, or more than the required amount of hardener liquid is used, then the hardening process will be accelerated. Conversely, if the temperature is cooler or less than the required amount of hardener is used, then the hardening process will be slowed down.

6. The flasking of the waxed-up denture and boiling away of the wax are both carried out in the normal way.

7. When deflasking, the teeth and the denture surfaces that have been insulated with VITAFOL H will come cleanly and

easily away from the plaster, meaning that time-consuming trimming can also be dispensed with.

VITAFOL H Hardener

Danger!

Warning!

Contains: tetraethyl silicate, ethyl silicate, bis(acetyloxy)diocetylstanann.

Flammable liquid and vapour.
Harmful if inhaled.
May cause respiratory irritation.
May cause damage to organs through prolonged or repeated exposure.
Causes serious eye irritation.

VITAFOL H Paste

Danger!

Warning!

Highly flammable liquid and vapour.
Causes serious eye irritation.
Work under an extraction unit.
Wear safety goggles / protective gloves / protective clothing.

Please refer to the safety data sheets for more detailed information! These can be downloaded at www.vita-zahnfabrik.com or ordered by fax under (+49) 7761-562-233.

The product was developed for use in the dental branch and must be processed according to the instructions for use. The information/illustrations contained herein do not constitute an assurance of properties and are not binding.

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VITA Zahnfabrik H. Rauter GmbH & Co.KG