VITA SPRAY-ON

Das System / The system / Le Système / El Sistema rociador / Sistema

VITA shade determination

VITA shade communication

VITA shade reproduction

VITA shade control


Edition 05.17

VITA – perfect match.

VITA

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Replacement parts list

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The following safety instructions should be observed when working with the VITA SPRAY-ON SYSTEM:

The liquids (INDICATOR and SPRAY-ON LIQUID), that are used for mixing together with the opaques for the spraying technique contain approx. 90% ethyl alcohol (ethanol).

Hazard description:

F highly flammable

Ethanol is an easily flammable, volatile liquid and can produce explosive vapor/air mixtures. If swallowed or inhaled in large amounts it may cause intoxication and symptoms of poisoning.

The following safety precautions should, therefore, be observed for storing solutions and mixtures containing ethanol and also for the spray technique:

- Keep the lid closed tightly and store in a cool and well ventilated room.
- Flammable liquids belonging to groups A1 and B (these include ethanol) may only be stored in 1 liter containers for use on the job (manual use).
- Protect against heat and keep away from possible sources of ignition.
- Observe – No smoking – restriction.
- Never spray into open flames or onto glowing objects.
- When working without a spray box, always wear safety glasses and gas filter mask (gas filter type A; equivalent to a spray painting mask); Spray against a moist towel or appropriate filter for absorbing the spray particles; always provide adequate ventilation.
- If the VITA SPRAY-Box is to be connected to an exhaust system or is to be used under an extractor, a high exchange of air must be guaranteed so that no dangerous solvent concentrations can be produced.

The VITA SPRAY-ON SYSTEM

Introduction

The VITA SPRAY-ON System offers the following advantages:

1. Easy to use.
2. Reduction of working time by up to more than 50%.
3. Saving of materials and costs, due to the storing of already mixed materials in closed retainers, ready for use when required.
4. Thin, well-proportioned layers of opaque porcelain.
5. More room to build up the other layers of porcelain, especially in areas that pose problems.
6. Optimum bond between opaque porcelain and metal alloy.
7. Depending on the technique employed, a variety of surface textures can be achieved, thereby allowing variation in the absorption, diffusion and reflection of light.

Installation

For installation and connection, see diagram page 11.

Adjusting the sprayer needle

The sprayer needle is very soft. Therefore, care must be taken that it is not bent. Successful spraying without clogging is dependant on the nozzle opening. Coarser-grained porcelains require a wider opening. So the needle should be pushed further back using the adjustment screw. When spraying fine-grained opaque porcelains, the needle should be positioned more to the front.

To move the needle forward, the adjusting screw and the locking nut situated at the rear of the spray gun are turned clockwise. Turning them anti-clockwise pulls the needle back.

Note: Should the spray gun clog up after a while, widen the nozzle opening so as to allow an unobstructed flow of the opaque porcelain. Just move the adjusting screw and the locking nut clockwise (that is to the left) by a half or a full turn, until the opaque porcelain flows freely.

Mixing the opaque porcelain

Fill the respective glass retainer up to a **maximum** of $\frac{1}{3}$ of its capacity. The vent hole of the adapter lid clogs up easily when the opaque porcelain is shaken. A blocked vent hole can lead to the development of a vacuum and impede proper spraying.

The more opaque porcelain is in the glass retainer, the more it tends to form a sediment at the bottom and become heavy bodied. For perfect spraying the opaque porcelain should have a watery, milk-like consistency.

The composition of the VITA SPRAY-ON LIQUID is such that it allows a light-bodied mixture. As it evaporates very quickly, more than one layer of opaquer can be applied without having to predry the coated framework in a furnace each time. After mixing the opaque porcelain powder with VITA SPRAY-ON LIQUID, the adapter lid is attached to the glass retainer and then connected to the spray gun.

Normal operating pressure is 0.7 - 1.0 bar. Attach a retainer containing a fresh opaque porcelain mixture. Hold the sprayer button depressed. If after two seconds no opaque porcelain is ejected from the nozzle, increase the pressure to between 1.1 and 1.4 bar. Should the opaque porcelain be sprayable, it ought to be carefully diluted until it can be sprayed at a pressure of 0.7 bar. If the opaque porcelain is still not sprayable at a pressure between 1.1 and 1.4 bar, the grain of the opaque porcelain is most probably too coarse for spraying with the needle in the current position. Retract the needle a little. Generally, you can adjust the pressure to suit yourself (Range: approx. 0.5 to 4.5 bar). The mixture is too thin, when spraying at a distance of 2 to 5 cm the opaque porcelain on the framework is runny. Should this be the case, simply add more opaque porcelain powder.

Spray at a distance of between 6 and 8 cm. **Example:** Spray opaque porcelain onto a larger coin from a distance of less than 6 cm. The opaque layer on the coin is very moist. If the distance is too small, the low pressure will lead to a wavy pattern on the surface. If you spray at a distance of more than 10 to 15 cm, the opaque porcelain will be almost dry when it hits the surface of the coin.

The opaque porcelain should display a slightly moist sheen, when it hits the surface of the framework. It should, however, not run. This allows the buildup of several layers of opaque porcelain. If the opaque porcelain hits the framework in too dry a state, the fired surface will be rougher. The rougher the surface is, the more light will be reflected in scattered form. The metal framework should in no case be heated before the opaque porcelain is applied.

When switching to a different shade of opaque porcelain, clean the spray gun by just spraying water until it is clear. Then continue spraying until no more water is ejected from the spray gun. You may now use a new shade of opaque porcelain. This procedure should also be followed when the spray gun is clogged up, whereby it might be necessary to retract the needle a little.

Application and firing

First, spray 1 to 3 thin (wash) layers onto the framework, and then place it in a furnace for firing. The wash layer can be applied using a single colour for all shades, as complete covering is not achieved. In doing so, you shorten this phase of work.

A further coating is applied until the framework displays a homogeneous shade and no metal shines through after the second firing. When spraying, always direct the spray gun from the occlusal plane to the gingival region. This is equivalent to an angle of 45° to the longitudinal axis and warrants no opaque porcelain being applied to the inner surface of the framework.

When applying opaque porcelain to a bridge, follow the above procedure correspondingly: Apply opaque porcelain to the entire bridge framework, with the exception of the base of the pontics.

For spraying the base surfaces, we recommend filling the abutment crowns with cotton wool. This prevents the opaque porcelain entering the inner surfaces of the framework. After taking out the cotton wool,

the inner surfaces can be cleaned with a **dry brush**.

Spraying Color porcelains

1. Mix 30% normal opaque porcelain and COLOR porcelains with VITA INDICATOR LIQUID to a more watery consistency than normal opaque porcelains.
2. Reduce the air pressure, so that spraying is just possible. The nozzle opening can be narrowed as the mixture is more liquid.
3. Spray lightly at a distance of approximately 1.5 cm and aim the spray gun downwards, away from the edge of the crown.
4. For a better control of the spraying process, self-made stencils can be employed.
5. After application of the opaque has been completed, the spray gun is cleaned by spraying while rinsing with water. The same applies for the adapter lid.

Note on maintenance

In the case of hard water it is recommended to rinse the pistol out occasionally with a calcium-dissolving agent.

It may be necessary from time to time to replace the Teflon needle guide and the Teflon seal for the spray nozzle head (B 109) if material does not flow smoothly.

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
With the unique VITA SYSTEM 3D-MASTER® all natural tooth shades are systematically determined and completely reproduced.

Please note: Our products should be used according to the working instructions. We cannot be held liable for damages resulting from incorrect handling or usage. The user is furthermore obliged to check the product before use with regard to its suitability for the intended area of applications. We cannot accept any liability if the product is used in conjunction with materials and equipment from other manufacturers which are not compatible or not authorized for use with our product. Furthermore, our liability for the correctness of this information is independent of the legal ground and, in as far as legally permissible, is limited to the invoiced value of the goods supplied excluding turnover tax. In particular, as far as legally permissible, we do not assume any liability for profit loss, for indirect damages, for consequential damages or for claims of third parties against the purchaser. Claims for damages based on fault liability (culpa in contrahendo, breach of contract, unlawful acts, etc.) can only be made in the case of intent or gross negligence.
Date of issue of these directions for use: 05.17

After the publication of these working instructions any previous versions become obsolete. The current version can be found at www.vita-zahnfabrik.com

VITA

VITA Zahnfabrik H. Rauter GmbH & Co.KG
Spitalgasse 3 · D-79713 Bad Säckingen · Germany
Tel. +49(0)7761/562-0 · Fax +49(0)7761/562-299
Hotline: Tel. +49(0)7761/562-222 · Fax +49(0)7761/562-446
www.vita-zahnfabrik.com · info@vita-zahnfabrik.com
 facebook.com/vita.zahnfabrik