VITA ENAMIC® STAINS KIT

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Produktbeschreibung

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Working instructions

Product description

The VITA ENAMIC STAINS KIT, which contains 6 different shade powders as well as accessories, is used for shade characterization (staining) of the surfaces of hybrid ceramic restorations made from VITA ENAMIC with its dual-network structure. The shades in these sets polymerize thanks to a combination of two curing systems that are adapted to suit one another. Autopolymerization is activated by mixing shade powder and liquid together. Light polymerization with a **wavelength of 350 – 500 nm** stimulates and completes final polymerization. This dual-curing system guarantees a long period of use, as well as ensuring that shade mixtures that are extremely opaque can set all the way through. Once the VITA ENAMIC STAINS have been applied, the surface on which shade characterization has been performed is sealed using VITA ENAMIC GLAZE and homogenized, increasing the durability of the shade layer in the moist oral environment.

Indication

System for shade customization of the surfaces of hybrid ceramic restorations made from VITA ENAMIC with its dual-network structure.

Contraindication

In the case of persons who are allergic to ingredients contained in the VITA ENAMIC STAINS KIT product components.

Technical data

Product	Components
VITA ENAMIC STAIN	cristobalite, dibenzoyl peroxide, dicyclohexyl phthalate
VITA ENAMIC STAINS LIQUID	methyl methacrylate, aromatic urethanacrylate
VITA ENAMIC GLAZE	methyl methacrylate, 2-Propenoic acid, reaction product with Pentaerythrite, Diphenyl(2,4,6-trimethylbenzoyl)phosphinoxide

Processing

Conditioning the surfaces to be characterized

For optimum coverage and retentive adhesion of the shade, the surface of the ENAMIC restoration to be characterized must be roughened and oil-free. Do not use on polished surfaces! The surface must be conditioned as follows:

- Sandblasting using Al₂O₃, at max. 50 µm and a blasting pressure of max. 1 bar or
- etching with a 5% hydrofluoric acid gel such as VITA CERAMICS ETCH as follows: Using a small disposable brush, apply VITA CERAMICS ETCH to the surfaces to be etched. Etching interval: 60 seconds Once the application time has come to an end, completely remove any residual acid from the etched surface by spraying down with copious amounts of water, by cleaning thoroughly using a steam jet device, or by cleaning in an oil-free ultra sonic bath using distilled water. Do not brush off, as this would lead to significant surface contamination.

Surfaces sandblasted using Al₂O₃ must also be cleaned thoroughly. The roughened surface must then be silanized, e.g. using VITASIL. The conditioned surface should no longer be touched.

Mixing the shades

Using a glass or acrylic mixing spatula, mix corresponding shade powder on the ceramic mixing palette with VITA ENAMIC STAINS LIQUID. The mixing ratio can be varied as required, as long as the mixture can still be easily spread. Varying the mixture allows the intensity of the hue to be adjusted, creating results ranging from watery/transparent through to opaque hues that offer full coverage.

Processing times

Period of use supported by the final shade mixture: approx. 10 minutes. The period of use can be extended by diluting with liquid. In the case of longer breaks between work, it is recommended that the mixing palette be covered using a lightproof lid in order to prevent the mixture from polymerizing prematurely. Shade mixtures that have dried out can be prepared using liquid without affecting the quality.

Applying the shade

Apply the shade using the brush provided (size no. 1) or using a disposable brush. If required, intermediate polymerization can be performed after each individual layer of shade is applied.

Important: Once staining has been completed and before sealing the surface using VITA ENAMIC GLAZE, the applied shade must be polymerized, as the shade layer could otherwise be inadvertently diluted as a result of subsequent sealing using VITA ENAMIC GLAZE.

Intermediate polymerization of the applied shade layer*

- Using polymerization lamps for dental clinic

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Manufacturer	Device	Programs	Duration
Ivoclar Vivadent	Bluephase C8	High (800 mW/cm²)	30 seconds
Ivoclar Vivadent	Astralis 10	HIP (1200 mW/cm²)	20 seconds
Ultradent	Valo LED	Xtra power (3200 mW/cm²)	6 seconds

- Using polymerization devices for dental laboratories

Manufacturer	Device	Duration
Dreve	Polylux PT	5 minutes
Hager & Werken	Speed Labolight	5 minutes
Heraeus Kulzer	Heraflash	2 x 180 seconds

* Due to statutory provisions regarding product liability, details cannot be provided on light polymerization devices from other manufacturers. The quality of polymerization depends on factors that include the luminous intensity and wavelength of the radiated light.

Important: The staining brushes are cleaned using VITA VM LC CLEANER. Before further use, brushes must be dried thoroughly, e.g. on tissue, as residual CLEANER on the brush can contaminate stain mixtures, preventing them from being correctly polymerized.

Sealing the applied shade using VITA ENAMIC GLAZE

The use of VITA ENAMIC GLAZE is mandatory as it is used for sealing and homogenizing surfaces on which shade characterization has been performed and increases the durability of the shade layer in the moist oral environment. Place a small drop of VITA ENAMIC GLAZE on the ceramic mixing palette and close the bottle again immediately.

Important: The sensitive photoinitiator contained in VITA ENAMIC GLAZE may already be activated by lighting at the work surface or by OR lighting. Using the disposable MicroBrush applicator provided, quickly apply a single coat of VITA ENAMIC GLAZE to the entire restoration to create a thin and even layer. Do not allow areas of excess glaze to form.

Final polymerization of the shade and VITA ENAMIC GLAZE*

The varnish can be set using any conventional light-curing device intended for dental applications **with a spectral range of 350 – 500 nm**. All coated surfaces must be fully polymerized.

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Manufacturer	Device	Programs	Duration ¹⁾
Ivoclar Vivadent	Bluephase C8	High (800 mW/cm²)	60 seconds
Ivoclar Vivadent	Astralis 10	HIP (1200 mW/cm²)	60 seconds
Ultradent	Valo LED	Xtra power (3200 mW/cm²)	30 seconds

- Using polymerization devices for dental laboratories

Manufacturer	Device	Duration ¹⁾
Dreve	Polylux PT	10 - 15 minutes
Hager & Werken	Speed Labolight	10 minutes
Heraeus Kulzer	Heraflash	3 x 180 seconds

* Due to statutory provisions regarding product liability, details cannot be provided on light polymerization devices from other manufacturers. The quality of polymerization depends on factors that include the luminous intensity and wavelength of the radiated light.

Important: ¹⁾The duration times indicated are minimum times. To prevent irritation of the mucosa, correct procedure must be observed during polymerization of VITA ENAMIC GLAZE. Once set, the surface should no longer be sticky. Please make sure that the luminous intensity and wavelength of the polymerization device used enables complete hardening.

Notes on risks

VITA ENAMIC STAIN

Contains: dibenzoyl peroxide, CAS: 94-36-0 May cause an allergic skin reaction. Contains: Cristobalite, CAS:14464-46-1 Causes damage to the lung through prolonged or repeated exposure. Route of exposure: Inhalative.

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Innovation MediTech GmbH, Max-Planck-Str. 31, 59423 Unna / Germany