
VITABLOCS® Guide 3D-MASTER®

Gebrauchsanleitung Working Instructions Mode d'emploi Instrucciones de uso Istruzioni di impiego

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Working Instructions for the VITABLOCS® Guide 3D-MASTER®

E

The VITABLOCS® Guide 3D-MASTER® with shade guide teeth made of original Mark II ceramic enables quick and accurate selection of the corresponding VITABLOCS®.

The shades of the VITABLOCS® Mark II, TriLuxe and TriLuxe forte are integral parts of the VITA SYSTEM 3D-MASTER®, which is the only tooth shade system on the market which takes into account all 3 dimensions of color and enables a systematic order principle to be applied to shade determination and shade reproduction: **value – chroma – hue**

VITABLOCS® are available in the following shades:

	Shade lightness (Value) →				
	light				dark
pale ↓ Intensity (Chroma)	0M1C	1M1C	2M1C	3M1C	–
	–	1M2C*	2M2C*	3M2C*	4M2C
	–	–	2M3C	3M3C	–
deep					

The shade samples of the VITABLOCS® Guide 3D-MASTER® are grouped linearly into lightness groups:

Group 0	Group 1	Group 2	Group 3	Group 4
0M1C	1M1C, 1M2C*	2M1C, 2M2C*, 2M3C	3M1C, 3M2C*, 3M3C	4M2C

VITABLOCS® TriLuxe and TriLuxe forte are only available in shades marked *.

Since the monochrome VITABLOCS® Mark II, and the individual layers of the VITABLOCS® TriLuxe und TriLuxe forte consist of a homogeneous mixture of translucent and dentine porcelains and are not layered consecutively like the shade sample teeth of the VITA Toothguide 3D-MASTER®, the corresponding block shade of the shade guide is approximated but not identical. This is expressed by the additional designation "C" after the shade designations for the VITABLOCS®, e.g. "2M2C".

The first "2" means **shade lightness group 2**, „M“ stands for the medium (neutral) **hue**, the second "2" stands for **shade intensity group 2**.

With the VITABLOCS® TriLuxe the middle layer (dentine or body layer) corresponds exactly to the appropriate shade of a monochrome VITABLOCS® Mark II. The incisal

or enamel layer is less intensively shaded, and hence more opaque than the middle layer. Thus the color gradient of the VITABLOCS® TriLuxe corresponds to that of natural teeth.

VITABLOCS® TriLuxe forte for even more esthetics...

VITABLOCS® TriLuxe forte are based on the clinically successful VITABLOCS® TriLuxe. The shade transition from the enamel to the cervical layer is even more finely nuanced, at the same time with greater emphasis of the chroma in the intensive, dentine or cervical layer. In connection with the increasing fluorescence towards the cervical area, this also ensures a more lifelike and natural shade effect which makes even better clinical results possible in the case of thin wall thicknesses.

Procedure

Using the shade sample and taking into account the lightness and intensity (chroma), determine the shade and take the corresponding block out of the assortment.

Tips for shade taking

- The shade guide should take place before preparation, since after preparation the tooth shade appears too white due to dehydration.
- If possible, determine tooth shades in daylight or under standardized daylight lamps (e.g. OSRAM LUMILUXE® DELUXE daylight 12-950) and not under OP lamps.
- Make your choice without hesitation; always accept the first decision as eye fatigue begins to set in after approx. 6-7 seconds.
- The entire surroundings should be kept in unobtrusive colors. That means no lipstick, no spectacles with tinted lenses and no brightly colored clothing (alternative: cover clothing with a grey cape).

Hygiene tips and instructions for care

The acrylic parts of your VITABLOCS® Guide 3D-MASTER® are made of a high-quality, easy-care material that is temperature resistant (up to 180 °C). The entire VITABLOCS® Guide 3D-MASTER® can be autoclaved at a temperature of up to max. 140 °C.

Disinfection

The VITABLOCS® Guide 3D-MASTER® can be surface disinfected. We recommend a customary, mild, spray disinfectant or disinfectant wipes. According to the manufacturers' product information these disinfectants are effective against HIV and HBV, and are available from most dental dealers. Please heed the respective manufacturer's instructions for use. Disinfectants containing phenol, compounds with phenylphenol groups or methyl ethyl ketone will damage the acrylic parts.

Note

The acrylic parts of the VITABLOCS® Guide 3D-MASTER® are neither monomer- nor UV-resistant, and must be protected from intense sunlight.

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