## VITA CAD-Temp® multiColor CTM-85/40

Kurzinformation

D

Stand: 02.16

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GB

Date of issue: 02.16

The VITA CAD-Temp® multiColor composite block comprised of acrylate polymer CTM-85/40 is used for the fabrication of long-term temporaries with a clinical wear time of up to 3 years\*. For information on the material indications, please see working instructions no. 1459. (Information on downloading is provided below).

In the software, block size CTM-85/40 (18 x 40 x 85 mm) is **only activated for milling using milling tools**. With regard to activation and usage of the new milling function as well as the capacity to process the block size, please note the applicable specifications and information provided by Sirona.

### Minimum software requirements:

- **CEREC Software ≥ 4.3**
- **inLab Software ≥ 4.2**

For inLab software 4.2, the inLab 4.2.5 service pack can be downloaded free of charge from the Sirona website at [www.sirona.com/de/inLabdownloads](http://www.sirona.com/de/inLabdownloads) in order to activate the new milling function and integrate the VITA CAD-Temp multiColor CTM-85/40 block.

### Basic hardware requirements for using the milling function:

- **CEREC MC XL Premium Package** (as of serial number 302001), **CEREC MC XL** (as of serial number 129001), or **inLab MC XL** (as of serial number 129001)
- Existing customers with milling units **that have an older serial number must replace a spindle motor before they can perform milling**. In the case of milling units with four motors, that means the left spindle motor in motor set 1; in the case of milling units with two motors, that means the left spindle motor. So that the spindle motors can be replaced, corresponding Sirona starter kits for milling are available from specialist dealers (REF 6448893 for design 2011 motors or REF 6451079 for closed motor brackets).
- **Existing customers with milling units in accordance with the serial numbers described above** already have devices equipped with new spindle motors and only require starter kits without a motor (REF 6454487) in order to begin using the new milling function.

Additional detailed instructions for use (working instructions, etc.) for VITA CAD-Temp are available from <https://www.vita-zahnfabrik.com/en/Download-centre-Product-information-612.html>

Thank you for placing your confidence in us – we wish you every success with VITA CAD-Temp!

\*Clinical study by the University of Tübingen.

Rx only

VITA Zahnfabrik has been certified in accordance with the Medical Device Directive and VITA CAD-Temp® bears the CE mark **CE** 0124.

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VITA Zahnfabrik H. Rauter GmbH & Co. KG