# VITA TITANKERAMIK

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For veneering of pure titanium  
and titanium alloys.  
Available in VITA SYSTEM 3D-MASTER and VITA classical A1 – D4 shades.

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## VITA TITANKERAMIK

is for veneering of pure titanium and titanium alloys ( $Ti_6Al_4V$ ).

### The special properties of Titan:

- titanium 1/degree of purity 99,5 - 99,8%
- Vickers hardness 180-250 (degree 1-4)
- thermal expansion coefficient  $9.6 \times 10^{-6} K^{-1}$
- biocompatible
- corrosion-resistant
- neutral taste / hypoallergenic
- high mechanical strength / lightweight
- low thermal conductivity

### Lattice transformation

The lattice transformation of titanium at 882 °C is of particular importance for the dental technician. Below this crystallisation temperature titanium occurs as hexagonal  $\alpha$ -titanium, and above this temperature as cubic body-centred  $\beta$ -titanium.

As well as making the material very brittle and impure, this transformation results in a 17% increase in volume.

### $\alpha$ -case

The surface of the titanium is hardened by diffusion of elements (oxygen, carbon, ...). This progression of hardness is caused by the reaction of the molten titanium with ingredients of the investment compound. The marginal layer becomes strongly embrittled and impure. The  $\alpha$ -case layer must be removed as a preparatory measure before veneering with dental ceramics.

### Passivation layer

Titanium has a strong tendency to spontaneously form a very fine oxide layer at room temperature. This protective layer prevents corrosion and is therefore responsible for the biocompatibility of the material. The passivation layer is not electron-conducting and self-regenerating.

**Adequate passivation must always be guaranteed after different working procedures.**

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
*When modelling crown and bridge frameworks it is always essential first and foremost to consider the low thermal conductivity of titanium and veneering ceramic. In the course of the entire working procedure both materials must be given the opportunity to give off heat at a constant rate. Providing corresponding cooling fins, which are modelled as interdental inlay-like fortifications has been proven to be beneficial. The use of special modelling materials is not required. When using modelling acrylics, a very low proportion of remaining monomer should be aimed at. The use of release spray is not required.*

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### **Correct processing**

Finishing is carried out with special, cross-cut burrs. The reduced wedge angle of this type of cutter leads to effective reduction of the material with effective cooling of the workpiece. Unsuitable use of processing materials may cause undesired oxide phases to form locally during finishing. In most cases the workpiece is characterised by a dark, discolored, slightly rough surface. No adequate adhesive bond can be achieved on these surfaces. The main priority here is to carefully remove shavings once more.

#### **Important:**

- slow speed (approx. 15.000 rpm)
- low application of pressure
- regularly clean the cutters using a steam cleaner
- the framework should be ground in one direction only
- sharp edges and overlapping should be avoided.
- clean, disposable corundum
- grain size 120-150 µm pressure of 2 bar
- angle of 45° to object

The passivating time of 5 minutes for the surface must then be adhered to.

Steam clean the surface.

Do not use acids or etching agents!

An oxidation firing is not required.

#### **TIP:**

Clean the cutters to prevent impurities

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The VITA BONDER enabled a breakthrough to be achieved with regard to increasing adhesive bond strength. The bonding of titanium oxides during firing procedures in the bonder matrix leads to a secure and durable adhesive bond. The titanium/ceramic adhesive bond is equal to conventional systems. There are different methods of applying the bonder:

- BONDER + OPAQUE LIQUID applied with glass spatula = difficult handling. Often results in incomplete coverage in the marginal areas, puddles may form.
- Bonder application using the VITA SPRAY-ON procedure = optimum coverage of the surface, easy handling
- Paste bonder. Optimal surface wetting. The paste bonder in the glass jar should be stirred with a plastic spatula before use. If, after longer periods of storage, the material has become too firm, its original consistency can be restored by adding certain amounts of Paste Opaque Liquid.

#### **Ceramic furnace:**

The ceramic furnace should be thoroughly tested in order to ensure the successful firing of VITA TITANKERAMIK.

#### **Requirements:**

- Ceramic furnaces used for firing different alloys are unsuitable.
- Thick oxide deposits may result in weakening of the bond and discoloration of the ceramic (discoloration of the firing tray)
- Furnaces subject to large temperature fluctuations are unsuitable
- Damaged firing trays as well as soiled furnace linings must first be exchanged.

The working parameters described by the manufacturer must be supplemented by the routine and experience of the dental technician. Temperature variations of as little as +/- 5°C can be clearly seen and may make all the difference between success and failure. The assessment of the surface glaze is the only certain way to ensure correct firing.

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#### **TIP:**

*Apply transparent ceramic to a titanium substructure already coated with opaque. If the edges are already broken after firing, the temperature is too high. If, however, the ceramic has a milky and inhomogeneous appearance, the temperature is too low. Raise or lower the temperature by 5°C at a time until the correct firing temperature is reached. Due to the low thermal conductivity of titanium/titanium ceramic materials each firing program begins with a standby temperature of 400°C.*

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**Important:**

- *The opaque must not be sprayed on too dry.*

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- *The recommended firing parameters must not be altered.*

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*The opaque can be applied with a brush or using the VITA SPRAY-ON procedure. The opaque application is always carried out in two firings. Both opaque firings must offer coverage.*

*1st opaque firing 790 °C 70 % coverage*

*2nd opaque firing 790 °C 100 % coverage*

***The classical wash opaque is not required here.***

**Firing recommendations:**

*All opaque firings should be carried out with full vacuum, incl. holding time. The surface of the fired opaque must not have an eggshell-like finish, but must be glossy. Slow cooling is not necessary in the case of these preparatory firings.*

**TIP:**

*Place the mixed opaque for the SPRAY-ON into the VITASONIC II (or another ultrasonic unit) for 30 seconds. For SPRAY-ON use only SPRAY-ON LIQUID.*

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*Previous working procedures, step-by-step*

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### *The main firings*

*The layering technique known from conventional VITA ceramics is still valid. Due to the higher translucency of the ceramics, it is necessary to apply the opaque dentine more thickly than usual.*

- *There is slightly more firing shrinkage than in conventional ceramics.*
- *Only provided original fluids should be used.*
- *VITA Modisol should be used as ceramic insulation.*
- *The ceramic powders should always be kept homogeneously moist.*
- *These can be re-wetted if necessary with distilled water.*
- *All firing programs have a standby temperature of 400 °C.*
- **All main firings and correction firings are carried out with full vacuum (that means heating time and holding time at end temperature). From the first dentin firings and all other subsequent firings it is recommended that all restorations should be slow cooled down to 400°C. Allow the restoration to remain inside the muffle during this cool down phase.**
- *No condensing is carried out during the entire layering process.*

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

Zervikal CO2 ochre  
Inzisal CO3 white

- *Individual effects can be added to the opaque-coated substructure using CO ceramics.*

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- *The opaque dentine is the shade-carrying base and covers the entire labial surface.*

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- *Individualisations can be carried out at this stage in order to increase the natural appearance and anatomical effect of the ceramic. The intensives IN1 - IN10 (MASTER KIT) can be used in pure or mixed form.*

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- *The VITA Interno ceramics offer a further possibility for individualisation (Int. 01 - Int. 12). **With VITA TITANKERAMIK do not use the Interno FLUID, but VITA MODELLING FLUID.***

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*The firing temperatures indicated are valid for the VITA VACUMAT furnaces. Important: when using a VITA VACUMAT 6000 **NO** operation with Quickstart.*

**From the opaque firing onwards all firing cycles should operate under full vacuum (incl. holding time). Slow cooling is recommended down to 400 °C.**

Brand / firing	Vt. °C	→ min.	↗ min.	↗ °C/min	Temp ca. °C	→ min.	↘ °C	→ min.	VAC °C/min
Pastenbonderbrand <i>Paste bonder firing</i>	400	6.00	6.00	67	830	1.00	–	–	7.00
Pulver Bonderbrand <i>Powder bonder firing</i>	400	2.00	6.00	67	830	1.00	–	–	7.00
Grundbrand <i>Opaque firing</i>	400	2.00	4.00	98	790	1.00	400**	0.00	5.00
Schultermasse <i>Shoulder firing</i>	400	6.00	7.00	53	770	1.00	400**	0.00	8.00
1. Dentinbrand <i>1<sup>st</sup> Dentine firing</i>	400	6.00	7.00	53	770	1.00	400**	0.00	8.00
2. Dentinbrand <i>2<sup>nd</sup> Dentine firing</i>	400	6.00	7.00	53	770	1.00	400**	0.00	8.00
Glanzbrand ohne Glasurmasse <i>Glaze firing without glaze</i>	400	0.00	4.00	93	770	1.00	400**	0.00	5.00
Malfarbenfixierbrand* <i>Stains firing*</i>	400	4.00	3.00	100	700	1.00	400**	0.00	–
Glanzbrand mit Glasurmasse Akz 25 <i>Glaze firing with AKZENT Plus glaze</i>	400	4.00	4.00	93	770	1.00	400**	0.00	5.00

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*When using dental ceramics, the firing result largely depends on the individual firing procedure of the user, i. e. among other aspects on the type of furnace, the location of the temperature sensor, the firing tray as well as the size of the workpiece during the firing cycles. Our application-technical recommendations (regardless whether they have been provided orally, in writing or in the form of practical instructions) are based on numerous own experiences and tests. The user, however, should consider this information only to provide basic values. If surface, transparency and degree of gloss should not correspond to the firing result that is achieved under optimal conditions, the firing procedure must be adjusted correspondingly. The crucial factors for the firing procedure are not the firing temperature displayed by the furnace but the appearance and the surface condition of the firing object after the firing process.*

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# CERAMICS OF THE VITA TITANKERAMIK MASTER KIT

## Information on application

### COLOR OPAQUE

CO 1	greyish-brown
CO 2	ochre
CO 3	white
CO 4	lilac
CO 5	grey
CO 6	pink
CO 7	light brown
CO 8	dark brown

### GINGIVA

G 1	light gingiva
G 2	dark gingiva

### MARGIN

The MARGIN porcelains (MAR) were developed in order to create an optimal aesthetic juncture between crown and tooth stump when the metal coping has been labially shortened in the area of the crown margin. The MARGIN porcelain applied can be hardened by applying heat. It is recommended to stabilize the shoulder by means of a hairdryer or using the heat emission from the firing chamber.

MARN Neutral MARGIN porcelain for mixing with MAR 1 - MAR 6

MAR 1	} The MARGIN porcelains are shade-matched to the VITA Toothguide 3D-MASTER® as well as the VITAPAN classical shade guide. For MARGIN classification see separate table.
MAR 2	
MAR 3	
MAR 4	
MAR 5	
MAR 6	

### CERVICAL

Applications: These materials are layered onto the dentine above the neck of the tooth extending into the approximal area to increase the illusion of depth.

CE 1	light orange	(light orange translucent)
CE 2	light yellow	(light-yellow translucent)
CE 3	gold	(gold-yellow translucent)
CE 4	clear	(neutral-translucent)
Creates the desired illusion of depth without influencing the color.		

### TRANSLUCENT

T 1	milky	(whitish translucent) for whitish incisal areas, margins and triangular protuberances
T 2	opal	(yellowish translucent) for ridges, wedges and covering the labial area
T 3	pearl	(orange translucent) for ridges, wedges and covering the labial area
T 4	clear	(neutral translucent)
T 5	polar	(light blue translucent) layered under the enamel material to create bluish incisal areas
T 6	blue	(blue translucent) for bluish wedges and ridges

### INTENSIVE

The INTENSIVES are used to apply small color effects in the dentine and enamel area. The choice of colors ranges from beige to yellow, orange, greyish blue and white.

IN 1	linen	(sandbeige)
IN 2	flesh	(yellow beige)
IN 3	corn silk	(corn yellow)
IN 4	pastel yellow	(light yellow)
IN 5	ochre	(ochre)
IN 6	gold earth	(orange)
IN 7	flame orange	(flame orange)
This very intensely colored porcelain is particularly useful for blending OPAQUE DENTINE, DENTINE and TRANSLUCENT.		
IN 8	brush	(brownish grey)
IN 9	cloud	(greyish blue)
A particularly fine-grained ceramic with a very intense shade of color. This ceramic is highly suitable for creating exceedingly thin, greyish-blue color effects to increase the illusion of depth in the incisal area.		
IN 10	snow	(white) For strong whitish effects. Particularly suitable for margins and triangular ridges in the occlusal area when mixed with DENTINE, ENAMEL or TRANSLUCENT materials.

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### Enamel classification for VITA SYSTEM 3D-MASTER®

1M1	1M2	2L1.5	2L2.5	2M1	2M2	2M3	2R1.5	2R2.5	3L1.5	3L2.5	3M1	3M2
EN2	EN2	EN2	EN2	EN2	EN2	EN2	EN2	EN2	EN2	EN2	EN2	EN2

3M3	3R1.5	3R2.5	4L1.5	4L2.5	4M1	4M2	4M3	4R1.5	4R2.5	5M1	5M2	5M3
EN2	EN2	EN2	EN4	EN4	EN4	EN4	EN4	EN4	EN4	EN4	EN4	EN4

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### Enamel classification for VITA classical

A1	A2	A3	A3,5	A4	B1	B2	B3	B4	C1	C2	C3	C4	D2	D3	D4
EN2	EN2	EN2	EN4	EN4	EN1	EN1	EN2	EN2	EN2	EN2	EN3	EN2	EN2	EN4	EN4

## Opakerzuordnung VITA SYSTEM 3D-MASTER®

### Opaque classification for VITA SYSTEM 3D-MASTER®

1M1	1M2	2L1.5	2L2.5	2M1	2M2	2M3	2R1.5	2R2.5	3L1.5	3L2.5	3M1	3M2
OP1	OP1	OP2	OP2	OP2	OP2	OP2	OP2	OP2	OP3	OP3	OP3	OP3

3M3	3R1.5	3R2.5	4L1.5	4L2.5	4M1	4M2	4M3	4R1.5	4R2.5	5M1	5M2	5M3
OP3	OP3	OP3	OP4	OP4	OP4	OP4	OP4	OP4	OP4	OP5	OP5	OP5

## Marginzuordnung VITA SYSTEM 3D-MASTER®

### Margin classification for VITA SYSTEM 3D-MASTER®

1M1	1M2	2L1.5	2L2.5	2M1	2M2	2M3	2R1.5	2R2.5	3L1.5	3L2.5	3M1	3M2
50% MAR1 + 50% MARN	MAR1	50% MAR2 + 50% MARN	MAR2	50% MAR2 + 50% MARN	MAR2	MAR2	50% MAR2 + 50% MARN	MAR2	50% MAR3 + 50% MARN	MAR3	50% MAR3 + 50% MARN	MAR3

3M3	3R1.5	3R2.5	4L1.5	4L2.5	4M1	4M2	4M3	4R1.5	4R2.5	5M1	5M2	5M3
MAR3	50% MAR3 + 50% MARN	MAR3	50% MAR4 + 50% MARN	MAR4	50% MAR4 + 50% MARN	MAR4	MAR4	50% MAR4 + 50% MARN	MAR4	50% MAR5 + 50% MARN	MAR5	MAR5

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### Margin classification for VITA classical

A1	A2	A3	A3,5	A4	B1	B2	B3	B4	C1	C2	C3	C4	D2	D3	D4
MAR2	MAR3	MAR3	MAR4	MAR5	MAR1	MAR2	MAR3	MAR4	25% MAR6 + 75% MARN	50% MAR6 + 50% MARN	MAR6	MAR6	MAR4	MAR4	50% MAR6 + 50% MARN

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*VITA TITANKERAMIK porcelains consist of a mixture of modified feldspar frits and minimal quantities of colour frits melted in to the shade-determining metal oxides. There is no known danger for the user if used according to manufacturer's instructions. In case of dust formation, it is necessary to use extraction unit or dust mask P2 for example (or grind while wet). Protective goggles should be worn when grinding the fired ceramic.*

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
The following products are subject to obligatory labeling:		
<b>VITA OPAQUE FLUID</b>	<p><b>Corrosive</b></p> <p>Causes severe burns. Store under lock and key and out of the reach of children. Do not eat and drink while working. In case of eye contact rinse immediately and thoroughly with water and consult physician. Do not allow to penetrate the sewage system; this product and its container must be disposed of according to the regulations for hazardous waste. Wear appropriate protective clothing when working. Wear protective gloves and protective goggles/face mask. In case of accident or unwellness consult physician immediately (if possible, show this label).</p>	
<b>VITA SPRAY-ON INDICATOR LIQUID und VITA SPRAY-ON LIQUID</b>	<p><b>Flammable</b></p> <p>Keep container tightly closed and store in a well-ventilated place. Keep away from sources of ignition – do not smoke. Do not allow to enter into the sewage system. This product and its container must be disposed of according to the regulations for hazardous waste.</p>	 

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*TITANKERAMIK is available in VITA SYSTEM 3D-MASTER and VITA classical A1 – D4 shades. Shade compatibility with all VITA SYSTEM 3D-MASTER and VITA classical A1 – D4 materials is ensured.*

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**VITA TITANKERAMIK, VITA AKZENT® Plus, VITA INTERNO®**

**Please note:** Our products must be used in accordance with the instructions for use. We accept no liability for any damage resulting from incorrect handling or usage. The user is furthermore obliged to check the product before use with regard to its suitability for the intended area of applications. We cannot accept any liability if the product is used in conjunction with materials and equipment from other manufacturers that are not compatible or not authorized for use with our product. Furthermore, our liability for the accuracy of this information is independent of the legal basis and, in as far as legally permissible, shall always be limited to the value as invoiced of the goods supplied, excluding value-added tax. In particular, as far as legally permissible, we do not assume any liability for loss of earnings, indirect damages, ensuing damages or for third-party claims against the purchaser. Claims for damages based on fault liability (culpa in contrahendo, breach of contract, unlawful acts, etc.) can only be made in the case of intent or gross negligence. The VITA Modulbox is not necessarily a component of the product. Date of issue of these working instructions: 07.17 After the publication of these working instructions any previous versions become obsolete. The current version can be found at [www.vita-zahnfabrik.com](http://www.vita-zahnfabrik.com)

The following materials are certified and bear the CE mark:  0124

**VITA TITANKERAMIK, VITA AKZENT® Plus, VITA INTERNO®**

# VITA

VITA Zahnfabrik H. Rauter GmbH & Co. KG  
 Spitalgasse 3 · D-79713 Bad Säckingen · Germany  
 Tel. +49(0)7761/562-0 · Fax +49(0)7761/562-299  
 Hotline: Tel. +49(0)7761/562-222 · Fax +49(0)7761/562-446  
[www.vita-zahnfabrik.com](http://www.vita-zahnfabrik.com) · [info@vita-zahnfabrik.com](mailto:info@vita-zahnfabrik.com)  
 [facebook.com/vita.zahnfabrik](https://www.facebook.com/vita.zahnfabrik)