

VITA SMART.FIRE

Instructions for performing a software update



Dear VITA SMART.FIRE customer,

We would like to inform you about the new software update for your VITA SMART.FIRE and ask you to install this update on your equipment, according to the following instructions:

if you have any questions or comments, please do not hesitate to contact us using the following details:




Tel.: 0049 77 61 / 562 – 101, -105 or -106

Fax: 0049 77 61 / 562 – 102

Email: instruments-service@vita-zahnfabrik.com

Best regards,
VITA Zahnfabrik
H. Rauter GmbH & Co. KG
Equipment servicing department

<p>1</p> <p>Export the update file onto an empty USB flash drive</p>	
<p>2</p> <p>Insert the USB flash drive into one of the USB ports</p>	
<p>3</p> <p>Press settings</p>	

<p>4 Press the update button</p>	 <p>A screenshot of a desktop environment with a light blue background. In the center, there are four application icons: a laptop, a pink icon with a gear, a blue icon with a gear and a white 'i', and a green icon with a white circular arrow. The green icon is circled in red. In the bottom right corner, there is a grey button with a white curved arrow pointing left.</p>
<p>5 You will be notified when the update is complete.</p>	 <p>A screenshot of a desktop environment titled 'VITA SMART.FIRE'. The desktop has a light blue background and several application icons. A grey dialog box is centered on the screen with the text 'Update completed' and a white checkmark icon. The icons on the desktop include 'VIT Stain', 'SUPRINITY PC Kombi/PASTE', 'SUPRINITY PC Stains-fixation', 'SUPRINITY PC GLAZE LT SPRAY', 'SUPRINITY PC GLAZE LT PASTE', 'YZ Stains-fixation', 'YZ GLAZE LT SPRAY', and 'YZ GLAZE LT PASTE'. A 'STANDARD' button is visible in the bottom right corner.</p>
<p>6 Now you will be asked if you would like to replace the current "STANDARD desktop" with the updated "STANDARD desktop"</p> <p><i>Note: If the current "STANDARD desktop" is replaced, then all of the changes you have made to this desktop will be lost.</i></p>	 <p>A screenshot of a desktop environment titled 'VITA SMART.FIRE'. The desktop has a light blue background and several application icons. A grey dialog box is centered on the screen with the text 'Desktop 'STANDARD' replace by the updated Desktop?' and two buttons: a white 'X' and a white checkmark. The icons on the desktop are the same as in the previous screenshot. A 'STANDARD' button is visible in the bottom right corner.</p>

Then you will be asked if you would like to replace current "STANDARD 2 desktop" with the updated "STANDARD 2 desktop"

7

Note: If the current "STANDARD 2 desktop" is replaced, then all of the changes you have made to this desktop will be lost

