VITA CERAMICS ETCH

DE

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VITA CERAMICS ETCH

EN

VITA CERAMICS ETCH is a hydrofluoric acid gel for etching ceramic restorations. It is used to create a retentive etched pattern on restorations made of silicate ceramic and hybrid ceramic in order to strengthen the adhesion between the ceramic surface and the adhesive composite.

- easy and precise application due to gel consistency
- easily distinguishable by its red colour

Composition

1 ml VITA CERAMICS ETCH contains 0.047 g hydrofluoric acid (HF)

Area of application

Silicate ceramic restorations (glass ceramic, feldspar ceramic) and hybrid ceramic (VITA ENAMIC).

Contraindications

VITA CERAMICS ETCH must not be used in the mouth or for etching tooth enamel.

Oxide ceramic (e.g. zirconium oxide), since restorations made of these ceramics are cemented using a modified adhesive technique without ceramic etching.

Safety advices

For extraoral use only! Contains hydrofluoric acid. Toxic if swallowed. Fatal in contact with skin. Causes severe skin burns and eye damage. Harmful if inhaled. Wear safety goggles / protective gloves / protective clothing. Store in a secure location. If swallowed, call the poisons information center immediately and refer to

the safety data sheet. In case of contact with clothing / skin, remove contaminated items of clothing immediately and rinse with large amounts of water. For specific measures, see the safety data sheet. In case of contact with the eyes, rinse for several minutes with water and consult a doctor / poison information center. This material and its container must be disposed of as hazardous waste.

Application

Apply VITA CERAMICS ETCH to the surfaces requiring etching with a small disposable brush or a plastic disposable applying needle.

When using the syringe with the application tip:

Screw a red application tip tightly onto the VITA CERAMICS ETCH syringe. The tip enables controlled extrusion from the syringe.

Prior to application on the ceramic restoration, first apply a small amount of VITA CERAMICS ETCH to the mixing pad to ensure even application and in order to become familiar with the flow behavior.

In the event of strong resistance to application, do not continue. Instead, attach a new tip and test again.

Dispose of used application tips immediately in accordance with statutory regulations in order to prevent contact.

Etching time

- Feldspar ceramic (e.g. VITABLOCS Mark II) 60 s
- Hybrid ceramic (VITA ENAMIC) 60 s
- Zirconia reinforced lithium silicate ceramic (VITA SUPRINITY PC) 20 s

Completely remove any acid remaining on the etched surface by intensive spraying with a lot of water or steam or by cleaning in an ultrasonic bath in distilled water. Do not brush off, since this would lead to severe imperfections on the surface!

Dry well before inserting, use silane coupling agent (e.g. VITA ADIVA C-PRIME) according to manufacturer's instructions.

Important

The ceramic restoration should be etched immediately before insertion or after the try-in, in order to avoid contamination of the etched ceramic surface. The etched surface must no longer be touched!

Storage

Store at temperatures between 2° and 28 °C. Protect from direct sunlight.

Waste disposal

The empty packing should be handed over to hazardous waste disposers.

Labelling

The expiry date (📅) and the batch no. (🏷️) are clearly marked on the containers.

Presentation

- VITA CERAMICS ETCH, dropper bottle, 6 ml
- VITA CERAMICS ETCH, disposable syringe, 3 ml, 20 plastic application needles

Toxic

Corrosive

VITA CERAMICS ETCH contains < 5 Vol. % hydrofluoric acid (HF)

VITA CERAMICS ETCH contains < 10 Vol. % sulphuric acid (H₂SO₄)

The safety datasheet can be downloaded at www.vita-zahnfabrik.com or requested by sending a fax to (+49) 7761-562-233.

The product was developed for dental use and must be processed according to the instructions for use. The information contained herein do not constitute an assurance of properties and are not binding. These instructions issued 01.2018.

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VITA Zahnfabrik H. Rauter GmbH & Co. KG