# VITA YZ<sup>®</sup> SOLUTIONS

---

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VITA YZ SOLUTIONS discs and blocks are blanks used for fabricating dental restorations.

**1. Indication:**

VITA YZ T and VITA YZ HT (type II, class 5; CTE 10.5x10<sup>-6</sup>/K, according to DIN EN ISO 6872): Fully anatomical crowns and bridges with up to 14 units\* in the anterior and posterior area; fully and partially veneered single tooth and bridge substructures\* with up to 14 units in the anterior and posterior area; single tooth and bridges\* with up to 14 units on direct screw-retained restorations in the anterior and posterior area; primary telescopes

VITA YZ ST\*\* (type II, class 5; CTE 10.3x10<sup>-6</sup>/K; according to DIN EN ISO 6872):

Fully anatomical crowns and bridges\* with up to 14 units\*\* in the anterior and posterior area; fully and partially veneered single tooth and bridge substructures\* with up to 14 units\*\* in the anterior and posterior area; single tooth and bridges\* with up to 14 units\*\* on direct screw-retained restorations in the anterior and posterior area; inlays, onlays, veneers, partial crowns

VITA YZ XT (type II, class 4a; WAK 10.0x10<sup>-6</sup>/K, according to DIN EN ISO 6872):

Fully anatomical single crowns and bridges with up to three units with one bridge pontic; fully and partially veneered single crown and bridge substructures with up to three units with one bridge pontic in the anterior and posterior area; inlays, onlays, veneers, partial crowns

\* Bridges and bridge superstructures with maximum of two adjoining bridge pontics

\*\* VITA YZ ST is limited in Canada to bridge indications with a max. of six units with a max. of two adjoining bridge pontics.

**2. Contraindication:**

VITA YZ T / VITA YZ HT / VITA YZ ST: for more than two adjacent bridge units; for two or more cantilever bridge units; parafunctions for veneered restorations, in particular for bruxism; for insufficient oral hygiene; for insufficient preparation results; for insufficient hard dental tissue; for patients suffering from allergies or sensitivities to the constituents; temporary bonding of VITA YZ ST

VITA YZ XT: for bridge restorations with more than three units; cantilever bridges; parafunctions for veneered restorations, in particular for bruxism; for insufficient oral hygiene; for insufficient preparation results; for insufficient hard dental tissue; for patients suffering from allergies or sensitivities to the constituents; for temporary bonding

**3. Design parameters**

All information refers to sintered restorations	Minimum wall thickness in mm	Connector cross-sections in mm <sup>2</sup>							
		VITA YZ T	VITA YZ HT	VITA YZ ST**	VITA YZ XT*	VITA YZ T	VITA YZ HT	VITA YZ ST	VITA YZ XT
Inlay, onlay, veneer	incisal	0.5	0.5	0.6	0.8	–	–	–	–
	occlusal	0.5	0.5	0.6	0.8	–	–	–	–
	circumferential	0.4	0.4	0.5	0.7	–	–	–	–
Anterior and posterior crowns (fully anatomical or substructure)	incisal	0.5	0.5	0.6	0.8	–	–	–	–
	occlusal	0.5	0.5	0.6	0.8	–	–	–	–
	circumferential	0.4	0.4	0.5	0.7	–	–	–	–
Fully anatomical anterior bridges and substructures with one pontic*	incisal	0.5	0.5	0.6	1.0	7	7	9	9
	circumferential	0.5	0.5	0.6	0.8	–	–	–	–
Fully anatomical posterior bridges and substructures with one pontic*	occlusal	0.6	0.6	0.7	1.2	9	9	12	12
	circumferential	0.5	0.5	0.6	1.0	–	–	–	–
Fully anatomical multi-unit anterior bridges and substructures with two pontics	incisal	0.6	0.6	0.8	–	9	9	12	–
	circumferential	0.5	0.5	0.6	–	–	–	–	–
Fully anatomical multi-unit posterior bridges and substructures with two pontics	occlusal	0.7	0.7	0.8	–	12	12	15	–
	circumferential	0.6	0.6	0.6	–	–	–	–	–
Cantilever bridges with an extension	incisal	0.7	0.7	0.8	–	12	12	15	–
	occlusal	0.7	0.7	0.8	–	–	–	–	–
	circumferential	0.5	0.5	0.6	–	–	–	–	–

\*) VITA YZ XT is limited to maximum of three units.

\*\*) In Canada VITA YZ ST is limited to bridges with six units.

**Please note:**

When designing the framework, pay attention to anatomically reduced tooth shapes. The cusp should be supported according to the anatomical progression. When preparing a ceramic veneer, a uniform layer thickness across the entire surface to be veneered must be ensured. The entire thickness of the ceramic layer, however, should not exceed 2 mm (the optimal layer thickness ranges from 0.7 to 1.2 mm).

When designing fully anatomical restorations, the minimum wall thicknesses for the respective material must be observed. The aim is to obtain a uniform material thickness. Clinical success cannot be reliably assured in the event that compliance with the given design guidelines and minimum wall thickness is not met.

**4. CAD process**

For seven units or more, a sinter support for the milled restoration should be designed for bridge frameworks from VITA YZ T / HT / ST. If no sinter support is attached to these bridge structures, it is possible for distortion to occur as a result of the sintering process, which may in turn lead to an imprecise fit in the finished restoration.

The instructions on the design of the sintering support from the individual software manufacturers are to be observed.

However, if this sinter support consists only of a leftover piece of the disc (sintering tongue), the volume should be reduced to ensure even heating during the sintering process.

**5. Enlargement factor**

All VITA YZ materials are porous before sintering. Sintering acts as a kind of "compression process," which results in the materials shrinking by approximately 20%. As a result, all restorations must be milled according to the specified enlargement factor to compensate for the shrinkage of the material during the sintering process. The corresponding blank imprint must be entered into the respective CAM software (as plain text on the disc or as a barcode on the block). Some software requires the factor (e.g., enlargement factor 1.2345) for the discs or the corresponding X, Y values (e.g., VGF: X, Y = 23.45) and for some providers, also the Z values (e.g., Z = 23.01). For blocks, the corresponding barcode is scanned into the CAM software.

**6. Milling the restorations**

VITA YZ processing must be performed with a validated CAD/CAM system. For information on the processing parameters, see

"Processing information for VITA YZ CNC machines." Further information is available on the internet at: <https://www.vita-zahnfabrik.com/systempartner>

- VITA SYSTEM SOLUTIONS with specific holder systems for: inLab System (Sirona Dental Systems GmbH); CEREC System (Sirona Dental Systems GmbH)
- VITA UNIVERSAL SOLUTIONS in universal disc geometries (Ø98.4 mm) for the following CAD/CAM systems: CORITEC series (imes-icore GmbH); ULTRASONIC series (DMG Mori AG); RXD line (Röders GmbH); inLab MC X5 (Sirona Dental Systems GmbH); S1/S2/K3/K4/K5 Impression (vhf manufacture AG)

Please note: restorations from VITA YZ ST and VITA YZ XT are not to be milled wet. They must be milled dry.

## 7. Reworking the ground/milled restoration

Restorations without sinter support: After completing the milling or grinding process, bridge restorations without sinter support as well as single crowns, are detached with a diamond burr milling cutter or a cross-cut tungsten carbide bur. Thickly milled margins must be reduced carefully. Fissures can be adjusted using a thin pointed tool.

Bridge restorations may not be separated using a diamond separating disc, since this may result in predetermined breaking points on the connectors. Restorations with sinter support: Bridges with eight or more units are not separated from the sintering frame or sintering tongue before the sintering process. All milling fibers must be removed before sintering.

During wet processing for restorations from VITA YZ T and VITA YZ HT (wet grinding), a cleaning firing should be carried out before sintering.

	Predry. °C	→ min.	↗ min.	↗ °C/ min.	T °C	→ min.	Vac. min.
YZ T cleaning firing	500	3.00	6.00	33	700	5.00	–
YZ HT cleaning firing	290	10.00	31.00	10	600	5.00	–

## 8. Characterizing frameworks and restorations

There are two options for characterizing VITA YZ SOLUTIONS: the brushing technique and the dipping technique.

Characterized frameworks must be dry before sintering. An infrared lamp (250 watts) or the corresponding pre-dry program in the VITA ZYRCOMAT 6000 MS can be used for drying.

For more information, refer to Instruction manual no. 10116 or [www.vita-zahnfabrik.com/cadcam](http://www.vita-zahnfabrik.com/cadcam)

## 9. Sintering in high-temperature furnaces

VITA YZ materials can be sintered in all high-temperature furnaces, which can be operated with the sintering parameters indicated above. The corresponding information of the furnace manufacturer must be observed. VITA does not grant a warranty or accept any liability for damage resulting from processing VITA YZ SOLUTIONS materials in furnaces from other manufacturers.

Storage during the sintering process: In order to ensure linear shrinkage during the sintering process and to ensure a high accuracy of fit for the VITA YZ materials, the restoration must be kept under optimal storage conditions. For storage on a bed of sintering balls, place anterior crowns and bridges on the labial or lingual surfaces, and place posterior crowns or bridges on the occlusal surface. Restorations with sinter support are placed centrally on a smooth, flat surface, either on an appropriate sinter support or directly on the VITA ZYRCOMAT 6000 MS firing base.

When using liquids, any contamination of the sintering chamber should be avoided. Cover the sintered objects for this purpose. Restorations made of VITA YZ T and VITA YZ HT can also be sintered in the high-speed mode. For further information, please refer to the VITA YZ SOLUTIONS Working Instructions.

Program name	Lift position %	T0 °C	↗ min.	↗ °C/ min.	T1 °C	→ min.	↘ °C	Lift position %
Pre-dry	50	25	7:21	17	150	30:00	–	50
YZ T Universal	100	25	88:32	17	1530	120:00	200	100
YZ HT Universal	100	25	83:49	17	1450	120:00	200	100
YZ ST Universal	100	25	188:08	8	1530	120:00	200	100
YZ XT Universal	100	25	356:15	4	1450	120:00	200	100

## 10. Reworking

The surface structure of ceramic materials is important for their load capacity. That is why reworking the sintered framework should generally be avoided.

However, if any reworking is required apart from grinding the holding pieces, the following basic factors need to be observed:

- Corrections need to be made with diamond instruments for wet grinding under water cooling and at a low pressure.
- Only fine grain diamonds in mint condition with red color coding (fine 27–76 µm) or lower (extra-fine, yellow 10–36 µm or ultrafine, white 4–14 µm) should be used
- When reworking, the minimum wall thicknesses and connector cross-sections should be observed, and notches and sharp edges should be avoided.
- Any reworking, especially near the bridge connectors, must be strictly avoided.

For further information, please refer to the VITA YZ SOLUTIONS Working Instructions.

## 11. Processing options

- Polishing: VITA SUPRINITY Polishing Set
  - Staining technique: VITA AKZENT Plus Stains for characterizing the shade and glazing of fully anatomical VITA YZ restorations.
  - Veneering: VITA VM 9 for individualizing reduced crowns and bridges.
  - Press-over technique: VITA PM 9 for pressing onto crown and bridge substructures.
- For further information on the application and on the corresponding firing temperatures, please refer to the respective working instructions of the material.

## 12. Bonding

VITA YZ SOLUTIONS can be luted with conventional, self-adhesive or adhesive methods. For more information on bonding, please refer to the manufacturer's information of the bonding material.

Adhesive bonding is required for inlays, onlays and veneers.

Note: Dental treatment and the integration of dental restorations entail the general risk of iatrogenic damage to hard tooth substance, pulp and/or oral soft tissue. The use of bonding systems and the integration of dental restorations involve the general risk of postoperative hypersensitivity.

In the event of non-compliance with the processing instructions of the products in use, the product characteristics can not be ensured so that product failure and irreversible damage to the natural hard tooth substance, pulp and/or oral soft tissue may result.

## 13. Clinical risks

Other, less likely, but later possible complications and risks in dental treatment are well known and are not associated with VITA materials. Patients should be informed by the dentist of the risk of later possible problems, such as fracture, chipping, debonding, occlusal roughness, crevices, over-contouring, marginal discrepancy (marginal gap), secondary caries, infection or other endodontic or periodontal problems.

## 14. General information

Occupational and health safety: When working with the product, wear suitable safety goggles/face protection and light respiratory protection.

## 15. Storage

VITA YZ materials should be stored in a dry and dark place and not exposed to direct sunlight.

## Please note

Our products must be used in accordance with the instructions for use.

We accept no liability for any damage resulting from incorrect handling or usage. The user is furthermore obliged to check the product before use with regard to its suitability for the intended area of applications. We cannot accept any liability if the product is used in conjunction with materials and equipment from other manufacturers that are not compatible or not authorized for use with our product and this results in damage. The VITA Modulbox is not necessarily a component of the product. Date of issue of this information: 10.18

After the publication of this information for use any previous versions become obsolete. The current version can be found at [www.vita-zahnfabrik.com](http://www.vita-zahnfabrik.com)

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VITA Zahnfabrik · H. Rauter GmbH & Co. KG