**VITA PM Einbettmasse**  
**VITA PM Investment**  
**VITA PM Revêtement**  
**VITA PM Revestimiento**  
**VITA PM Rivestimento**

Gebrauchsanweisung  
Instructions for use  
Mode d'emploi  
Instrucciones de uso  
Istruzioni per l'uso

Date of issue: 10.19

**VITA – perfect match.**

**VITA**

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
## Physical data\*

\* in accordance with DIN EN ISO 9694, manufacturer data

Test parameter	VITA PM Investment	
Beginning of setting	[min]	approx. 9
Compressive strength	[MPa]	approx. 4.0
Flowability	[mm]	130 – 140
Concentration of liquid	[%]	75

## Safety information:

Caution!  
Contains quartz  
and cristobalite.

The following needs to be observed during processing:

- Causes damage to organs after prolonged or repeated exposure.
- Use a local dust extraction device.
- In case of insufficient ventilation, wear suitable respiratory equipment.
- Do not inhale dust.
- May cause mechanical irritation of the skin and mucosa of the eye and respiratory tract.
- Do not mix the powder with other materials.
- Dispose of the contents/container in accordance with local regulations.

The following needs to be observed when processing the liquid:

- Caution! Mixing liquids are alkaline!
- If mixing liquid comes into contact with skin, rinse IMMEDIATELY with copious amount of water!
- If mixing liquid comes into contact with the eyes: rinse with copious amount of water while keeping the eyelid open. In case of concerns: consult an ophthalmologist.
- Do not mix the liquid with other liquids.
- Protect the liquid from freezing!

When heating above 200 °C, the following needs to be observed:

- Formation of ammonia vapors.
- In the case of local contact, irritation and possibly caustic effect on the skin and mucosa (eyes, respiratory tract and gastrointestinal tract after swallowing) may occur.

Do not open the furnace during preheating as this would cause aggressive combustion. The preheat furnace should only be filled to half of its maximum capacity.

This investment material contains quartz powder! Please take suitable measures to avoid inhalation, such as wearing a dust respirator.

The figures given by our research and development department are guidelines and can be individually adjusted.

Should you have any further questions, please call our product specialists  
Phone: ++49/7761/562-222.

The product was developed for dental use and must be processed according to the instructions for use. The information contained herein do not constitute an assurance of properties and are not binding. These instructions issued: 10.19


Graphit-free, phosphate-bonded precision investing material for “speed” preheating especially for VITA PM press ceramic

	<b>Storage</b>	<b>keep the powder in cool, dry place</b>	<b>Do not store the mixing liquid at less than 5°C (41°F), as it is sensitive to frost</b>		
	Working temperature	Approx. 22°C (71.6°F) (room temperature)			
Mixing ratio	100 g Powder – 22 ml Liquid	<b>Mixing ratio for sachets:</b>			
		Powder (g)	Liquid (ml)		
		100	22		
		200	44		
Liquid concentration (examples)	<b>Liquid concentration for 1000 ml liquid</b>				
	Concentration	Liquid (ml)	Distilled Water (ml)		
	90 %	900	100		
	80 %	800	200		
	70 %	700	300		
	60 %	600	400		
	50 %	500	500		
	40 %	400	600		
Concentrations		Sachet 1 x 100 g		Sachet 2 x 100 g (200 g)	
	Concentration (%)	Liquid (ml)	Dist. Water (ml)	Liquid (ml)	Dist. Water (ml)
	100	22	–	44	–
	95	21	1	42	2
	90	20	2	40	4
	85	19	3	38	6
	80	18	4	36	8
	75	17	5	34	10
	70	15	7	30	14
	65	14	8	28	16
	60	13	9	26	18
	55	12	10	24	20
50	11	11	22	22	
45	10	12	20	24	
40	9	13	18	26	

Expansion control	<b>Expansion control (%)</b>		
		Liquid (ml)	Dist. Water (ml)
	<b>Overpressing technique (crowns, bridges):</b>	75 %	25 %
	<b>Substructure-free staining and layering technique:</b>		
	Anterior and premolar crowns	75 - 80 %	25 - 20 %
	Molar crowns	75 %	25 %
	Veneers	70 %	30 %
	Inlays 1-layer and 2-layer	50 - 60 %	50 - 40 %
	MOD-Inlays	75 %	25 %
	Onlays	85 - max. 90 %	15 - 10 %

For anterior and premolar crowns attention has to be paid to the fact that better fitting can be obtained with slim and lean preparations at the higher expansion rate of 80%.  
Expansion figures given above are guidelines and have to be adjusted accordingly i.e. due to different preparation models, press pellets and press temperature.

Stirring		Use a spatula to stir the investment material by hand until the powder has been wetted thoroughly.
Mixing in vacuum	60 seconds	Place in vacuum for 15 seconds, but without switching the paddles on. Mix for 60 seconds. Check the vacuum/paddles frequently for correct functioning. Inadequate vacuum results in inaccurate fit and bubbles on the casting.
Working time	Approx. 6 minutes at approx. 22°C (71.6°F) (room temperature)	The working time depends on the room temperature. Warmth curtails the working time.
Investing		Fill the ring with investment: only use a vibrator if it is required to improve the flow. Avoid high vibration! This will lead to the formation of bubbles and breakdown of the mixture.
Pre-drying time	30 minutes after start of mix 1. Remove level gauge and sprue base after 20 minutes 2. Leave muffle to cool for 10 minutes	
Placing the ring		1. Straighten bottom of ring (plaster knife/sand paper) 2. After 30 minutes, place the ring in the 850°C preheated furnace
Preheating	850°C	
Hold-times		Hold-time of ring once preheating temperature (850°C) is attained again 100 g ring: at least 50 minutes 200 g ring: at least 75 minutes From three 100 g or 200 g muffles in the preheating furnace, the hold time must be increased by 15 minutes in each case.
Investing		Investment to be blasted with 50 µm glass beads only. Before reaching pressed objects, remove remaining investment material with a pressure of <b>max. 2 bar</b> from same.

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