# VITA AMBRIA® Calibration of the press furnace

## ➤ Setting the optimal press temperature

We recommend investing in a wax grid together with the restoration for the first pressing.

After that, you can use these pictures to evaluate your press results and adjust the press temperature accordingly.

- 20 °C

The grid has not flowed out, the restoration looks opaque.  
The pressing temperature is far too low.  
Procedure: Raise the temperature in 10°C-steps.

- 15 °C

The grid has a smooth surface, but parts of the grid have not completely flowed out.  
The pressing temperature is still too low.  
Procedure: Raise the temperature in 5°C-steps.

- 5 °C

The grid has flowed out almost completely.  
The pressing temperature is minimally low.  
Procedure: Raise the temperature slightly (approx. 5°C).

+/- 0 °C

The grid has completely flowed out, the surface is smooth and shows no reaction layer.  
The pressing temperature is optimal.

+ 5 °C

The grid has flowed out, the surface shows a light reaction layer.  
The pressing temperature is a bit too high.  
Procedure: Minimally lower the temperature (approx. 5°C).

+ 15 °C

The grid has flowed out, the surface is whitish and partly also porous.  
The reaction layer is clearly visible.  
The press temperature is too high.  
Procedure: Lower the temperature in 5°C-steps.

+ 20 °C

The reaction layer is very pronounced, there are small cracks and holes in the edge area.  
The restoration can appear transparent / colorless.  
The press temperature is too high.  
Procedure: Lower the temperature in 10°C-steps.

Find more languages online:  
[www.vita-zahnfabrik.com/ambria](http://www.vita-zahnfabrik.com/ambria)

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