

## Safety Data Sheet

according to UN GHS (ST/SG/AC.10/11/Rev.9)

### VITA VIONIC BOND I

Revision date: 17.01.2024

Product code: 288

Page 1 of 8

## SECTION 1: Identification of the substance/mixture and of the company/undertaking

### 1.1. Product identifier

VITA VIONIC BOND I

### 1.2. Relevant identified uses of the substance or mixture and uses advised against

#### Use of the substance/mixture

Use as laboratory reagent

### 1.3. Details of the supplier of the safety data sheet

|                         |                                       |                             |
|-------------------------|---------------------------------------|-----------------------------|
| Company name:           | VITA Zahnfabrik H.Rauter GmbH & Co.KG |                             |
| Street:                 | Spitalgasse 3                         |                             |
| Place:                  | D-79713 Bad Säckingen                 |                             |
| Post-office box:        | 1338                                  |                             |
|                         | D-79704 Bad Säckingen                 |                             |
| Telephone:              | +49(0)7761-562-0                      | Telefax: +49(0)7761-562-299 |
| E-mail:                 | info@vita-zahnfabrik.com              |                             |
| Contact person:         | regulatory affairs                    |                             |
| E-mail:                 | info@vita-zahnfabrik.com              |                             |
| Internet:               | www.vita-zahnfabrik.com               |                             |
| Responsible Department: | Regulatory Affairs                    |                             |

### Further Information

medical device

## SECTION 2: Hazards identification

### 2.1. Classification of the substance or mixture

#### UN GHS (ST/SG/AC.10/11/Rev.9)

Flammable liquids: Flam. Liq. 2  
Skin corrosion/irritation: Skin Irrit. 2  
Respiratory or skin sensitisation: Skin Sens. 1  
Specific target organ toxicity - single exposure: STOT SE 3 (respiratory tract irritation)

### 2.2. Label elements

#### UN GHS (ST/SG/AC.10/11/Rev.9)

#### Hazard components for labelling

2,2'-ethylenedioxydiethyl dimethacrylate  
dibenzoyl peroxide; benzoyl peroxide  
methyl methacrylate

Signal word: Danger

Pictograms:



#### Hazard statements

H225 Highly flammable liquid and vapour.  
H315 Causes skin irritation.  
H317 May cause an allergic skin reaction.  
H335 May cause respiratory irritation.

#### Precautionary statements

P403+P233 Store in a well-ventilated place. Keep container tightly closed.

## Safety Data Sheet

according to UN GHS (ST/SG/AC.10/11/Rev.9)

### VITA VIONIC BOND I

Revision date: 17.01.2024

Product code: 288

Page 2 of 8

|      |  |
|------|--|
| P235 | Keep cool.   |
| P233 | Keep container tightly closed.   |
| P210 | Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking. |
| P261 | Avoid breathing dust/fume/gas/mist/vapours/spray.  |
| P280 | Wear protective gloves/protective clothing/eye protection/face protection/hearing protection.  |

#### 2.3. Other hazards

No information available.

### SECTION 3: Composition/information on ingredients

#### 3.2. Mixtures

##### Hazardous components

| CAS No   | Chemical name   | Quantity    |
|----------|---|-------------|
|          | Classification (UN GHS (ST/SG/AC.10/11/Rev.9))                            |             |
| 80-62-6  | methyl methacrylate   | 75 - < 80 % |
|          | Flam. Liq. 2, Skin Irrit. 2, Skin Sens. 1, STOT SE 3; H225 H315 H317 H335 |             |
| 109-16-0 | triethylene glycol dimethacrylate   | 5 - < 10 %  |
|          | Skin Irrit. 2, Eye Irrit. 2, Skin Sens. 1, STOT SE 3; H315 H319 H317 H335 |             |
| 94-36-0  | dibenzoyl peroxide; benzoyl peroxide                                      | 1 - < 5 %   |
|          | Org. Perox. B, Eye Irrit. 2A, Skin Sens. 1; H241 H319 H317                |             |

### SECTION 4: First aid measures

#### 4.1. Description of first aid measures

##### After inhalation

Provide fresh air. Medical treatment necessary.

##### After contact with skin

After contact with skin, wash immediately with plenty of water and soap.

##### After contact with eyes

After contact with the eyes, rinse with water with the eyelids open for a sufficient length of time, then consult an ophthalmologist immediately.

##### After ingestion

Rinse mouth immediately and drink plenty of water.

#### 4.2. Most important symptoms and effects, both acute and delayed

No information available.

#### 4.3. Indication of any immediate medical attention and special treatment needed

Treat symptomatically.

### SECTION 5: Firefighting measures

#### 5.1. Extinguishing media

##### Suitable extinguishing media

Carbon dioxide (CO<sub>2</sub>), Foam, Extinguishing powder.

##### Unsuitable extinguishing media

Water.

#### 5.2. Special hazards arising from the substance or mixture

Highly flammable. Vapours can form explosive mixtures with air.

#### 5.3. Advice for firefighters

## Safety Data Sheet

according to UN GHS (ST/SG/AC.10/11/Rev.9)

### VITA VIONIC BOND I

Revision date: 17.01.2024

Product code: 288

Page 3 of 8

Wear a self-contained breathing apparatus and chemical protective clothing. Full protection suit.

#### Additional information

Use water spray jet to protect personnel and to cool endangered containers. Suppress gases/vapours/mists with water spray jet. Collect contaminated fire extinguishing water separately. Do not allow entering drains or surface water.

### SECTION 6: Accidental release measures

#### 6.1. Personal precautions, protective equipment and emergency procedures

##### General advice

Remove all sources of ignition. Do not breathe gas/fumes/vapour/spray. Avoid contact with skin, eyes and clothes. Use personal protection equipment.

#### 6.2. Environmental precautions

Do not allow uncontrolled discharge of product into the environment. Danger of explosion

#### 6.3. Methods and material for containment and cleaning up

##### Other information

Absorb with liquid-binding material (sand, diatomaceous earth, acid- or universal binding agents). Treat the recovered material as prescribed in the section on waste disposal.

#### 6.4. Reference to other sections

Safe handling: see section 7

Personal protection equipment: see section 8

Disposal: see section 13

### SECTION 7: Handling and storage

#### 7.1. Precautions for safe handling

##### Advice on safe handling

If handled uncovered, arrangements with local exhaust ventilation have to be used. Do not breathe gas/fumes/vapour/spray.

##### Advice on protection against fire and explosion

Keep away from sources of ignition - No smoking. Take precautionary measures against static discharges. Vapours can form explosive mixtures with air.

##### Advice on general occupational hygiene

Remove contaminated, saturated clothing immediately. Draw up and observe skin protection programme. Wash hands and face before breaks and after work and take a shower if necessary. When using do not eat or drink.

#### 7.2. Conditions for safe storage, including any incompatibilities

##### Requirements for storage rooms and vessels

Keep container tightly closed. Keep in a cool, well-ventilated place. Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking.

##### Hints on joint storage

Do not store together with: Oxidising agent. Pyrophoric or self-heating substances.

### SECTION 8: Exposure controls/personal protection

#### 8.1. Control parameters

## Safety Data Sheet

according to UN GHS (ST/SG/AC.10/11/Rev.9)

### VITA VIONIC BOND I

Revision date: 17.01.2024

Product code: 288

Page 4 of 8

#### Exposure limits (EH40)

| CAS No   | Substance                  | ppm | mg/m <sup>3</sup> | fibres/ml | Category      | Origin |
|----------|----------------------------|-----|-------------------|-----------|---------------|--------|
| 128-37-0 | 2,6-Di-tert-butyl-p-cresol | -   | 10                |           | TWA (8 h)     | WEL    |
| 94-36-0  | Dibenzoyl peroxide         | -   | 5                 |           | TWA (8 h)     | WEL    |
| 80-62-6  | Methyl methacrylate        | 50  | 208               |           | TWA (8 h)     | WEL    |
|          |                            | 100 | 416               |           | STEL (15 min) | WEL    |

#### 8.2. Exposure controls



##### Appropriate engineering controls

If handled uncovered, arrangements with local exhaust ventilation have to be used. Do not breathe gas/fumes/vapour/spray.

##### Individual protection measures, such as personal protective equipment

##### Eye/face protection

Wear eye/face protection.

##### Hand protection

When handling with chemical substances, protective gloves must be worn with the CE-label including the four control digits. The quality of the protective gloves resistant to chemicals must be chosen as a function of the specific working place concentration and quantity of hazardous substances. For special purposes, it is recommended to check the resistance to chemicals of the protective gloves mentioned above together with the supplier of these gloves. Recommended glove articles KCL Butoject Butyl caoutchouc (butyl rubber)  
Breakthrough time: 60 min

##### Skin protection

Wear suitable protective clothing.

##### Respiratory protection

Technical ventilation of workplace Provide adequate ventilation as well as local exhaustion at critical locations.

### SECTION 9: Physical and chemical properties

#### 9.1. Information on basic physical and chemical properties

|   |                |
|---|----------------|
| Physical state:   | Liquid         |
| Colour:   | colourless     |
| Odour:  | stinging       |
| Melting point/freezing point:                             | not determined |
| Boiling point or initial boiling point and boiling range: | 101 °C         |
| Flammability:   | not applicable |
| Lower explosion limits:                                   | 2,1 vol. %     |
| Upper explosion limits:                                   | 12,5 vol. %    |
| Flash point:  | 10 °C          |
| Auto-ignition temperature:                                | 430 °C         |
| Decomposition temperature:                                | not determined |
| pH-Value:   | not determined |
| Water solubility:   | No             |

## Safety Data Sheet

according to UN GHS (ST/SG/AC.10/11/Rev.9)

### VITA VIONIC BOND I

Revision date: 17.01.2024

Product code: 288

Page 5 of 8

Solubility in other solvents  
not determined

Partition coefficient n-octanol/water:

not determined

Vapour pressure:  
(at 50 °C)

<=1100 hPa

Density:

not determined

Relative vapour density:

not determined

#### 9.2. Other information

##### Information with regard to physical hazard classes

Self-ignition temperature

Solid:

not applicable

Gas:

not applicable

Oxidizing properties

Not oxidising.

##### Other safety characteristics

Evaporation rate:

not determined

Solid content:

1,0 %

### SECTION 10: Stability and reactivity

#### 10.1. Reactivity

Highly flammable.

#### 10.2. Chemical stability

The product is stable under storage at normal ambient temperatures.

#### 10.3. Possibility of hazardous reactions

No known hazardous reactions.

#### 10.4. Conditions to avoid

Keep away from sources of heat (e.g. hot surfaces), sparks and open flames. Vapours can form explosive mixtures with air.

#### 10.5. Incompatible materials

No information available.

#### 10.6. Hazardous decomposition products

No known hazardous decomposition products.

### SECTION 11: Toxicological information

#### 11.1. Information on toxicological effects

##### Acute toxicity

Based on available data, the classification criteria are not met.

##### ATEmix calculated

ATE (oral) > 5000 mg/kg; ATE (dermal) > 5000 mg/kg; ATE (inhalation vapour) > 50 mg/l; ATE (inhalation dust/mist) > 12.5 mg/l

| CAS No  | Chemical name       |                   |         |        |        |
|---------|---------------------|-------------------|---------|--------|--------|
|         | Exposure route      | Dose              | Species | Source | Method |
| 80-62-6 | methyl methacrylate |                   |         |        |        |
|         | dermal              | LD50 > 5000 mg/kg |         |        |        |

## Safety Data Sheet

according to UN GHS (ST/SG/AC.10/11/Rev.9)

### VITA VIONIC BOND I

Revision date: 17.01.2024

Product code: 288

Page 6 of 8

#### Irritation and corrosivity

Causes skin irritation.

Serious eye damage/eye irritation: Based on available data, the classification criteria are not met.

#### Sensitising effects

May cause an allergic skin reaction. (methyl methacrylate; triethylene glycol dimethacrylate; dibenzoyl peroxide; benzoyl peroxide)

#### Carcinogenic/mutagenic/toxic effects for reproduction

Based on available data, the classification criteria are not met.

#### STOT-single exposure

May cause respiratory irritation. (methyl methacrylate)

#### STOT-repeated exposure

Based on available data, the classification criteria are not met.

#### Aspiration hazard

Based on available data, the classification criteria are not met.

#### Additional information on tests

The mixture is classified as hazardous according to regulation (EC) No 1272/2008 [CLP].

### SECTION 12: Ecological information

#### 12.1. Toxicity

The product is not: Ecotoxic.

#### 12.2. Persistence and degradability

The product has not been tested.

#### 12.3. Bioaccumulative potential

The product has not been tested.

#### 12.4. Mobility in soil

The product has not been tested.

#### 12.6. Endocrine disrupting properties

This product does not contain a substance that has endocrine disrupting properties with respect to non-target organisms as no components meets the criteria.

#### 12.7. Other adverse effects

No information available.

#### Further information

Avoid release to the environment.

### SECTION 13: Disposal considerations

#### 13.1. Waste treatment methods

##### Disposal recommendations

Do not allow to enter into surface water or drains. Dispose of waste according to applicable legislation.

##### Contaminated packaging

Handle contaminated packages in the same way as the substance itself.

### SECTION 14: Transport information

#### Marine transport (IMDG)

|  |  |
|--|--|
| <b>14.1. UN number or ID number:</b>     | UN 1993  |
| <b>14.2. UN proper shipping name:</b>    | FLAMMABLE LIQUID, N.O.S. (methyl methacrylate) |
| <b>14.3. Transport hazard class(es):</b> | 3  |
| <b>14.4. Packing group:</b>              | II   |
| Hazard label:                            | 3  |

## Safety Data Sheet

according to UN GHS (ST/SG/AC.10/11/Rev.9)

### VITA VIONIC BOND I

Revision date: 17.01.2024

Product code: 288

Page 7 of 8



Special Provisions: 274  
 Limited quantity: 1 L  
 Excepted quantity: E2  
 EmS: F-E, S-E  
 Segregation group: ammonium compounds

#### Air transport (ICAO-TI/IATA-DGR)

**14.1. UN number or ID number:** UN 1993  
**14.2. UN proper shipping name:** FLAMMABLE LIQUID, N.O.S. (methyl methacrylate)  
**14.3. Transport hazard class(es):** 3  
**14.4. Packing group:** II  
 Hazard label: 3



Special Provisions: A3  
 Limited quantity Passenger: 1 L  
 Passenger LQ: Y341  
 Excepted quantity: E2  
 IATA-packing instructions - Passenger: 353  
 IATA-max. quantity - Passenger: 5 L  
 IATA-packing instructions - Cargo: 364  
 IATA-max. quantity - Cargo: 60 L

#### 14.5. Environmental hazards

ENVIRONMENTALLY HAZARDOUS: No

#### 14.6. Special precautions for user

Warning: Combustible liquid.

#### 14.7. Maritime transport in bulk according to IMO instruments

not applicable

### SECTION 15: Regulatory information

#### National regulatory information

Employment restrictions: Observe restrictions to employment for juveniles according to the 'juvenile work protection guideline' (94/33/EC).  
 Skin resorption/Sensitization: Causes allergic hypersensitivity reactions.

### SECTION 16: Other information

#### Changes

This data sheet contains changes from the previous version in section(s): 2.

## Safety Data Sheet

according to UN GHS (ST/SG/AC.10/11/Rev.9)

### VITA VIONIC BOND I

Revision date: 17.01.2024

Product code: 288

Page 8 of 8

#### Abbreviations and acronyms

ADR: Accord européen sur le transport des marchandises dangereuses par Route  
(European Agreement concerning the International Carriage of Dangerous Goods by Road)  
IMDG: International Maritime Code for Dangerous Goods  
IATA: International Air Transport Association  
GHS: Globally Harmonized System of Classification and Labelling of Chemicals  
EINECS: European Inventory of Existing Commercial Chemical Substances  
ELINCS: European List of Notified Chemical Substances  
CAS: Chemical Abstracts Service  
LC50: Lethal concentration, 50%  
LD50: Lethal dose, 50%

#### Further Information

The above information describes exclusively the safety requirements of the product and is based on our present-day knowledge. The information is intended to give you advice about the safe handling of the product named in this safety data sheet, for storage, processing, transport and disposal. The information cannot be transferred to other products. In the case of mixing the product with other products or in the case of processing, the information on this safety data sheet is not necessarily valid for the new made-up material.

---

*(The data for the hazardous ingredients were taken respectively from the last version of the sub-contractor's safety data sheet.)*