


## VITASONIC II

Ultrasonic mixer for the VITA In-Ceram® slip technique



**VITA**

### Product description

High-frequency ultrasonic mixer specially developed for the VITA In-Ceram slip technique

### Area of application

To mix VITA In-Ceram powder slip

## VITASONIC II

### Technical functions

Timer allows both short and continuous operation mode

### Supplied with

#### Prod. no. DSON220

Mixer, with: 1 mains cable,  
1 lid with insert, 1 operating manual

### Printed materials

– Operating manual no. 481M

### Technical data

Device width	220 mm
Device depth	95 mm
Device height	170 mm
Container width	190 mm
Container depth	85 mm
Container height	60 mm
Capacity	0.8 l
Weight	2.2 kg
Casing	Stainless steel
Power supply	230 volt AC, 50 Hz or 100/110 volt AC, 50/60 Hz
Power consumption	max. 80 watt
Frequency	25 KHz

### Features / Benefits

- Lid with special insert for the VITA In-Ceram slip technique
- Can also be used as an ultrasonic cleaning device if lid is removed
- Robust stainless steel casing
- Timer allows both short and continuous operation mode