

VITA ENAMIC® multiColor – shade gradient included

The multi-chromatic CAD/CAM blanks for ideal anterior esthetics



VITA – perfect match

VITA

Your benefits with VITA ENAMIC multiColor:

- Maximum time efficiency and esthetics
- Minimally invasive restorations
- Reliable results
- Simplified processing and placement precision offer maximum cost effectiveness
- Available in 5 VITA SYSTEM 3D-MASTER shades (1M1-HT, 1M2-HT, 2M2-HT, 3M2-HT, 4M2-HT)
- Available in 2 sizes (EMC-14 und EMC-16)

VITA ENAMIC® multiColor – shade gradient included

The multi-chromatic CAD/CAM blanks for ideal anterior esthetics

APPLICATION OPTIONS



Benefits



Enables cost-effective restorations

Time efficiency. Truly natural esthetic results. Minimally invasive restorations. Long-lasting results. There are numerous requirements for successful anterior restorations. With VITA ENAMIC multiColor, they are perfectly fulfilled in a single block.

Maximum time efficiency and esthetics

VITA ENAMIC multiColor's integrated, natural shade gradient includes six finely nuanced layers with a smooth transition from cervical to incisal edge. Firing is completely omitted. Staining or characterizing are no longer required. Just polish – that's it!

Minimally invasive restorations

The dual ceramic polymer network structure of VITA ENAMIC multiColor can be milled very smoothly and as thin as 0.3 mm, enabling particularly low wall thickness and thinly tapering edges. Suitable in cases of limited space availability.

Reliable results

VITA ENAMIC is a material with integrated buffer and crack-stop function. Like a shock absorber, the hybrid ceramic absorbs all forces acting on it.

Adhesive bonding

Thanks to the dominant ceramic network, VITA ENAMIC can be etched and silanized in the usual way with hydrofluoric acid. As a result, excellent bond strength values are achieved during adhesive bonding.

VITA