

# Safety Data Sheet

according to WHMIS

## VITA VM OPAQUE FLUID

Revision date:

Product code: 168

Page 1 of 7

### 1. Identification

#### Product identifier

VITA VM OPAQUE FLUID

#### Relevant identified uses of the substance or mixture and uses advised against

##### Use of the substance/mixture

Use as laboratory reagent

#### Details of the supplier of the safety data sheet

Company name:	VITA Zahnfabrik H.Rauter GmbH & Co.KG	
Street:	Spitalgasse 3	
Place:	D-79713 Bad Säckingen	
Post-office box:	1338	
	D-79704 Bad Säckingen	
Telephone:	+49(0)7761-562-0	Telefax: +49(0)7761-562-299
e-mail:	info@vita-zahnfabrik.com	
Contact person:	regulatory affairs	
e-mail:	info@vita-zahnfabrik.com	
Internet:	www.vita-zahnfabrik.com	
Responsible Department:	Regulatory Affairs	

**Emergency telephone number:** +49-(0)761-19240

#### Further Information

medical device

### 2. Hazard identification

#### Classification of the substance or mixture

##### WHMIS 2015

Skin corrosion/irritation: Skin Corr. 1

Serious eye damage/eye irritation: Eye Dam. 1

#### Label elements

##### WHMIS 2015

Signal word: Danger

Pictograms:



#### Hazard statements

Causes severe skin burns and eye damage.

#### Precautionary statements

Do not breathe dust/fume/gas/mist/vapours/spray.

Wear protective gloves/protective clothing/eye protection/face protection.

IF SWALLOWED: Rinse mouth. Do NOT induce vomiting.

IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water or shower.

IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do.

Continue rinsing.

#### Other hazards

No information available.

### 3. Composition/information on ingredients

**Safety Data Sheet**

according to WHMIS

**VITA VM OPAQUE FLUID**

Revision date:

Product code: 168

Page 2 of 7

**Mixtures****Hazardous components**

CAS No	Chemical name	Quantity
64-17-5	ethyl alcohol	1 - < 5% (*)
7790-69-4	Lithium nitrate	0.1 - < 1% (*)

(\*) The actual concentration is withheld as a trade secret.

**4. First-aid measures****Description of first aid measures****General information**

First aider: Pay attention to self-protection! Remove affected person from the danger area and lay down.

**After inhalation**

Provide fresh air. Medical treatment necessary.

**After contact with skin**

After contact with skin, wash immediately with plenty of water and soap. Take off immediately all contaminated clothing and wash it before reuse. If skin irritation occurs: Get medical advice/attention.

**After contact with eyes**

In case of contact with eyes flush immediately with plenty of flowing water for 10 to 15 minutes holding eyelids apart and consult an ophthalmologist.

**After ingestion**

Observe risk of aspiration if vomiting occurs. Rinse mouth immediately and drink plenty of water. Do NOT induce vomiting. Adverse human health effects and symptoms: Gastric perforation. Call a physician immediately. Do not allow a neutralisation agent to be drunk.

**Most important symptoms and effects, whether acute or delayed**

No information available.

**Indication of immediate medical attention and special treatment needed**

Treat symptomatically.

**5. Fire-fighting measures****Extinguishing media****Suitable extinguishing media**

Co-ordinate fire-fighting measures to the fire surroundings.

**Specific hazards arising from the hazardous product**

Non-flammable.

**Special protective equipment and precautions for fire-fighters**

Wear a self-contained breathing apparatus and chemical protective clothing. Full protection suit.

**Additional information**

Suppress gases/vapours/mists with water spray jet. Collect contaminated fire extinguishing water separately. Do not allow entering drains or surface water.

**6. Accidental release measures****Personal precautions, protective equipment and emergency procedures****General measures**

Provide adequate ventilation. Do not breathe gas/fumes/vapour/spray. Avoid contact with skin, eyes and clothes. Use personal protection equipment.

**Environmental precautions**

Do not allow to enter into surface water or drains.

**Safety Data Sheet**

according to WHMIS

**VITA VM OPAQUE FLUID**

Revision date:

Product code: 168

Page 3 of 7

**Methods and material for containment and cleaning up****Other information**

Absorb with liquid-binding material (sand, diatomaceous earth, acid- or universal binding agents). Treat the recovered material as prescribed in the section on waste disposal.

**Reference to other sections**

Safe handling: see section 7

Personal protection equipment: see section 8

Disposal: see section 13

**7. Handling and storage****Precautions for safe handling****Advice on safe handling**

If handled uncovered, arrangements with local exhaust ventilation have to be used. Do not breathe gas/fumes/vapour/spray.

**Advice on protection against fire and explosion**

No special fire protection measures are necessary.

**Conditions for safe storage, including any incompatibilities****Requirements for storage rooms and vessels**

Keep container tightly closed. Keep locked up. Store in a place accessible by authorized persons only. Provide adequate ventilation as well as local exhaust at critical locations.

**Hints on joint storage**

No special measures are necessary.

**8. Exposure controls/Personal protection****Control parameters****Exposure controls****Appropriate engineering controls**

If handled uncovered, arrangements with local exhaust ventilation have to be used. Do not breathe gas/fumes/vapour/spray.

**Protective and hygiene measures**

Remove contaminated, saturated clothing immediately. Draw up and observe skin protection programme.

Wash hands and face before breaks and after work and take a shower if necessary. When using do not eat, drink, smoke, sniff.

**Eye/face protection**

Suitable eye protection: goggles.

**Hand protection**

When handling with chemical substances, protective gloves must be worn with the CE-label including the four control digits. The quality of the protective gloves resistant to chemicals must be chosen as a function of the specific working place concentration and quantity of hazardous substances. For special purposes, it is recommended to check the resistance to chemicals of the protective gloves mentioned above together with the supplier of these gloves. Recommended glove articles KCL Dermatril P Breakthrough time: 480 NBR (Nitrile rubber)

**Skin protection**

Use of protective clothing.

**Safety Data Sheet**

according to WHMIS

**VITA VM OPAQUE FLUID**

Revision date:

Product code: 168

Page 4 of 7

**Respiratory protection**

Provide adequate ventilation as well as local exhaust at critical locations. Technical ventilation of workplace

**9. Physical and chemical properties****Information on basic physical and chemical properties**

Physical state:	Liquid	
Colour:	colourless	
Odour:	characteristic	
pH-Value:		12,7

**Changes in the physical state**

Melting point:		not determined
Boiling point or initial boiling point and boiling range:		100 °C
Flash point:		not determined

**Flammability**

Solid/liquid:		not applicable
Gas:		not applicable

**Explosive properties**

The product is not: Explosive.

Lower explosive limits:		not determined
Upper explosive limits:		not determined

**Self-ignition temperature**

Solid:		not applicable
Gas:		not applicable

Decomposition temperature:		not determined
----------------------------	--	----------------

**Oxidizing properties**

Not oxidising.

Vapour pressure: (at 50 °C)		<=1100 hPa
--------------------------------	--	------------

Density:		1,00000 g/cm <sup>3</sup>
----------	--	---------------------------

**Solubility in other solvents**

not determined

Partition coefficient n-octanol/water:		not determined
--	--	----------------

Relative vapour density:		not determined
--------------------------	--	----------------

Evaporation rate:		not determined
-------------------	--	----------------

**Other information**

Solid content:		0,0 %
----------------	--	-------

**10. Stability and reactivity****Reactivity**

Possibility of hazardous reactions.

**Chemical stability**

The product is stable under storage at normal ambient temperatures.

**Possibility of hazardous reactions**

Exothermic reaction with: Acid, Peroxides, Oxidizing agent.

**Safety Data Sheet**

according to WHMIS

**VITA VM OPAQUE FLUID**

Revision date:

Product code: 168

Page 5 of 7

**Conditions to avoid**

none

**Incompatible materials**

Keep away from: Acid, Oxidizing agent, Peroxides.

**Hazardous decomposition products**

No known hazardous decomposition products.

**11. Toxicological information****Information on toxicological effects****Acute toxicity**

Based on available data, the classification criteria are not met.

CAS No	Chemical name				
	Route of exposure	Dose	Species	Source	Method
64-17-5	ethyl alcohol				
	dermal	LD50 mg/kg	7060		

**Irritation and corrosivity**

Causes severe skin burns and eye damage.

Causes serious eye damage.

**Sensitizing effects**

Based on available data, the classification criteria are not met.

**Carcinogenic/mutagenic/toxic effects for reproduction**

Based on available data, the classification criteria are not met.

**STOT-single exposure**

Based on available data, the classification criteria are not met.

**STOT-repeated exposure**

Based on available data, the classification criteria are not met.

**Aspiration hazard**

Based on available data, the classification criteria are not met.

**Additional information on tests**

The mixture is classified as hazardous according to regulation (EC) No 1272/2008 [CLP].

**12. Ecological information****Ecotoxicity**

The product is not: Ecotoxic.

**Persistence and degradability**

The product has not been tested.

**Bioaccumulative potential**

The product has not been tested.

**Mobility in soil**

The product has not been tested.

**Other adverse effects**

No information available.

**Further information**

Avoid release to the environment.

**13. Disposal considerations**

**Safety Data Sheet**

according to WHMIS

**VITA VM OPAQUE FLUID**

Revision date:

Product code: 168

Page 6 of 7


**Waste treatment methods****Disposal recommendations**

Do not allow to enter into surface water or drains. Dispose of waste according to applicable legislation.


**Contaminated packaging**

Wash with plenty of water. Completely emptied packages can be recycled.

**14. Transport information****Marine transport (IMDG)**

<b><u>UN number:</u></b>	UN 1824
<b><u>United Nations proper shipping name:</u></b>	SODIUM HYDROXIDE SOLUTION
<b><u>Transport hazard class(es):</u></b>	8
<b><u>Packing group:</u></b>	III
Hazard label:	8
	
Special Provisions:	223
Limited quantity:	5 L
Excepted quantity:	E1
EmS:	F-A, S-B
Segregation group:	18 - alkalis

**Air transport (ICAO-TI/IATA-DGR)**

<b><u>UN number:</u></b>	UN 1824
<b><u>United Nations proper shipping name:</u></b>	SODIUM HYDROXIDE SOLUTION
<b><u>Transport hazard class(es):</u></b>	8
<b><u>Packing group:</u></b>	III
Hazard label:	8
	
Special Provisions:	A3 A803
Limited quantity Passenger:	1 L
Passenger LQ:	Y841
Excepted quantity:	E1
IATA-packing instructions - Passenger:	852
IATA-max. quantity - Passenger:	5 L
IATA-packing instructions - Cargo:	856
IATA-max. quantity - Cargo:	60 L

**Environmental hazards**

ENVIRONMENTALLY HAZARDOUS: No

**15. Regulatory information****Canadian regulations****16. Other information**

**Safety Data Sheet**

according to WHMIS

**VITA VM OPAQUE FLUID**

Revision date:

Product code: 168

Page 7 of 7

**Abbreviations and acronyms**

ADR: Accord européen sur le transport des marchandises dangereuses par Route  
(European Agreement concerning the International Carriage of Dangerous Goods by Road )  
IMDG: International Maritime Code for Dangerous Goods  
IATA: International Air Transport Association  
GHS: Globally Harmonized System of Classification and Labelling of Chemicals  
EINECS: European Inventory of Existing Commercial Chemical Substances  
ELINCS: European List of Notified Chemical Substances  
CAS: Chemical Abstracts Service  
LC50: Lethal concentration, 50%  
LD50: Lethal dose, 50%  
CLP: Classification, labelling and Packaging  
REACH: Registration, Evaluation and Authorization of Chemicals  
GHS: Globally Harmonised System of Classification, Labelling and Packaging of Chemicals  
UN: United Nations  
DNEL: Derived No Effect Level  
DMEL: Derived Minimal Effect Level  
PNEC: Predicted No Effect Concentration  
ATE: Acute toxicity estimate  
LL50: Lethal loading, 50%  
EL50: Effect loading, 50%  
EC50: Effective Concentration 50%  
ErC50: Effective Concentration 50%, growth rate  
NOEC: No Observed Effect Concentration  
BCF: Bio-concentration factor  
PBT: persistent, bioaccumulative, toxic  
vPvB: very persistent, very bioaccumulative  
RID: Regulations concerning the international carriage of dangerous goods by rail  
ADN: European Agreement concerning the International Carriage of Dangerous Goods by Inland Waterways  
(Accord européen relatif au transport international des marchandises dangereuses par voies de navigation  
intérieures)  
EmS: Emergency Schedules  
MFAG: Medical First Aid Guide  
ICAO: International Civil Aviation Organization  
MARPOL: International Convention for the Prevention of Marine Pollution from Ships  
IBC: Intermediate Bulk Container  
SVHC: Substance of Very High Concern  
For abbreviations and acronyms, see table at <http://abbrev.esdscom.eu>

**Further Information**

The information is based on the present level of our knowledge. It does not, however, give assurance of product properties and establishes no contract legal rights. The receiver of our product is singularly responsible for adhering to existing laws and regulations.

---

*(The data for the hazardous ingredients were taken respectively from the last version of the sub-contractor's safety data sheet.)*