

**Safety Data Sheet**

according to 29 CFR 1910.1200(g)

**VITA AKZENT LC CLEANER**

Revision date:

Product code: 328

Page 1 of 8

**1. Identification****Product identifier**

VITA AKZENT LC CLEANER

**Recommended use of the chemical and restrictions on use****Use of the substance/mixture**

Use as laboratory reagent

**Details of the supplier of the safety data sheet**

Company name:	VITA Zahnfabrik H.Rauter GmbH & Co.KG	
Street:	Spitalgasse 3	
Place:	D-79713 Bad Säckingen	
Post-office box:	1338	
	D-79704 Bad Säckingen	
Telephone:	+49(0)7761-562-0	Telefax: +49(0)7761-562-299
e-mail:	info@vita-zahnfabrik.com	
Contact person:	regulatory affairs	
e-mail:	info@vita-zahnfabrik.com	
Internet:	www.vita-zahnfabrik.com	
Responsible Department:	Regulatory Affairs	
<b><u>Emergency phone number:</u></b>	+49-(0)761-19240	

**2. Hazard(s) identification****Classification of the chemical****29 CFR Part 1910.1200**

Flammable liquids: Flam. Liq. 2

Serious eye damage/eye irritation: Eye Irrit. 2A

**Label elements****29 CFR Part 1910.1200****Signal word:** Danger**Pictograms:****Hazard statements**

Highly flammable liquid and vapor

Causes serious eye irritation

**Precautionary statements**

Keep away from heat/sparks/open flames/hot surfaces. - No smoking.

Keep container tightly closed.

Ground/bond container and receiving equipment.

Wear protective gloves/protective clothing/eye protection/face protection.

**Hazards not otherwise classified**

No information available.

**3. Composition/information on ingredients****Mixtures**

**Safety Data Sheet**

according to 29 CFR 1910.1200(g)

**VITA AKZENT LC CLEANER**

Revision date:

Product code: 328

Page 2 of 8

**Hazardous components**

CAS No	Components	Quantity
64-17-5	ethyl alcohol	90 %
78-93-3	butanone; ethyl methyl ketone	1 %

**4. First-aid measures****Description of first aid measures****After inhalation**

Provide fresh air.

**After contact with skin**

Wash with plenty of water. Take off contaminated clothing and wash it before reuse.

**After contact with eyes**

After contact with the eyes, rinse with water with the eyelids open for a sufficient length of time, then consult an ophthalmologist immediately.

**After ingestion**

Rinse mouth immediately and drink plenty of water.

**Most important symptoms and effects, both acute and delayed**

No information available.

**Indication of any immediate medical attention and special treatment needed**

Treat symptomatically.

**5. Fire-fighting measures****Extinguishing media****Suitable extinguishing media**Carbon dioxide (CO<sub>2</sub>), Foam, Extinguishing powder.**Unsuitable extinguishing media**

Water.

**Specific hazards arising from the chemical**

Highly flammable. Vapors may form explosive mixtures with air.

**Special protective equipment and precautions for fire-fighters**

In case of fire: Wear self-contained breathing apparatus.

**Additional information**

Use water spray/stream to protect personnel and to cool endangered containers. Suppress gases/vapors/mists with water spray jet. Collect contaminated fire extinguishing water separately. Do not allow entering drains or surface water.

**6. Accidental release measures****Personal precautions, protective equipment and emergency procedures****General measures**

Remove all sources of ignition. Do not breathe gas/fume/vapor/spray. Avoid contact with skin, eyes and clothes. Use personal protection equipment.

**Environmental precautions**

Do not allow uncontrolled discharge of product into the environment. Explosion risk.

**Methods and material for containment and cleaning up****Other information**

Absorb with liquid-binding material (sand, diatomaceous earth, acid- or universal binding agents). Treat the recovered material as prescribed in the section on waste disposal.

**Safety Data Sheet**

according to 29 CFR 1910.1200(g)

**VITA AKZENT LC CLEANER**

Revision date:

Product code: 328

Page 3 of 8

**Reference to other sections**

Safe handling: see section 7

Personal protection equipment (PPE): see section 8

Disposal: see section 13

**7. Handling and storage****Precautions for safe handling****Advice on safe handling**

No special measures are necessary.

**Advice on protection against fire and explosion**

Keep away from sources of ignition - No smoking. Take precautionary measures against static discharges.

Vapors may form explosive mixtures with air.

**Conditions for safe storage, including any incompatibilities****Requirements for storage rooms and vessels**

Keep container tightly closed. Keep in a cool, well-ventilated place. Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking.

**Hints on joint storage**

Do not store together with: Oxidizing agent. Pyrophoric or self-heating substances.

**8. Exposure controls/personal protection****Control parameters****Exposure limits**

CAS No.	Substance	ppm	mg/m <sup>3</sup>	f/cc	Category	Origin
78-93-3	2-Butanone (Methyl ethyl ketone)	200	590		TWA (8 h)	PEL
78-93-3	2-Butanone	200	590		TWA (8 h)	REL
		300	885		STEL (15 min)	REL
64-17-5	Ethyl alcohol (Ethanol)	1000	1900		TWA (8 h)	PEL
64-17-5	Ethyl alcohol	1000	1900		TWA (8 h)	REL

**Exposure controls****Protective and hygiene measures**

Remove contaminated, saturated clothing immediately. Draw up and observe skin protection programme.

Wash hands and face before breaks and after work and take a shower if necessary. When using do not eat, drink, smoke, sniff.

**Eye/face protection**

Wear eye/face protection.

**Hand protection**

When handling with chemical substances, protective gloves must be worn with the CE-label including the four control digits. The quality of the protective gloves resistant to chemicals must be chosen as a function of the specific working place concentration and quantity of hazardous substances. For special purposes, it is recommended to check the resistance to chemicals of the protective gloves mentioned above together with the supplier of these gloves. Recommended glove articles KCL Dermatrill P Breakthrough time: 30 min NBR (Nitrile rubber)

**Safety Data Sheet**

according to 29 CFR 1910.1200(g)

**VITA AKZENT LC CLEANER**

Revision date:

Product code: 328

Page 4 of 8

**Skin protection**

Flame-retardant protective clothing. Wear anti-static footwear and clothing. Wear suitable protective clothing.

**Respiratory protection**

In case of inadequate ventilation wear respiratory protection. Provide adequate ventilation as well as local exhaust at critical locations. Technical ventilation of workplace

**9. Physical and chemical properties****Information on basic physical and chemical properties**

Physical state:	Liquid	
Color:	colorless	
pH-Value:		not determined

**Changes in the physical state**

Melting point/freezing point:		not determined
Boiling point or initial boiling point and boiling range:		78 °C
Flash point:		13 °C

**Flammability**

Solid/liquid:		not applicable
Gas:		not applicable

**Explosive properties**

The product is not: Explosive.

Lower explosion limits:		3,5 vol. %
Upper explosion limits:		15 vol. %
Auto-ignition temperature:		425 °C

**Self-ignition temperature**

Solid:		not applicable
Gas:		not applicable

Decomposition temperature:		not determined
----------------------------	--	----------------

**Oxidizing properties**

Not oxidising.

Vapor pressure: (at 50 °C)		<=1100 hPa
-------------------------------	--	------------

Density:		0,80000 g/cm <sup>3</sup>
----------	--	---------------------------

**Solubility in other solvents**

not determined

Partition coefficient n-octanol/water:		not determined
----------------------------------------	--	----------------

Viscosity / dynamic: (at 20 °C)		1,2 mPa·s
------------------------------------	--	-----------

Relative vapour density:		not determined
--------------------------	--	----------------

Evaporation rate:		not determined
-------------------	--	----------------

**Other information**

Solid content:		0,0 %
----------------	--	-------

**10. Stability and reactivity****Reactivity**

Highly flammable.

**Safety Data Sheet**

according to 29 CFR 1910.1200(g)

**VITA AKZENT LC CLEANER**

Revision date:

Product code: 328

Page 5 of 8

**Chemical stability**

The product is stable under storage at normal ambient temperatures.

**Possibility of hazardous reactions**

No known hazardous reactions.

**Conditions to avoid**

Keep away from sources of heat (e.g. hot surfaces), sparks and open flames. Vapors may form explosive mixtures with air.

**Incompatible materials**

No information available.

**Hazardous decomposition products**

No known hazardous decomposition products.

**11. Toxicological information****Information on toxicological effects****Acute toxicity**

Based on available data, the classification criteria are not met.

CAS No	Components				
	Exposure route	Dose	Species	Source	Method
64-17-5	ethyl alcohol				
	dermal	LD50 mg/kg	7060		
78-93-3	butanone; ethyl methyl ketone				
	dermal	LD50 mg/kg	5000		

**Irritation and corrosivity**

Causes serious eye irritation

Skin corrosion/irritation: Based on available data, the classification criteria are not met.

**Sensitizing effects**

Based on available data, the classification criteria are not met.

**Carcinogenic/mutagenic/toxic effects for reproduction**

Based on available data, the classification criteria are not met.

**Specific target organ toxicity (STOT) - single exposure**

Based on available data, the classification criteria are not met.

**Specific target organ toxicity (STOT) - repeated exposure**

Based on available data, the classification criteria are not met.

Carcinogenicity (OSHA):

No ingredient of this mixture is listed.

Carcinogenicity (IARC):

Ethanol in alcoholic beverages (CAS 64-17-5) is listed in group 1.

Carcinogenicity (NTP):

No ingredient of this mixture is listed.

**Aspiration hazard**

Based on available data, the classification criteria are not met.

**Additional information on tests**

The mixture is classified as hazardous according to regulation (EC) No 1272/2008 [CLP].

**12. Ecological information****Ecotoxicity**

The product is not: Ecotoxic.

**Persistence and degradability**

**Safety Data Sheet**

according to 29 CFR 1910.1200(g)

**VITA AKZENT LC CLEANER**

Revision date:

Product code: 328

Page 6 of 8

The product has not been tested.

**Bioaccumulative potential**

The product has not been tested.

**Mobility in soil**

The product has not been tested.

**Other adverse effects**

No information available.

**Further information**

Avoid release to the environment.

**13. Disposal considerations****Waste treatment methods****Disposal recommendations**

Do not allow to enter into surface water or drains. Dispose of waste according to applicable legislation.

**Contaminated packaging**

Wash with plenty of water. Completely emptied packages can be recycled. Handle contaminated packages in the same way as the substance itself. Waste codes/waste designations according to EWC/AVV

**14. Transport information****Marine transport (IMDG)**

**UN number:** UN 1170  
**UN proper shipping name:** ETHANOL (ETHYL ALCOHOL)  
**Transport hazard class(es):** 3  
**Packing group:** II  
 Hazard label: 3



Special Provisions: 144  
 Limited quantity: 1 L  
 Excepted quantity: E2  
 EmS: F-E, S-D

**Air transport (ICAO-TI/IATA-DGR)**

**UN number:** UN 1170  
**UN proper shipping name:** ETHYL ALCOHOL  
**Transport hazard class(es):** 3  
**Packing group:** II  
 Hazard label: 3



Special Provisions: A3 A58 A180  
 Limited quantity Passenger: 1 L  
 Passenger LQ: Y341  
 Excepted quantity: E2  
 IATA-packing instructions - Passenger: 353  
 IATA-max. quantity - Passenger: 5 L

**Safety Data Sheet**

according to 29 CFR 1910.1200(g)

**VITA AKZENT LC CLEANER**

Revision date:

Product code: 328

Page 7 of 8

IATA-packing instructions - Cargo: 364

IATA-max. quantity - Cargo: 60 L

**Environmental hazards**

ENVIRONMENTALLY HAZARDOUS: No

**Special precautions for user**

Warning: Combustible liquid.

**Transport in bulk according to Annex II of MARPOL 73/78 and the IBC Code**

not applicable

**15. Regulatory information****U.S. Regulations****National regulatory information**

SARA Section 304 CERCLA:

Methyl ethyl ketone (78-93-3): Reportable quantity = 5,000 (2270) lbs. (kg)

SARA Section 311/312 Hazards:

ethyl alcohol (64-17-5): Fire hazard, Immediate (acute) health hazard

Methyl ethyl ketone (78-93-3): Fire hazard, Immediate (acute) health hazard

Clean Air Act Section 112(b):

Methyl ethyl ketone (78-93-3)

**State Regulations****Safe Drinking Water and Toxic Enforcement Act of 1986 (Proposition 65, State of California)**

This product can not expose you to chemicals known to the State of California to cause cancer, birth defects or other reproductive harm.

**16. Other information**

Revision date:

Revision No: 1

**Abbreviations and acronyms**ADR: Accord européen sur le transport des marchandises dangereuses par Route  
(European Agreement concerning the International Carriage of Dangerous Goods by Road)

IMDG: International Maritime Code for Dangerous Goods

IATA: International Air Transport Association

GHS: Globally Harmonized System of Classification and Labelling of Chemicals

EINECS: European Inventory of Existing Commercial Chemical Substances

ELINCS: European List of Notified Chemical Substances

CAS: Chemical Abstracts Service

LC50: Lethal concentration, 50%

LD50: Lethal dose, 50%

CLP: Classification, labelling and Packaging

REACH: Registration, Evaluation and Authorization of Chemicals

GHS: Globally Harmonised System of Classification, Labelling and Packaging of Chemicals

UN: United Nations

DNEL: Derived No Effect Level

DMEL: Derived Minimal Effect Level

PNEC: Predicted No Effect Concentration

ATE: Acute toxicity estimate

LL50: Lethal loading, 50%

EL50: Effect loading, 50%

EC50: Effective Concentration 50%

ErC50: Effective Concentration 50%, growth rate

NOEC: No Observed Effect Concentration

**Safety Data Sheet**

according to 29 CFR 1910.1200(g)

**VITA AKZENT LC CLEANER**

Revision date:

Product code: 328

Page 8 of 8

BCF: Bio-concentration factor

PBT: persistent, bioaccumulative, toxic

vPvB: very persistent, very bioaccumulative

RID: Regulations concerning the international carriage of dangerous goods by rail

ADN: European Agreement concerning the International Carriage of Dangerous Goods by Inland Waterways  
(Accord européen relatif au transport international des marchandises dangereuses par voies de navigation intérieures)

EmS: Emergency Schedules

MFAG: Medical First Aid Guide

ICAO: International Civil Aviation Organization

MARPOL: International Convention for the Prevention of Marine Pollution from Ships

IBC: Intermediate Bulk Container

SVHC: Substance of Very High Concern

For abbreviations and acronyms, see table at <http://abbrev.esdscom.eu>**Other data**

The information is based on the present level of our knowledge. It does not, however, give assurance of product properties and establishes no contract legal rights. The receiver of our product is singularly responsible for adhering to existing laws and regulations.

---

*(The data for the hazardous ingredients were taken respectively from the last version of the sub-contractor's safety data sheet.)*